

Cimilre Breast Pump Comparison

cimilre®



Cimilre S6 / S6+



Cimilre P1



Cimilre S5+



Cimilre E1

	Cimilre S6 / S6+	Cimilre P1	Cimilre S5+	Cimilre E1
Weight / Profile	1.71 / 1.9 lbs.	0.55 lbs.	0.9 lbs.	0.37 lbs.
Maximum Suction	-285 mmHg	-285 mmHg	-285 mmHg	-285mmHg
Dual-Motor Technology	No	No	Yes	No
Display	Backlit LCD	Backlit LCD	Backlit LCD	LED
Included breast shield sizes	21 & 24 mm*	21 & 24 mm*	21 & 24 mm*	21 & 24 mm*
Free alternative breast shields	16, 19, 28, 32 mm	16, 19, 28, 32 mm	16, 19, 28, 32 mm	16, 19, 28, 32 mm
Personalize with decorative skin	Yes	Yes	Yes	No
Battery (Rechargeable lithium-ion)	No / Yes	Yes	Yes	Yes
Night Light	3 levels of light	No	No	No
Warranty (Pump)	2-Year	2-Year	2-Year	2-year
Massage Mode	Suction Levels	1-5	1-5	1-5
	Cycle Speeds	Adjustable (60-80 cpm)	Fixed (49-70 cpm)	Fixed (50-86 cpm)
Expression Mode	Suction Levels	1-12	1-10	1-10 levels
	Cycle Speeds	Adjustable (30-46 cpm)	Fixed (25-56 cpm)	Fixed (24-60 cpm)

*Breast shield sizes included may vary during our transition from 24mm and 28mm to 21mm and 24mm.



Stay connected
@cimilrebreastpumps

cimilrebreastpumps.com

Contact us

sales@cimilrebreastpumps.com
help@cimilrebreastpumps.com