

## DiaFil Flow™

Light-Cured Flowable Composite

**Buy 1 Econo Pack,  
Get 1 DiaEtch Econo Pack  
FREE!**



- Perfect adaptation and handling
- Excellent color adaptability
- Pit & fissure sealant
- Excellent tensile strength and low shrinkage
- Highly radiopaque
- Class III, IV and V direct restorations
- Base/Liner under posterior restorations in Class I, II

### Econo Pack

Contains: 4 x 2g flowable syringe & 40 tips.  
Shades: A1, A2, A3, A3.5, B2, B3, C2, A20, A30.

## DiaTemp™ & DiaTemp™ FLOW

Light Cured Temporary Restorative Material

**Buy 2 of any,  
Get 1 Dia-ProSeal™ Intro Kit  
FREE!**

### FEATURES

- 1 DiaTemp: high-viscosity
- 2 DiaTemp Flow: low-viscosity
- 3 Great for implant access holes
- 4 Elastic & easy to handle
- 5 Simple and fast removal in one piece
- 6 Strong, durable, & odorless
- 7 Excellent adhesion



[A2003-1101] **BLUE** DiaTemp  
[A2003-1201] **YELLOW** DiaTemp  
Contains: 3 x 3 g syringe

[A2003-2101] **BLUE** DiaTemp Flow  
[A2003-2201] **YELLOW** DiaTemp Flow  
Contains: 5 x 1.8 g syringe and 10 tips

## DiaPlus™

One-Step Bonding Agent (5th gen)

**BUY 2, GET 1 FREE!**



- 1 Direct & indirect applications (5<sup>th</sup> gen)
- 2 Low film thickness
- 3 No adhesive show-through
- 4 Virtually no post-op sensitivity
- 5 Ethanol-based
- 6 Compatible with dual- & self-cured resin
- 7 Exceptional & reliable bond strength

[A2001-2101] DiaPlus Bond (5ml)

## Dia-Ultrasonic®

Ultrasonic Endo Activator

**COMING SOON!**

[4009-1001] Regular Kit

### FEATURES

- 1 Dual-frequency: 30 KHz and 45 KHz
- 2 Flexible and strong NiTi tip
- 3 Easy to use and compact
- 4 Fast charging: 3.5 hours for a full charge
- 5 Long lasting: up to 500 cycles
- 6 Powerful and effective irrigation



## DIA-PROSEAL™

Epoxy Resin-based Root Canal Sealer



**Buy 2 Regular Kits,  
Get 1 Dia-Prep Pro EDTA FREE!**

### FEATURES

- 1 Excellent sealing ability
- 2 Low shrinkage & solubility
- 3 Outstanding flowability
- 4 Highly antimicrobial with Ca(OH)<sub>2</sub>
- 5 Fast preparation w/ dual syringe
- 6 Setting time: 3-4 hours

[1003-201] Regular Kit  
16 g dual-syringe (Base 6.5 g + Catalyst 9.5 g),  
1 Mixing Pad and 1 Spatula