



DIRECT SYSTEM

VENEER

edelweiss
DENTISTRY

shaping the future of dentistry




Claudio Novelli
D.M.D PhD Private-Dozent
Clinical and Scientific Director

edelweiss VENEERs LONG

«I started using VENEERs LONG to restore my periodontal patients with long teeth due to gingival recession. However I soon realized that also many not-periodontal patients have teeth with increased length-to-width ratio (i.e. l/w ratio > 1,4:1) and are well suited for edelweiss VENEERs LONG. Now I like to always have VENEERs LONG in my edelweiss VENEERs stock because it makes the system more flexible and allows me to treat a wider range of esthetic patients. Regular edelweiss VENEERs are designed with ideal length-to-width ratio (i.e. 1,2:1) to fit the large majority of patients. However, sometimes patients have teeth that are not ideal in size and it is not unusual to see longer teeth with increased length-to-width ratio. This is a common clinical finding in periodontal patients where soft tissue recession produces root exposure and increased tooth length. edelweiss VENEERs LONG are designed to optimally fit periodontal patients and allow to cover the exposed roots and to restore natural tooth appearance with a single-appointment, minimally-invasive restoration.»

For more information visit this link:
<https://www.edelweissdentistry.com/press>



EDELWEISS VENEERs LONG



vs. EDELWEISS VENEERs



ADVANTAGES

- One appointment
- Minimally invasive
- Natural look
- Long-lasting
- Biomechanical
- Biocompatible
- Antibacterial
- Economical

VENEERs LONG increase the flexibility of the edelweiss system because now a larger variety of veneer sizes is available. Until now edelweiss VENEERs were available with regular length-to-width ratio 1,2:1. Now also VENEERs LONG with increased length-to-width ratio 1,4:1 are available and dentists can choose the veneer proportion that best fits their patients.

CASE STUDY - OPEN BITE CLOSER

with edelweiss VENEERS LONG



CASE STUDY - PERIODONTAL FUNCTION & ESTHETICS

with edelweiss VENEERS LONG



SIZES

REF 17036-37	edelweiss VENEER 11 LONG	S, M, L	REF 17065-67	edelweiss VENEER 11/21 LONG	S, M, L
REF 17038-40	edelweiss VENEER 12 LONG	S, M, L	REF 17068-70	edelweiss VENEER 12/22 LONG	S, M, L
REF 17041-43	edelweiss VENEER 13 LONG	S, M, L	REF 17071-73	edelweiss VENEER 13-23 LONG	S, M, L
REF 17050-52	edelweiss VENEER 21 LONG	S, M, L			
REF 17053-55	edelweiss VENEER 22 LONG	S, M, L			
REF 17056-58	edelweiss VENEER 23 LONG	S, M, L			

edelweiss
DENTISTRY

edelweiss dentistry products gmbh • Austria
office@edelweissdentistry.com

www.edelweissdentistry.com

#822021 / 2022-04



edelweiss dentistry ® and the edelweiss Logo ® are registered trademarks of edelweiss dentistry products gmbh • Austria

© 2022 edelweiss dentistry products gmbh • Austria. All material contained in this brochure are the property of edelweiss dentistry products gmbh • Austria. No part of this publication may be reproduced, distributed, or transmitted in any form or by any means, without the prior written permission of the edelweiss dentistry products gmbh • Austria. All Rights Reserved.