


Technical Aspects

DIRECT SYSTEM PEDIATRIC CROWN



Size S - Upper	mesial - distal A	vestibular - oral B	incisal - cervical C
E & F / 51 & 61	5.7	5.0	4.9
D & G / 52 & 62	5.1	4.5	4.1
C & H / 53 & 63	7.0	6.5	5.0
B & I / 54 & 64	7.1	8.2	3.8
A & J / 55 & 65	9.3	9.1	3.9

Size S - Lower	mesial - distal A	vestibular - oral B	incisal - cervical C
P & O / 81 & 71	4.7	4.2	3.8
Q & N / 82 & 72	4.8	4.5	4.2
R & M / 83 & 73	6.3	6.0	4.6
S & L / 84 & 74	8.3	6.6	4.3
T & K / 85 & 75	9.8	8.4	4.2

Size M - Upper	mesial - distal A	vestibular - oral B	incisal - cervical C
E & F / 51 & 61	6.4	5.6	5.4
D & G / 52 & 62	5.7	5.0	4.6
C & H / 53 & 63	7.7	7.4	5.5
B & I / 54 & 64	7.9	9.0	4.3
A & J / 55 & 65	10.3	10.1	4.4

Size M - Lower	mesial - distal A	vestibular - oral B	incisal - cervical C
P & O / 81 & 71	5.2	4.7	4.2
Q & N / 82 & 72	5.3	5.0	4.6
R & M / 83 & 73	7.0	6.5	5.0
S & L / 84 & 74	9.2	7.3	4.3
T & K / 85 & 75	10.8	9.4	4.5

Size L - Upper	mesial - distal A	vestibular - oral B	incisal - cervical C
E & F / 51 & 61	7.0	6.1	6.1
D & G / 52 & 62	6.3	5.5	5.0
C & H / 53 & 63	8.6	8.2	6.1
B & I / 54 & 64	8.7	9.9	4.8
A & J / 55 & 65	11.4	11.1	4.9

Size L - Lower	mesial - distal A	vestibular - oral B	incisal - cervical C
P & O / 81 & 71	5.7	5.2	4.6
Q & N / 82 & 72	5.8	5.5	5.1
R & M / 83 & 73	7.7	7.2	5.6
S & L / 84 & 74	10.3	8.1	5.3
T & K / 85 & 75	11.9	10.3	5.1

Size XL - Upper	mesial - distal A	vestibular - oral B	incisal - cervical C
E & F / 51 & 61	7.7	6.7	6.7
D & G / 52 & 62	6.9	6.1	5.5
C & H / 53 & 63	9.5	9.0	6.7
B & I / 54 & 64	9.6	10.9	5.3
A & J / 55 & 65	12.5	12.2	5.4

Size XL - Lower	mesial - distal A	vestibular - oral B	incisal - cervical C
P & O / 81 & 71	6.3	5.7	5.1
Q & N / 82 & 72	6.4	6.1	5.6
R & M / 83 & 73	8.5	7.9	6.2
S & L / 84 & 74	11.3	8.9	5.8
T & K / 85 & 75	13.1	11.3	5.6