

NucleoVac Vacuum Regulator


1 Contents

NucleoVac Vacuum Regulator	
REF	740641
740641	1

2 How to connect the NucleoVac Vacuum Regulator

Compatible Tubing: 8 mm (inner diameter)

NucleoVac Vacuum Regulator



*Vacuum can cause implosions. Please make sure to use appropriate vessels. When working with vacuum use caution and wear safety glasses near a vacuum manifold under pressure.

3 Service Bioanalysis

Technical support and customer service:

MACHEREY-NAGEL Germany
 Tel.: +49 (0) 24 21 969 270
 e-mail: TECH-BIO@mn-net.com