

Receiver Columns 35 µm

1 Contents

Receiver Columns 35 µm			
REF	740524.10	740524.50	740524.250
Receiver Columns 35 µm	10	50	250

2 Product description

Receiver Columns are mini spin columns with an inserted hydrophobic frit of 35 µm pore size. They can be used for general filtration purposes as well as for retaining chromatographic resins (e.g., Sephadex® G25, G50, or Sephacryl® S200). Receiver Columns 35 µm are delivered with capped outlets, inserted into a Collection Tube (2 mL). Receiver Columns 35 µm are for use with suitable bench-top centrifuges.

3 Specifications

Column material:	Polypropylene
Filter frit material:	Polyethylene
Filter frit pore size:	Approx. 35 µm
Filter frit diameter:	6.8 mm
Effective filtration area:	36 mm ²
Max. volume:	0.7 mL
Outlet connection:	Fits in female luer

4 Product use restriction / warranty

Attention! Receiver Columns 35 µm are developed, designed, and sold for research purpose only.

All goods have to be checked immediately upon receipt. Damage claims are only acceptable in writing within 8 working days of receipt of all the goods. In case of legitimate claims the buyer can only require replacement of the goods. If replacement is not possible, the buyer has the right to choose between alternative products with same value or refund. The buyer cannot claim further compensation. All returns must first be authorized by us in writing. We offer a oneyear warranty from date of delivery that our products will conform to applicable specifications set forth in the product specifications if not sold to persons set forth in § 13 of German Civil Code. In such a case, the provisions of the German Civil Code shall be valid.

This warranty does not cover defects in or damage to the products which are due to improper installation or maintenance, misuse, neglect or any use other than ordinary commercial application.

Any discharge from liability will be void if a defect results from a heavily negligent or intentional breach of contract or if the buyer will be bodily injured as a result from a negligent or intentional breach.

The seller shall not be liable for the goods being fit for a particular purpose unless otherwise agreed upon, to which the buyer intends to put them. This warranty is strictly exclusive. Any further damage compensation is impossible..

Trademarks:

Sephadex® and Sephacryl® are registered trademarks of GE Healthcare companies

MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany
Tel.: +49 24 21 969-0 · Fax: +49 24 21 969-199 · tech-bio@mn-net.com · www.mn-net.com