

Regular Cut Offs for Drug Tests

Abbreviation	Description	Cut-off Level
AMP	Amphetamine	1000ng/mL
BAR	Barbiturates	300ng/mL
BUP	Buprenorphine	10ng/mL
BZO	Benzodiazepines	300ng/mL
COC	Cocaine	300ng/mL
mAMP	Methamphetamine	1000ng/mL
MDMA	MDMA Ecstasy	500ng/mL
MOP	Morphine	300ng/mL
MTD	Methadone	300ng/mL
OPI	Opiate	2000ng/mL
OXY	Oxycodone	100ng/mL
PCP	Phencyclidine	25ng/mL
PPX	Propoxyphene	300ng/mL
TCA	Tricyclic Antidepressants	1000ng/mL
THC	Marijuana/Cannabinoids	50ng/mL
TRA	Tramadol	200ng/mL
EDDP	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine	300ng/mL
K2	JWH-018 / JWH-073	50ng/mL / 25ng/mL

Drug Test Available

510K Cleared, CLIA waived and CE Marked:

AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OPI, OXY, PCP, PPX, TCA, THC

For Forensic Use only: TRA, EDDP, K2

Build-in Adulteration parameters offering visual check for urine integrity:

pH, Specific Gravity, Creatinine, Nitrite, Oxidant, Glutaraldehyde, Bleach, Pyridinium Chlorchromate

Special Cut Offs for Drug Tests

Abbreviation	Description	Cut-off Level
AMP	Amphetamine	300ng/mL
		500ng/mL
BZO	Benzodiazepines	200ng/mL
COC	Cocaine	100ng/mL
		150ng/mL
mAMP	Methamphetamine	300ng/mL
		500ng/mL
MOP	Morphine	100ng/mL
MTD	Methadone	200ng/mL
OXY	Oxycodone	300ng/mL
THC	Marijuana/Cannabinoids	20ng/mL
		40ng/mL
EDDP	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine	100ng/mL
TRA	Tramadol	100ng/mL
		1000ng/mL

Please contact us for ordering requirements of Special Cut Offs.

Drug Test Available

510K Cleared and CE Marked:

AMP (300ng/mL), BZO(200ng/mL), COC(100ng/mL), mAMP(500ng/mL), mAMP(300ng/mL), MOP(100ng/mL), MTD(200ng/mL), THC(40ng/mL)

For Forensic Use only:

AMP (500ng/mL), COC(150ng/mL), OXY(300ng/mL), EDDP(100ng/mL), THC(20ng/mL), TRA(100ng/mL, 1000ng/mL)

Build-in Adulteration parameters offering visual check for urine integrity:

pH, Specific Gravity, Creatinine, Nitrite, Oxidant, Glutaraldehyde, Bleach, Pyridinium Chlorochromate