

ANALYTICAL LABORATORIES (SINGAPORE) PTE LTD

Analytical Chemists; Environmental and Materials Testing
8 Kaki Bukit Place, Singapore 416186 Tel: 67460886 Fax: 67463830 Email: admin@analabs.com.sg
CO. REG NO. 197302347G GST REG NO. M2-0017430-5





REPORT

Lab No

: AC/MC/8306/15

Company Name

: Biocomm Private Limited

Date Received

: 02/12/2015

Date Reported: 07/12/2015

Sample Description

: One sample of mouthrinse

Date Tested: 02/12/2015-05/12/2015

The sample consisted of 2 x 100ml plastic bottles of mouthrinse marked as shown below.

On microbiological examination in accordance with the United States Pharmacopela 37, the following result was obtained:

Pearlle White Fluorinze Anti-bacterial Fluoride Mouthrinse (Membrane Filtration Method)

B/N: PFN 0142
Date: 01/12/2015

Total Aerobic Microbial Count (Membrane Filtration Method)

(SCDA, 35°C/3 days)

<10 CFU/ml</p>

Notes: 1) CFU - Colony Forming Units

2) SCDA - Soybean-Casein Digest Agar

3) < - Less than

Cheong Teck Geok Testing Officer

Phang Ken Aun Manager

ANALYTICAL LABORATORIES (S) PTE, LTD.

Manufactured by Biocomm Pte Ltd 53 Ubi Ave 1 #05-55 Paya Ubi Industrial Park Singapore 408934. www.biocomm.sg

/jaz



ANALYTICAL LABORATORIES (SINGAPORE) PTE LTD

Analytical Chemists; Environmental and Materials Testing 8 Kaki Bukit Place, Singapore 416186 Tel: 67460886 Fax: 67463830 Email: admin@analabs.com.sg CO. REG NO. 197302347G GST REG NO. M2-0017430-5





REPORT

Lab No

: AC/MC/8375/15

Company Name

: Biocomm Private Limited

Date Received

04/12/2015

Date Reported: 07/12/2015

Sample Description

: One sample of mouthrinse

Date Tested: 04/12/2015-07/12/2015

The sample consisted of 2 x 100ml plastic bottles of mouthrinse marked as shown below.

On microbiological examination in accordance with the United States Pharmacopeia 37, the following result was obtained:

Pearlie White Fluorinze Anti-bacterlal Fluoride Mouthrinse (100ml)	Total Aerobic Microbial Count (Membrane Filtration Method) (SCDA, 35°C/3 days)
B/N: PFN 0143 Date: 03/12/2015	Not detected (<10 CFU/ml)

Notes: 1) CFU - Colony Forming Units

2) SCDA - Soybean-Casein Digest Agar
3) < -1 esether

Cheong Teck Geok Testing Officer

Phang Ken Aun Manager

ANALYTICAL LABORATORIES (S) PTE. LTD.

Manufactured by Biocomm Pte Ltd 53 Ubl Ave 1 #05-55 Paya Ubi Industrial Park Singapore 408934. www.biocomm.sg

/jaz