TTAPE ACHILLES TENDONITIS



BEFORE YOU START

FOR BEST RESULTS



APPLY ONE HOUR **BEFORE ACTIVITY**



TRIM HAIR



CLEAN SKIN WITH RUBBING ALCOHOL

YOU WILL NEED

1FULL STRIP

2 HALF STRIPS







STRIPS ONE:

Tear backing paper two inches from end to create starting anchors.

STRIP TWO &

THREE: Tear backing paper at the center of strips to begin application.



TO CREATE A HALF STRIP:

Fold a full strip in half and cut folded edge rounding the corners.

ANCHOR	STRETCH	ANCHOR
(L :		+)
		!
\leftarrow		
- :		



WATCH STEP-BY-STEP **VIDEO INSTRUCTIONS.**

HOW TO APPLY

STRIP ONE



Flex the foot to create a stretch on the back of the ankle. Place the anchor on the bottom of the foot with no stretch.



Apply tape with moderate stretch under the heel and along the achilles tendon.



Lay down the second anchor with no stretch using the last two inches of tape.

STRIP TWO



Apply moderate stretch and adhere a half strip across the first strip below the ankle.



Lay down the last two inches at each end with no stretch.

STRIP THREE



Repeat the same technique with the second half strip above the ankle. Smooth over the full application to ensure tape is completely applied to the skin.