

# BIOMIC® V3

## Clinical Microbiology System

### Open system to automate and standardize microbiology tests

- Instant digital reading & interpretation
- 2-way LIS/LIMS interface
- CLSI and EUCAST guidelines
- Save high resolution test images
- Resistance monitoring & antibiograms
- 24 inch touch screen monitor
- No routine maintenance



### Optional Testing Modules:




#### Disk Diffusion

- BD, Oxoid, etc
- Etest® & MIC Test Strips




#### Organism Identification

- API® (bioMerieux)
- RapID™ (Thermo)
- Crystal™ (BD)
- Liofilchem® ID Systems



#### Agar Dilution

- Standard Methods
- Custom Methods



#### Broth Microdilution\*

- Custom Plates
- Commercial Plates




#### Colony Count

- Colony Count Plates
- Chromogenic Agar Plates



#### MIC Strip

- Etest® (bioMerieux)
- MTS® (Liofilchem)



#### Urine Screen

- Custom Methods
- Commercial Methods

\*Broth Microdilution - Note in USA: For Research Use Only. Not for use in diagnostic procedures.