

PRODUCT CATALOG

DISPOZ

CLINICAL PLASTICWARES



FC-BIOS Sdn Bhd
(188200-H)
Excellence in Laboratory & Process Supplies

FC-BIOS SDN BHD (HQ)
31 & 33, Jalan SS15/4B,
47500 Subang Jaya, Selangor
Tel: +603-5635 1559
Fax: +603-5635 1306

FC-BIOS SDN BHD (PENANG)
462A, Jalan Permatang Rawa, Bandar Perda,
14000 Bukit Mertajam, Pulau Pinang, Malaysia
Tel: +604-5380 559
Fax: +604-5391 306

CONTENTS

Centrifuge Tube	3
Microtube	4
Test Tube	5
Transfer Pipette	6
Pipette Tips	7
Swab	8
Inoculating Loops	9
	10
Petri Dish	10
L-Type Cell Spreader	11
Embedding Cassette	12
Flange Cup Cap	13
Sample Cup	14
Sample Container	15

- Manufactured from high quality polypropylene
- Sterilized with Ethylene Oxide (EO)
- Temperature range of -80°C to 121°C
- Withstand centrifugal forces to 9500 xg RCF (15ml) and 10,500 xg RCF (50ml)
- Can be autoclaved and will resist temperatures up to 121°C
- Graduations on clear tubes are printed in bold easy-to-read black ink
- Packed in individual bags of 50pcs (15ml) & 25pcs (50ml) for easy distribution
- White writing area provided
- Medical Device Authority (MDA) approved
Registered under Section 5(1) of Act 737
MDA Reg. No.: IVD81814793218A



15ml Centrifuge Tube



50ml Centrifuge Tube

Cat No.	Description	Size	Pack	Case
LB1S00155010	Centrifuge Tubes- Bulk Packing	15ml	50	1000
LC1S00502505	Centrifuge Tubes- Bulk Packing	50ml	25	500

- Manufactured from high quality polypropylene
- Volume of 1.5ml
- Eppendorf type
- Graduated at 0.1mL, 1.0mL, 1.5mL
- Temperature range of -80°C to 121°C
- Can be autoclaved and will resist temperatures up to 121°C
- Withstand centrifugal forces to 17,500 xg RCF
- Medical Device Authority (MDA) approved

**Registered under Section 5(1) of Act 737
MDA Reg. No.: IVD57173796118A**



Cat No.	Description	Size	Pack	Case
LB1N00150510	Microtube Eppendorf Type	1.5ml	500	10,000



Test Tube - PP

- Manufactured from high quality polypropylene
- Available with and without stopper
- Temperature range of -80°C to 121°C
- Can be autoclaved and will resist temperatures up to 121°C
- Leak-proof
- Stoppers easily removed with single hand
- Medical Device Authority (MDA) approved



Test Tube - PS


- Manufactured from high quality polystyrene
- Available with and without stopper
- High clarity
- Temperature range of -30°C to 60°C
- Non-autoclavable
- Leak-proof
- Stoppers easily removed with single hand
- Medical Device Authority (MDA) approved

Registered under Section 5(1) of Act 737
MDA Reg. No.: IVD88704932218A


Cat No.	Description	Size	Diameter	Bag	Case
TA1N00050505	Test Tube – PS	5ml	12 x 75mm	500	5000
TA2N00050505	Test Tube – PP	5ml	12 x 75mm	500	5000

- Manufactured from high quality low density polypropylene
- Length – 160mm
- Graduated at 1mL, 2mL, 3mL
- Ideal for sampling and decanting infectious or toxic liquids
- Safe, convenient and non-toxic
- Flexible stem can be bent to draw liquid from narrow mouth test tubes
- Non-sterile or Sterile (EO) available
- High clarity plastic with low liquid retention properties
- Medical Device Authority (MDA) approved

Registered under Section 5(1) of Act 737
MDA Reg. No.: IVD35120796018A
MDA Reg. No.: IVD22198781218A



1ml Transfer Pipette



3ml Transfer Pipette



3ml Sterile Disposable Transfer Pipette

Cat No.	Description	Size	Pack	Case
LA2S0033502	Sterile Disposable Transfer Pipette	3ml	500	2000
LA1N00310505	Disposable Transfer Pipette	1ml	500	5000
LA2N00330505	Disposable Transfer Pipette	3ml	500	5000

- Manufactured from high quality polypropylene
- Can be autoclaved and will resist temperatures up to 121°C

- Temperature range of -80°C to 121°C
- Medical Device Authority (MDA) approved

Registered under Section 5(1) of Act 737
MDA Reg. No.: IVD11227802618A
MDA Reg. No.: IVD38448802518A



Pipette Tips - Gilson Type



Pipette Tips – Eppendorf Type

Cat No.	Description	Size	Type	Pack	Case
LE1N00011030	Pipette Tips - Yellow	200µL	Eppendorf	1000	30,000
LE3N00100510	Pipette Tips - Blue	1000µL	Eppendorf	500	10,000
LE2N00201030	Pipette Tips - Yellow	200µL	Gilson	1000	30,000
LE4N00100510	Pipette Tips - Blue	1000µL	Gilson	500	10,000



**Swab Applicator in Pouch
(individually wrapped)**

- Sterilized with Ethylene Oxide (EO)
- Length - 150mm
- Cotton bud diameter - 5mm
- Medical Device Authority (MDA) approved

**Registered under Section 5(1) of Act 737
MDA Reg. No.: IVD16846796918A**



Rayon Transport Swab


- Sterilized with Ethylene Oxide (EO)
- Label provided
- Length - 150mm
- Cotton bud diameter - 5mm
- Plain swab (no media)
- Medical Device Authority (MDA) approved

**Registered under Section 5(1) of Act 737
MDA Reg. No.: GMD76695796418A**

Cat No.	Description	Diameter	Pack	Case
MA4S00150505	Swab Applicator in Pouch	5mm	500	5000
MB1S20160102	Rayon Transport Swab	5mm	100	2000

- Manufactured from high quality polypropylene
- Volume of 1ul or 10ul
- Disposable loops with needle type for collection and inoculation by streaking or puncturing method
- Sterilized by Ethylene Oxide (EO)
- Medical Device Authority (MDA) approved

Registered under Section 5(1) of Act 737
MDA Reg. No.: IVDA10906122-94464



Cat No.	Description	Size	Pack	Case
ML2S00010106	Inoculating loops with needle tip	1ul	10	6000
ML1S00100106	Inoculating loops with needle tip	10ul	10	6000



Disposable Petri Dish

- Manufactured from non-cytotoxic virgin polystyrene
- Manufactured in an ISO clean room
- Excellent optical clarity
- Stacking ring for easy storage and handling
- Sterilized with EO gas

Cat No.	Description	Diameter	Pack	Case
DZ02-MP1S00900105	Disposable Petri Dish	90 x 15 mm	20	500




Cell Spreader L-Type

- Manufactured from non-cytotoxic virgin polystyrene
- Manufactured in an ISO clean room
- Sterilized with EO gas

Cat No.	Description	Diameter	Pack	Case
DJ01-MC1S00005105	Cell Spreader L-Type	12mm	500	5000

- Manufactured from high quality polypropylene
- Available in yellow & white
- Medical Device Authority (MDA) approved

Registered under Section 5(1) of Act 737
MDA Reg. No.: IVD15011796818A



Cat No.	Description	Color	Bag	Case
HA1N00002525	Embedding Cassette	White	250	2500
HA2N00002525	Embedding Cassette	Yellow	250	2500

- Manufactured from high quality polypropylene
- Stoppers easily removed with single hand
- Available in red & blue
- Medical Device Authority (MDA) approved

Registered under Section 5(1) of Act 737
MDA Reg. No.: IVDA3844319-30093



Cat No.	Description	Diameter	Pack	Case
FA1N00120520	Flange Plug Cap	12mm	500	20000

- Manufactured from high quality polystyrene
- Volume of 3mL
- Compatible with Hitachi 7150 series & Hitachi 7060 series
- Medical Device Authority (MDA) approved
Registered under Section 5(1) of Act 737
MDA Reg. No.: IVDA3441019-30469



Cat No.	Description	Size	Pack	Case
DA1N00030550	Sample Cups	3ml	500	5000



Specimen Container

- Manufactured from high quality polypropylene
- Sterilized with Ethylene Oxide (EO)
- Leak proof
- Label provided
- Dimensions - 40mm (D) x 65mm (H)
- Volume of 60mL
- Smooth opening and closing of cap
- Medical Device Authority (MDA) approved

Registered under Section 5(1) of Act 737
MDA Reg. No.: IVDA6824922-89519



Stool Container with Spoon

- Manufactured from high quality polystyrene
- Sterilized with Ethylene Oxide (EO)
- Dimensions - 25mm (D) x 80mm (H)
- Label provided
- Volume of 30mL
- Medical Device Authority (MDA) approved

Registered under Section 5(1) of Act 737
MDA Reg. No.: IVD49522793318A

Cat No.	Description	Size	Pack	Case
CA2S10600205	Specimen Container (with label)	60ml	20	500
CA4S00600205	Specimen Container (without label)	60ml	20	500
CB1S10305005	Stool Container with Spoon	30ml	50	500

SELANGOR - HQ

31 & 33, Jalan SS15/4B, 47500 Subang Jaya,
Selangor Darul Ehsan, Malaysia.

PENANG - BRANCH

462A, Jalan Permatang Rawa, Bandar Perda, 14000 Bukit Mertajam,
Pulau Pinang, Malaysia.



+6019-220 8238



+603-5635 1559/+604-538 0559



+603-5635 1306/+604-539 1306



enquiry@fcbios.com.my



www.fcbios.com.my