

colitag™

The revolutionary
16-48 hour water test



EPA
APPROVED



GET YOUR FREE
SAMPLE TODAY!



Colitag™ 16-48 Hour Testing Window

Colitag™ is a presence-absence (P/A) and MPN enzyme substrate test that features a 16-48 hour read window for total coliform and E. coli determinations and detects 1 MPN of total coliform or E. coli in a 100mL water sample. Colitag™ is relied upon in a growing market of over 24 countries.



Run It When You're Ready, Read It When You're Ready

Colitag™ offers unrivalled flexibility, allowing you to use one test instead of three to read your samples at any time within 16-48 hours. Increase lab efficiency and get the job done on your own time.



Improve Your Laboratory Efficiency

Colitag's bright blue fluorescence indicating a positive E. coli sample eliminates the need for time-consuming confirmations sometimes required by other enzyme substrate media with lower intensity fluorescence. This teamed with the flexible reading window allows you to streamline your lab processes.

Comparison

The Colitag™ offers up to 8x more flexibility compared to the industry standard enabling better workflow planning & management, simultaneously improving efficiency of the laboratory.



Colitag




Brand C



Brand C-18

Minimum Incubation	Shortest	Longest	Moderate
Incubation Period	16 to 48 hours	24 to 28 hours	18 to 22 hours
Shelf-life	18 months	12 months	15 months
Price	\$	\$\$	\$\$\$

Valid Results for 32 hours*



Colitag™ package

Package A

Presence-absent method

» Colitag water test kit


- N13-9850 Colitag water test kits, 20 tests/kit
- N13-9851 Colitag water test kits, 100 tests/kit

» Colitag Comparator

- N13-9857 Colitag comparator, 100 mL

» Sampling Bottles

- N13-9859 Sample containers, screwtop, with sodium thiosulfate, 100/pack
- N13-9856 Sample containers, screwtop, without sodium thiosulfate, 100/pack



Package B

Presence-absent method

» Colitag water test kit

- N13-9850 Colitag water test kits, 20 tests/kit
- N13-9851 Colitag water test kits, 100 tests/kit

» Colitag Comparator

- N13-9857 Colitag comparator, 100 mL

» Sampling Bag

- N02-B01402 Whirl-Pak® Stand Up Thio-Bags®, 100 mL, 100/pack

