

COVID-19 HiViral™ Transport Swab Kit

Specially used for collection, transportation and preservation of human nasopharyngeal or oropharyngeal virus samples from the collection site to the testing laboratory.

What's included in the kit?

- 3.0ml viral transport medium in a self-standing conical bottom tube.
- Sterile nasopharyngeal nylon flocked swab with breakpoint
- Sterile oropharyngeal viscose swab with breakpoint.
- Biohazard bag.



ISO 13485 and ISO 9001 certified, FDA & MDA Approved.



Patient friendly swab with ergonomic design and soft texture for maximum yield.



Support a wide range of viruses including SARS-CoV-2 (2019-nCoV)

Item Code	Item Description	Packing
H05-MS2760S-50NO	HiViral Transport Kit (non-inactivated type) - in bulk packing	50 tests/pack
H05-MS5760L-50NO	HiViral Transport Kit (inactivated type) - in bulk packing	50 tests/pack



MDA Reg. No.: IVDA7300421-73898



Whatsapp us to request quotation & get a FREE sample
019-220 8238 (Mr. Farris)