

BULK BAG PIPETTE TIPS (NON-STERILE)

Category	Cat. No.	Volume & Type	Unit/Pk	Pk/Cs
Non-filter tips	T521000	10ul graduated	1000	20
	T521010-Y	200ul bevelled yellow	1000	20
	T521014-Y	200ul graduated yellow	1000	20
	T521015	300ul graduated	1000	20
	T521016	1000ul graduated	500	10
	T521016-B	1000ul graduated blue	500	10
	T521017	1000ul wide bore	500	10
	T521031	5ml graduated	250	10
Filter Tips	T521033	10ml graduated	100	10
	T527100	10ul graduated	1000	20
	T527104	200ul graduated	1000	10
MAXIPENSE™ Ultra Low Retention Tips (non-filter tips)	T527106	1000ul graduated	500	10
	T523100	10ul graduated	1000	20
	T523101	200ul bevelled	1000	10
	T523102	200ul graduated	1000	10
	T523103	300ul graduated	1000	10
	T523105	1000ul graduated	500	10



Clean room manufactured | Certified DNase, RNase & pyrogen free
100% pure virgin polypropylene | Plasma treated low retention tips
ISO 13485 & ISO 9001 certified

RACKED PIPETTE TIPS (STERILE)

Category	Cat. No.	Volume & Type	Unit/Pk	Pk/Cs
Non-filter tips	T522100	10ul graduated	960 (10 racks)	5
	T522101-Y	200ul graduated yellow	960 (10 racks)	5
	T522102	300ul graduated	960 (10 racks)	5
	T522103-B	1000ul graduated blue	960 (10 racks)	5
MAXIPENSE™ Ultra Low Retention Tips (filter tips)	T529100	10ul graduated	960 (10 racks)	5
	T529101	20ul graduated	960 (10 racks)	5
	T529103	100ul graduated	960 (10 racks)	5
	T529104	200ul graduated	960 (10 racks)	5
	T529106	1000ul graduated	960 (10 racks)	5

Compatible with leading pipette brands





Pipette Tip Compatibility Chart



Brand	Model	Volume Range	Channels	10µl		200µl		300µl		1000µl	
				NF	F	NF	F	NF	F	NF	F
Biohit	mLine	0.1-3µl, 0.5-10µl	1	✓	✓						
		2-20µl, 10-100µl, 20-200µl	1			✓	✓				
		100-1000µl	1							✓	✓
		0.5-10µl	8/12	✓	✓						
		5-100µl	8/12			✓	✓				
		30-300µl	8/12					✓	✓		
	Proline	0.1-2.5µl, 0.5-10µl	1	✓	✓						
		2-20µl, 5-50µl, 10-100µl, 20-200µl	1			✓	✓				
		100-1000µl	1							✓	✓
		0.5-10µl	8/12	✓	✓						
		5-50µl	8/12			✓	✓				
		50-300µl	8/12					✓	✓		
Brand	Transferpette	0.1-1µl, 0.5-10µl	1	✓	✓						
		2-20µl micro	1								
		2-20µl, 5-50µl, 10-100µl, 20-200µl	1			✓	✓				
		100-1000µl	1							✓	✓
		0.5-10µl	8/12	✓	✓						
		2-20µl	8/12			✓	✓				
		2.5-25µl, 5-50µl, 10-100µl, 20-200µl	8/12			✓	✓				
		30-300µl	8/12					✓	✓		
	Transferpette S	0.1-1µl, 0.1-2.5µl, 0.5-10µl	1	✓	✓						
		2-20µl, 10-100µl, 20-200µl	1			✓	✓				
		100-1000µl	1							✓	✓
		0.5-10µl	8/12	✓	✓						
		5-50µl, 10-100µl, 20-200µl	8/12			✓	✓				
		30-300µl	8/12					✓	✓		
Eppendorf	Research Plus	0.1-2.5µl, 0.5-10µl	1	✓	✓						
		2-20µl micro	1								
		2-20µl, 10-100µl, 20-200µl	1			✓	✓				
		30-300µl	1					✓	✓		
		100-1000µl	1							✓	✓
		0.5-10µl	8/12	✓	✓						
		10-100µl	8/12			✓	✓				
		30-300µl	8/12					✓	✓		
	Reference	0.1-2.5µl, 0.5-10µl	1	✓	✓						
		2-20µl micro	1								
		2-20µl, 10-100µl, 50-200µl	1			✓	✓				
		100-1000µl	1							✓	✓

NF: non-filter tip F: filter tip

Brand	Model	Volume Range	Channels	10µl		200µl		300µl		1000µl	
				NF	F	NF	F	NF	F	NF	F
Gilson	Pipetman Classic, G, L: P2	0.2-2µl, 1-10µl	1	✓	✓						
	Pipetman Classic, G, L: P20	2-20µl, 20-100µl, 50-200µl	1			✓	✓				
	Pipetman Classic, G, L: P1000	200-1000µl	1							✓	✓
	Pipetman P10L	0.5-10µl, 2-20µl, 20-200µl	8/12			✓	✓				
	Pipetman P300L	20-300µl	8/12					✓	✓		
Nichiryo	Nichipet Premium, EX, EX Plus	0.1-2µl	1	✓	✓						
		2-20µl, 10-100µl, 20-200µl	1			✓	✓				
		100-1000µl	1							✓	✓
	Nichipet 7000	0.5-10µl	8/12	✓	✓						
		5-50µl, 40-200µl	8/12			✓	✓				
		50-300µl	8/12					✓	✓		
Rainin	Pipet-Lite SL-2XLS	0.1-2µl, 0.5-10µl	1	✓	✓						
	Pipet-Lite SL-20XLS	2-20µl, 10-100µl, 20-200µl	1			✓	✓				
	Pipet-Lite SL-300XLS	20-300µl	1					✓	✓		
	Pipet-Lite SL-1000XLS	100-1000µl	1							✓	✓
	Classic PR-2	0.1-2µl, 0.5-10µl	1	✓	✓						
	Classic PR-20	2-20µl, 10-100µl, 20-200µl	1			✓	✓				
	Classic PR-1000	100-1000µl	1							✓	✓
Socorex	Calibra® 822	0.2-2µl, 1-10µl micro	1	✓	✓						
		1-10µl, 2-20µl, 10-100µl, 20-200µl	1			✓	✓				
		100-1000µl	1							✓	✓
		1-10µl	8	✓	✓						
		10-100µl, 20-200µl	8/12			✓	✓				
	Acura 825	0.1-2µl, 0.5-10µl	1	✓	✓						
		1-10µl, 2-20µl, 5-50µl	1			✓	✓				
		100-1000µl	1							✓	✓
		0.5-10µl	8/12	✓	✓						
		5-50µl, 20-200µl	8/12			✓	✓				
		40-350µl	8/12					✓			
Thermo/Finn	Finnpette 4500 /F1 /F2	0.2-2µl micro, 0.5-5µl micro	1	✓		✓		✓			
		0.5-5µl, 2-20µl, 5-50µl, 10-100µl, 20-200µl	1			✓	✓	✓	✓		
		100-1000µl	1							✓	✓
		0.5-10µl	8/12	✓		✓		✓			
		5-50µl	8/12			✓	✓	✓	✓		
		50-300µl	8/12			✓	✓	✓	✓		

NF: non-filter tip F: filter tip