

Simple is
always the
best.

TOMY

SX-500 | SX-700 | FLS-1000



Watch
product video



YouTube

<https://bit.ly/3yH37RZ>

Test it today.




WhatsApp us for
FREE Trial & Demo
Mr. Finn
019-220 8238



Autoclave Specifications

MODEL/SPECIFICATIONS		T01-SX-500	T01-SX-700	T01-FLS-1000
Operating Temperature Range	Sterilizing	105 - 135°C (0.019 - 0.212MPa)		105 - 123°C
	Heating	45 - 104°C (0 - 0.015MPa)		
	Warming	45 - 95°C		
Temperature Range	Display	Digital		
Pressure Gauge	Display	Analog		
Chamber Dimension		φ 325 x 733mm	φ 370 x 774mm	φ 450 x 678mm
Chamber Capacity		Internal volume: 58L	Internal volume: 79L	Internal volume: 100L
Body Dimension		410W, 477D, 970H	470W, 528D, 1003H	680W, 600D, 1022H
Weight		60kg	72kg	108kg

Autoclave Height Design




Autoclave Accessories



Autoclave Long Basket

Model	Item Code	Dimensions
SX-500	ACA-315T	Φ 300 x 557 mm
SX-700	ACA-700T	Φ 345 x 561 mm
FLS-1000	ACA-1000TP	Φ 415 x 456 mm



Support Stand (Qty: 1)

Model	Item Code
SX-500	SBSS-325
SX-700	SBSS-345



Autoclave Bucket for Waste

Model	Item Code	Dimensions
SX-500	ACA-315K	Φ 300 x 280 mm
SX-700	ACA-700K	Φ 345 x 280 mm
FLS-1000	ACA-1000K	Φ 415 x 239 mm



Sterilizing Container, ST-Zero

Model	Item Code	Dimensions
ST-Zero	151518	106W, 350D, 80H

*Can accommodate 25 x 25ml disposable pipettes.



Standard Autoclave Basket

Model	Item Code	Dimensions
SX-500	ACA-315B	Φ 300 x 182 mm
SX-700	ACA-700B	Φ 345 x 182 mm
FLS-1000	ACA-1000B	Φ 415 x 233 mm



Built-In Printer

Model	Item Code
FLS-1000	T01-FLS1000PRINT



Durham Tube Basket

Model	Item Code	Dimensions
SX-500	STD-500	128, 16 mm dia.
SX-700	STD-700	172, 16 mm dia.



External Sensor

Model	Item Code
FLS-1000	T01-203590013



Chemical Indicator

Model	Item Code
ProChem SSW	R01-CI-SSW
Autoclave Tape	A72-YST-1950



Autoclave Deodorizer

Model	Item Code
HiMedia Fresh	H05-LA008A-10X50NO
Deodorizing Pearls	H05-LA008A-5X50NO



Biological Indicator

Model	Item Code
EzTest	R01-EZS/5 or R01-EZS/6
Magna Amp	R01-MA/6 or R01-MA/6-2



Autoclave Disposable Bag

Model	Item Code
Autoclavable Bags, 12"x16"	S28-AB-1216P
Autoclavable Bags, 18"x30"	S28-AB-1830P



Autoclave Cleaning Kit

Model	Item Code
Cleaning Kit	F07-CLEANING KIT