



# Rack-tube Compatibility Chart

**1 Know your tube external diameter & rack mouth size.**

**2 Check with the compatibility list below.**

Rack mouth diameter should be at least 1-2 mm bigger than test tube external diameter.



	Mouth Size (mm)	Cat No.	Centrifuge Tube/ Tube Size (ml)						Other Tube Mouth Size (mm)					
			0.5	1.5	2.0	5.0	15	50	≤ 12	16	18	20	25	30
<b>POLYWIRE™ Test Tube / Half Test Tube Rack</b> 	13	T202110-W/ T201010-W							✓					
	16	T202111-W/ T201011-W							✓					
	20	T202112-W/ T201012-W				✓	✓		✓	✓	✓			
	25	T202113-W/ T201013-W				✓	✓			✓	✓	✓		
	30	T202114-W/ T201014-W				✓	✓	✓		✓	✓	✓	✓	
<b>Test Tube Stand</b> 	20	T202030-W					✓			✓	✓			
	25	T202040-W					✓			✓	✓	✓		
	30	T202050-W					✓	✓			✓	✓	✓	
<b>COMBILOK™ 4 Way Tube Rack</b> 	Multi	T202081-A	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓	✓
<b>Universal Combi 4 Way Tube Rack</b> 	Multi	T202080-A				✓	✓	✓	✓	✓	✓	✓	✓	✓