


# Chamber capacity



Through its unique top-loading design, compact size and large capacity, TOMY autoclave stackable basket system allows up to 3 baskets (2 baskets for FLS-1000) in the chamber to maximize loads and reduced number of cycle needed to sterilize your items.



## Chamber capacity comparison

Schott bottles volume	SX-500 (58ℓ)	Brand-H (50ℓ)	SX-700 (79ℓ)	Brand-H (85ℓ)	FLS-1000 (100ℓ)	Brand-H (110ℓ)
10 mL	153	126	204	168	202	168
25 mL	153	126	204	168	202	168
50 mL	93	72	126	102	124	102
100 mL	63	48	84	68	82	68
150 mL	48	39	66	54	64	54
250 mL	39	27	51	42	50	42
500 mL	24	21	33	28	34	28
750 mL	14*	15	18*	22	26	22
1000 mL	10*	8	16*	18	24	18
2000 mL	6*	2	8*	10	6	10

\*with autoclave support stand