



We provide all PPE Products

Isolation Gowns

Essential protection for the frontline



New York

Toronto

Montreal Head Office

4105-B Matte Blvd. Brossard,
Quebec J4Y 2P4 Canada

Contact Us for a Quote

 +1-289-339-4440

 info@med-supply.net

 www.med-supply.net



About Us.....	p. 3
Medical and Protective Gowns.....	p. 7
AAMI Level Gowns.....	p. 12
AAMI Level 1 Gowns.....	p. 15
AAMI Level 2 Gowns.....	p. 18
AAMI Level 3 Gowns.....	p. 25
Protective Gowns.....	p. 39



About Us

Med Supply specializes in the importation and distribution of products from international markets in a B2B setting. Our Canadian division has extensive experience in sourcing, from four different continents, and with teams on the ground in each jurisdiction. Med Supply offers a door-to-door full service experience with an emphasis on cost efficiency.

Within our Medical Division, Med Supply is committed to providing excellent service to, and high-quality products for, frontline healthcare professionals, hospitals, clinics, patients, the food service industry, among others.



Health
Canada

Santé
Canada

id: 155050

Licence Number

12569

Numéro de la licence

**Medical Device
Establishment Licence**

**Licence d'établissement
pour les instruments médicaux**

A2Z WHOLESALE INC.

4105 BOULEVARD MATTE
SUITE B
BROSSARD, QUEBEC
CANADA
J4Y 2P4

FDA | U.S. Food and Drug Administration



Registration Confirmation

Facility: A2Z WHOLESALE INC. Brossard, Quebec, CANADA

You have successfully entered your facility registration and device listing information. You should print a copy of this page for your records. Listing numbers appear below for the products manufactured, developed, or processed at this facility.

As a manufacturer, specification developer, or single-use device reprocessor, you are required to pay an annual fee for medical device facility registration.

You will receive another e-mail providing you with your registration number in approximately 30 to 90 days. Until your registration number is assigned, reference your Owner/Operator number in any correspondence with the Center for Devices and Radiological Health.

Your registration will be valid through Dec 31, 2020. An e-mail will be sent to the Owner/Operator and the Official Correspondent 90 days before the facility is required to re-register for Fiscal Year 2020 with instructions on how and when to re-register.

Note: Registering your device facility and listing your devices does not, in any way, constitute FDA approval of your facility or your devices.

Should you have any questions, please send an e-mail to reglist@cdrh.fda.gov (mailto:reglist@cdrh.fda.gov).

The Owner/Operator Number for this Registration is: 10074953

Facility Information

Initial Importer: N

Facility Name: A2Z WHOLESALE INC

Address: 4105-B Matte Blvd
Brossard, Quebec, J4Y 2P4, CANADA

Foreign Trade Zone: N

6/3/2020

FDA Registration for A2Z WHOLESALE INC

Owner/Operator Information

Owner/Operator Number: 10074953

Contact Name: Ishan Singh

Company: A2Z WHOLESALE INC

Address: 4105-B Matte Blvd
Brossard, QUEBEC, J4Y 2P4, CANADA

Telephone: 1 - 514 - 3836660

Fax: -

E-mail: ricky@shumaker.ca

Official Correspondent Information

Contact Name: David Lennarz

Company: Registrar Corp

Address: 144 Research Drive
Hampton, VIRGINIA, 23666, UNITED STATES

Telephone: 1 - 757 - 2240177

Fax: 1 - 757 - 2240179

E-mail: david.lennarz@registrarcorp.com

DUNS Number: 139242874

United States Agent Information

Contact Name: David Lennarz

Contact Title: Mr

Business Name: Registrar Corp

Address: 144 Research Drive
Hampton, Virginia, 23666, UNITED STATES

Phone: 757 - 2240177

Fax: 757 - 2240179

DUNS Number: 139242874

E-mail: david.lennarz@registrarcorp.com

6/3/2020

FDA Registration for A2Z WHOLESALE INC

Device Listings

Listing Number	Premarket Submission Number/Type	Product Code(s)	Device Name(s)	Activities	Importers
D406670	Exempt	OEA	Non-surgical isolation gown	Foreign Exporter	

Date of Initial Registration: Wed Jun 03 08:28:02 EDT 2020

- Federal Corporation Number: 924712-2
- Business Number: 803045780
- Ontario Corporation Number: 003112591
- NEQ: 1174250168
- Health Canada NPN: 80099806
- CA Number: CA64578
- GST Number: 803045780
- QST Number: 1226258414



MEDICAL AND PROTECTIVE GOWNS



About Us

Med Supply range of medical and protective gowns protect the clothing of health care providers. They also protect visitors and patients because they prevent the transfer of microorganisms and body fluids in patient isolation situations.

Pathogens are unable to penetrate the material, which protects both the wearer and the patient.

Gowns are one part of an overall infection-control strategy, and can therefore benefit many other industries, such as food and service, aviation, tourist, hotel and hospitality, etc.

1. Polyethylene “PE”

Isolation gowns in PE material are lightweight, cost-effective and highly protective. The single layer allows breathability and flexibility of movement.

2. Polypropylene “PP” Spunbond Non-woven

Polypropylene material is spun and then directly dispersed into a web to create a spunbond non-woven a fabric-like material. Single-layer, lightweight gowns offer high degree of breathability and comfort with adequate protection.

3. PP + PE Laminated Non-woven

Combining the best attributes of polyethylene and polypropylene ensures a strong protection while maximizing comfort. Isolation gowns in this material are multi-layer and medium weighted.

4. SMS Lightweight

Multi-layer Spunbond-Meltblown-Spunbond is a 3-ply lightweight material that offers moderate protection while ensuring comfort and breathability.

5. SMS Medium weight

Multi-layer Spunbond-Meltblown-Spunbond is a 3-ply medium weight material that offers strong protection while ensuring comfort and breathability.

6. PPM – PP Moki-Tex

PP Moki-Tex is a proprietary, innovative and industry-first material that is exclusive to A2Z.

PPM offers a high level of comfort and breathability, a moderate protection, and still remains very affordable.

A. Open Back



B. Half Back Closed



Neck Style- Over the head



C. Full Cover



D. Apron Style



Neck Style- Neck tie



Sleeve and Cuff Options



Thumbloop Sleeve




Elastic Cuff Sleeve




Nylon Knit Cuff Sleeve


Gowns' Colour Options



Blue




Yellow



White

AAMI LEVEL GOWNS

ANSI/AAMI PB70 STANDARDS BARRIER PROTECTION



AATCC-42: water resistance impact penetration test determines the ability of a material to resist water penetration under spray impact.

AATCC-127: water resistance hydrostatic pressure test determines the ability of a material to resist water penetration under constant contact with increasing pressure.

In accordance with ANSI/AAMI standards, all of our AAMI Level gowns are tested to meet the highest level of protection for which the gown is rated.

As shown in the figure, all areas of the medical gown except bindings, cuffs, and hems are considered critical zones of protection and are tested as per the AAMI PB70 standard. All seams must have the same liquid barrier protection as the rest of the gown.

The critical zone comprises at least areas A and B. The back of the gown, shown as area D, may be non protective.

We offer a complete range of AAMI Level isolation gowns for both low risk and high risk environments. We work with internationally-accredited laboratories to conduct the necessary tests to ensure our gowns meet the highest level of protection for which each gown is rated.

Level 1

Model 6C

- AATCC 42 - Water penetration ≤ 4.5 g
- Used for minimal risk situations, such as basic care, standard hospital medical unit and cover gowns for visitors

Level 2

Model 1BY; Model 4C

- AATCC 42 - Water penetration ≤ 1.0 g
- AATCC 127 - Hydrostatic pressure ≥ 20 cm water column
- Used for low risk situations, such as blood drawing procedures, suturing, ICU and use in pathology labs
- Provides a barrier to larger amounts of fluid penetration through splatter and some fluid exposure through soaking

Level 3

Model 1A; Model 1B; Model 1CY;
Model 3C; Model 5C

- AATCC 42 - Water penetration ≤ 1.0 g
- AATCC 127 - Hydrostatic pressure ≥ 50 cm water column
- Used in moderate risk solutions, such as inserting an IV line, arterial blood draw, emergency room, trauma, etc.
- Provides a barrier to larger amounts of fluid penetration through splatter and more fluid exposure through soaking than Level 2

Level 4

- F1671 (Viral): No penetration at 2 psi
- Used in high-risk situations, such as during long, fluid intense surgical procedures

AAMI Level 1 Gowns

Technical Specs

AAMI Level 1



- Classic gown with front and back cover
- Fabric-like non-woven first-to-market Moki-tex material
- Neck tie
- Back waist-tie
- Proprietary Moki-Tex technology
- Lightweight
- AATCC 42 - Water penetration ≤ 4.5 g
- Sleeve options: Elastic cuff or Nylon cuff
- Non sterile
- Colour options: Blue or Yellow
- Size: Universal

SKU	Size/Description	Qty.
Classic gown, flatpack, choice of sleeve and colour		
6C-IREYU1	Elastic cuff, yellow	5/pk, 100/cs
6C-IRNYU1	Nylon cuff, yellow	5/pk, 100/cs
6C-IREBU1	Elastic cuff, blue	5/pk, 100/cs
6C-IRNBU1	Nylon cuff, blue	5/pk, 100/cs



Model 6C – Universal Size

medsupply



AAMI Level 2 Gowns

Technical Specs

AAMI Level 2


- Impervious
- Easy donning and quick removal
- Partial open-back for breathability
- Back waist-tie
- PE material
- Lightweight
- AATCC 42 - Water penetration ≤ 1.0 g
- AATCC 127 - Hydrostatic pressure ≥ 20 cm water column
- Thumbloop sleeve
- Non sterile
- Size: Universal & X-Large

SKU	Size/Description	Qty.
Over-the-head, thumbloop, flatpack, blue		
1BY-OTBU2	Universal	5/pk, 200/cs
1BY-OTBX2	X-Large	5/pk, 200/cs




Model 1BY – Universal Size

medsupply



Model 1BY – X-Large Size


medsupply



Technical Specs

AAMI Level 2



- Classic gown with front and back cover
- Multi-layer protection with maximum breathability
- Neck tie
- Back waist-tie
- SMS material
- Lightweight
- AATCC 42 - Water penetration ≤ 1.0 g
- AATCC 127 - Hydrostatic pressure ≥ 20 cm water column
- Sleeve options: Elastic cuff or Nylon cuff
- Non sterile
- Colour options: Blue or Yellow
- Size: Universal



SKU	Size/Description	Qty.
SMS gown, flatpack, choice of sleeve and colour		
6C-IREYU1	Elastic cuff, yellow	4/pk, 80/cs
6C-IRNYU	Nylon cuff, yellow	4/pk, 80/cs
6C-IREBU1	Elastic cuff, blue	4/pk, 80/cs
6C-IRNBU1	Nylon cuff, blue	4/pk, 80/cs

Model 4C – Universal Size

medsupply



AAMI Level 3 Gowns

Technical Specs

AAMI Level 3


- Impervious
- Easy donning and quick removal
- Open-back for breathability
- Back waist-tie
- PE material
- Lightweight
- AATCC 42 - Water penetration ≤ 1.0 g
- AATCC 127 - Hydrostatic pressure ≥ 50 cm water column
- Thumbloop sleeve
- Non sterile
- Colour options: Blue or Yellow
- Generously sized for greater coverage and flexibility

SKU	Size/Description	Qty.
Over-the-head, thumbloop, flatpack, blue or yellow		
1A-IOTYU3	Thumbloop, yellow	10/pk. 400/cs
1A-IOTBU3	Thumbloop, blue	10/pk. 400/cs



Model 1A – Universal Size

medsupply



Technical Specs

AAMI Level 3


- Impervious
- Easy donning and quick removal
- Partial open-back for breathability
- Back waist-tie
- PE material
- Lightweight
- AATCC 42 - Water penetration ≤ 1.0 g
- AATCC 127 - Hydrostatic pressure ≥ 50 cm water column
- Thumbloop sleeve
- Non sterile
- Colour options: Blue or Yellow
- Generously sized for greater coverage and flexibility

SKU	Size/Description	Qty.
Over-the-head, thumbloop, flatpack, blue or yellow		
1B-IOTYU3	Thumbloop, yellow	10/pk, 300/cs
1B-IOTBU3	Thumbloop, blue	10/pk, 300/cs



Model 1B – Universal Size

medsupply



Model 1B – Test Reports

[illegible][illegible][illegible][illegible]

SGS

Order No. _____ Date Received _____
MONTMARTRE

Cashier Name _____



(*) Size of Sample (*)


Technical Specs

AAMI Level 3

- Impervious
- Easy donning and quick removal
- Full cover front and back
- Back waist-tie and front waist tie
- PE material
- Lightweight
- AATCC 42 - Water penetration ≤ 1.0 g
- AATCC 127 - Hydrostatic pressure ≥ 50 cm water column
- Thumbloop sleeve
- Non sterile
- Generously sized for greater coverage and flexibility

SKU	Size/Description	Qty.
Over-the-head, thumbloop, flatpack, blue		
1CY-OTBU3	Thumbloop, blue	5/pk, 200/cs






Technical Specs

AAMI Level 3

- Classic gown with front and back cover
- Fabric-like feel on the inside and impervious laminated protection on the outside
- Neck tie
- Back waist-tie
- PP + PE coated material
- Medium weight
- AATCC 42 - Water penetration ≤ 1.0 g
- AATCC 127 - Hydrostatic pressure ≥ 50 cm water column
- Sleeve options: Elastic cuff or Nylon cuff
- Non sterile
- Colour options: Blue or Yellow or White
- Size: Universal





SKU	Size/Description	Qty.
Classic gown, flatpack, choice of sleeve and colour		
3C-IREYU3	Elastic cuff, yellow	4/pk, 80/cs
3C-IRNYU3	Nylon cuff, yellow	4/pk, 80/cs
3C-IREBU3	Elastic cuff, blue	4/pk, 80/cs
3C-IRNBU3	Nylon cuff, blue	4/pk, 80/cs
3C-IREWU3	Elastic cuff, white	4/pk, 80/cs
3C-IRNWU3	Nylon cuff, white	4/pk, 80/cs



Model 3C – Universal Size


medsupply



Technical Specs

AAMI Level 3

- Classic gown with front and back cover
- Multi-layer protection with maximum breathability
- Neck tie
- Back waist-tie
- SMS material
- Medium weight
- AATCC 42 - Water penetration ≤ 1.0 g
- AATCC 127 - Hydrostatic pressure ≥ 50 cm water column
- Sleeve options: Elastic cuff or Nylon cuff
- Non sterile
- Colour options: Blue or Yellow
- Size: Universal





SKU	Size/Description	Qty.
SMS gown, flatpack, choice of sleeve and colour		
5C-IREYU3	Elastic cuff, yellow	5/pk, 50/cs
5C-IRNYU3	Nylon cuff, yellow	5/pk, 50/cs
5C-IREBU3	Elastic cuff, blue	5/pk, 50/cs
5C-IRNBU3	Nylon cuff, blue	5/pk, 50/cs



Model 5C – Universal Size

medsupply



Protective Gowns



Safety starts well before the moment a clinician enters the patient's room. Safety is equally important across various industries, such as the food service industry, the maintenance and cleaning industry as well as other close-contact service industries.

In situations where standard precaution is required, or for contact precaution – whether for visitors, patients or service providers – protective gowns are a key component.

Med Supply range of protective gowns are well suited for patient care, visitors, nursing contact precaution, food service and other close-contact service industries where standard precautions are needed.

Technical Specs



- Classic gown with front and back cover
- Fabric-like non-woven spunbond
- Neck tie
- Back waist-tie
- PP material
- Lightweight
- Sleeve options: Elastic cuff or Nylon cuff
- Non sterile
- Colour options: Blue or Yellow
- Size: Universal

SKU	Size/Description	Qty.
Classic gown, flatpack, choice of sleeve and colour		
2C-IREYU0	Elastic cuff, yellow	5/pk, 100/cs
2C-IRNYU0	Nylon cuff, yellow	5/pk, 100/cs
2C-IREBU0	Elastic cuff, blue	5/pk, 100/cs
2C-IRNBU0	Nylon cuff, blue	5/pk, 100/cs



Model 2C – Universal Size

medsupply



Technical Specs

- Classic apron style coverall protection
- Easy donning and quick removal
- Open-back for breathability
- Back waist-tie
- PE material
- Lightweight
- Sleeveless design
- Non sterile
- Colour options: Blue or yellow
- Universal size

SKU	Size/Description	Qty.
Over-the-head, flatpack, blue or yellow		
1D-Y	Yellow	10/pk, 400/cs
1D-B	Blue	10/pk, 400/cs



Our Happy Clients





We provide all PPE Products

New York

Toronto

Montreal Head Office

4105-B Matte Blvd. Brossard,
Quebec J4Y 2P4 Canada

Contact Us for a Quote

 +1-289-339-4440

 info@med-supply.net

 www.med-supply.net