

TIGER MULTI

Polishin Compound



YAMAHACHI JAPAN

*Packing

120g



400g



*Usage

Polishing for Titanium Alloy, Pure Titanium, Chrome-Cobalt Alloy, Hard and Soft Metal, Acrylics, Composite

*Direction for Use

1. Smudge Tiger Multi to a polishing brush.
2. Polish necessary parts.

*Examples of Use

For Metal Crown Polishing



After Polishing



For Composite Resin Polishing



After Polishing



Unique Dental Supply Inc.
Quality • Service • Innovation

1-888-532-0554
orders@udscanada.com
www.udscanada.com