



Thermo Scientific P10DS Gradient Gel Caster
P10DS Gradient Gel Caster Instructions


1. Load the glass cassettes with the simple clamping mechanism




2. Pour an agarose plug in the specially designed agarose well



3. Cast the acrylamide solution in the cassettes and once solidified, load your samples.



4. Run the gel (s)



© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

>> For more information, visit: www.thermoscientific.com/owlsci