Dynamic Tilt Treatment Table



Dynamic Tilt Table Sliding Top Tilt Table

The Dynamic Tilt Table can help improve patient recovery rates whilst reducing risks to carers. This is achieved through progressive partial weight bearing exercises using the sliding top on the tilt table which can significantly strengthen de-conditioned lower extremities. Repeated exercise on the Dynamic Tilt Table allows physical progress to be assessed in an objective manner.

Patient handling causes many injuries to health care workers and the use of a Dynamic Tilt Table overcomes the need to move patients out of bed on to the floor for weight bearing exercise critical for their recovery.

The Dynamic Tilt Table is built to handle patients up to 350kg. Moving the table to a patient's bedside reduces patient discomfort and the costs of transfers. The optional electronic braking system makes moving the table effortless. Simply transfer the patient to the table, raise the backrest slightly for comfort and resistance, remove the insert, release the lock and your patient is ready to exercise with the level of resistance controlled by the angle of the table. Upholstery Colour: BLUE And Control

Code:Description:Colour:ABTILTDynamic Tilt TableBlue

WHITELEY HEALTHCARE

Phone: 1300 303 755 Fax: 1300 303 788 Email: sales@whiteleyhealthcare.com.au Mail: PO Box 576, Hornsby NSW 1630 Australia

