



您的自測結果必須通過「醫療診斷實驗室」測試來進行進一步驗證。
Your self-test result must be further verified by having a confirmatory laboratory test in medical lab.

注意事項 Warnings and Precautions

輪狀病毒、腺病毒及諾如病毒的快速測試均使用獨立測試條進行，三者之間的測試結果皆不相關。
Rapid test of Rotavirus, Adenovirus and Norovirus are conducted on individual test paper and there is no correlation between the test results.

此測試的結果只可用作初步篩選或在醫生指引下監察病情，用戶應先閱讀說明書來了解測試結果背後的含意。
The results of this test are for preliminary screening only or for monitoring under a doctor's guidance. Users should read the instruction manual to understand the implications of the test results.

此為定性測試，橫線的顏色深淺不影響判斷測試結果屬陽性或陰性。
This is a qualitative test, and the shade of the line does not affect whether the test result is positive or negative.

儲存條件: 2-30°C
Storage Condition: 2-30°C

使用條件: 18-28°C
Usage Condition: 18-28°C

