

Improve Better Patient Outcomes with Longer Lasting Restorations.

Increase office productivity and save time and material costs



OFFER VALID: APRIL 1, 2023 - JUNE 30, 2023

Specialty Carbides, Diamonds, Polishers and All Procedural Kits

BUY 4 GET 1 AT NO CHARGE**

CARBIDES: GREAT WHITE® ULTRA, ENDOGUIDE® AND ENDO SAFE END
DIAMONDS: PIRANHA®, REVELATION®, TDA AND GREAT WHITE®Z DIAMONDS
POLISHERS: JAZZ® POLISHERS

CAN MIX AND MATCH SHAPES TO EQUAL OR LESSER VALUE.
 NO CHARGE PRODUCT TO EQUAL OR LESSER VALUE.
 VALID FOR 3, 5, AND 25 PACKS OFFER EXCLUDES STANDARD CARBIDES, GREAT WHITE® GOLD AND TRIMMING AND FINISHING BURS.



STERILE R AVAILABLE IN ROLLS OF 25

Piranha® - Reduce sterilization costs & increase office productivity with a fresh diamond for each patient

Jazz® Polishers - Create faster, more beautiful polishes and eliminate up to 3 steps when finishing and polishing



All Lab Products

BUY 4 GET 1 AT NO CHARGE**

CHOOSE FROM: GREAT WHITE® LAB CUTTERS, JAZZ® LAB POLISHERS, EXPRESS LINE AND EXPRESS LINE-H DIAMOND DISCS, DURADI™ AND LAZER SHARP™ CAD/CAM BURS

Can mix and match shapes to equal or lesser value - No charge product to equal or lesser value

©2023 SS White® All Rights Reserved. SS White®, Express Line®, Great White® and Jazz® are registered trademarks of SS White Dental. DuraDi™ and Lazer Sharp™ are trademarks of SS White Dental. *Market research and test studies available upon request.**No charge goods must be of equal or lesser value. All offers are limited to stock on hand and are not valid on prior purchases. Offers may not be combined with any other sale, promotion, discount, code, coupon and/or offer. Sales promotion choice must be indicated on invoice sent to SS White® on applicable order. SS White® reserves the right to discontinue a promotion at any time. SS White's 100% guarantee is on all of our products. If you are not completely satisfied with your purchase, SS White® will issue a full refund of the purchase price through your distributor. To claim no charge goods: Mail a copy of your invoice by July 31, 2023 to: SS White Burs, Inc., Attn: Promotions Department, 1145 Towbin Ave., Lakewood, NJ 08701 or email freegoods@sswhitedental.com. Please note that these are end-user doctor promotions and must be redeemed by the qualifying practice. Please allow 4-6 weeks for free delivery. SS White® does its best to ensure that its photos/illustrations are as true to color as possible. However, due to inconsistencies of various monitors, lighting sources, digital photography, and manufacturing lot variations, SS White cannot guarantee that the color that appears on your screen/printed material accurately portrays the true color of the received product. Images are intended as a guide only.



Better Patient Outcomes
 Improved Efficiency
 Faster Practice Growth

1145 Towbin Avenue Lakewood, New Jersey 08701
 Visit us on the web at www.sswhitedental.com



Now, Save Even More Time with the Expanded SS White® Lineup of Zirconia Cutting Great White® Z Diamonds

Perfect for crown and bridge adjustment and removal, lab fabrication and creating endodontic access through aluminum oxide, lithium disilicate and zirconia copings.

STERILE R
AVAILABLE IN ROLLS OF 25

NEW PRODUCT



Description	368-023F	368-023M	379-023F	379-023M	801-010F	801-012F	801-014F	801-016M	801-018F	801-018M
Part Number	18165	18166	18164	18167	18168	18169	18162	18172	18163	18173
Sterile Part Number	-	-	10901 (25 Pk)	-	-	-	-	-	10902 (25 Pk)	-
Grit	F	M	F	M	F	F	F	M	F	M
Head Diameter (mm)	2.3	2.3	2.3	2.3	1.0	1.2	1.4	1.6	1.8	1.8
Head Length (mm)	5.0	5.0	4.2	4.3	1.0	1.2	1.3	1.6	1.7	1.7
LOA-Overall Length (mm)	19.0	19.0	19.0	19.0	19.0	19.0	19.0	19.0	19.0	19.0
Pack Size	5/Pk	5/Pk	5/Pk	5/Pk	5/Pk	5/Pk	5/Pk	5/Pk	5/Pk	5/Pk

NEW PRODUCT



Description	850-018F	850-018M	856-018F	856-018M	856-021F	856-021M	862-012F	862-012M	881-012F	881-016F
Part Number	18174	18175	18161	18176	18177	18178	18179	18181	18182	18183
Sterile Part Number	-	-	10900 (25 Pk)	-	-	-	-	-	-	-
Grit	F	M	F	M	F	M	F	M	F	F
Head Diameter (mm)	1.8	1.8	1.8	1.8	2.1	2.1	1.2	1.2	1.2	1.6
Head Length (mm)	10.0	10.0	8.0	8.0	8.0	8.0	8.0	8.0	8.0	8.0
LOA-Overall Length (mm)	23.0	23.0	22.0	22.0	22.0	22.0	22.0	22.0	22.0	22.0
Pack Size	5/Pk	5/Pk	5/Pk	5/Pk	5/Pk	5/Pk	5/Pk	5/Pk	5/Pk	5/Pk

The proprietary bonding process with fine grit synthetic diamond particles creates more consistent particle coverage, durability and cutting efficiency, while providing significant savings on material cost and chair time.


GREAT WHITE® Z DIAMONDS

**BUY 4
GET 1
AT NO CHARGE****


Endodontic Access Through Zirconia




Zirconia Crown Removal



Finish Laboratory Copings



Crown Adjustments



100% Satisfaction Guaranteed*

©2023 SS White® All Rights Reserved. SS White®, Express Line®, Great White® and Jazz® are registered trademarks of SS White Dental. DuraDI™ and Lazer Sharp™ are trademarks of SS White Dental. *Market research and test studies available upon request.**No charge goods must be of equal or lesser value. All offers are limited to stock on hand and are not valid on prior purchases. Offers may not be combined with any other sale, promotion, discount, code, coupon and/or offer. Sales promotion choice must be indicated on invoice sent to SS White® on applicable order. SS White® reserves the right to discontinue a promotion at any time. †SS White's 100% guarantee is on all of our products. If you are not completely satisfied with your purchase, SS White® will issue a full refund of the purchase price through your distributor. To claim no charge goods: Mail a copy of your invoice by July 31, 2023 to: SS White Burs, Inc., Attn: Promotions Department, 1145 Towbin Ave., Lakewood, NJ 08701 or email freegoods@sswhitedental.com. Please note that these are end-user doctor promotions and must be redeemed by the qualifying practice. Please allow 4-6 weeks for free delivery. SS White® does its best to ensure that its photos/illustrations are as true to color as possible. However, due to inconsistencies of various monitors, lighting sources, digital photography, and manufacturing lot variations, SS White cannot guarantee that the color that appears on your screen/printed material accurately portrays the true color of the received product. Images are intended as a guide only.

SS WHITE Dental®

Better Patient Outcomes
Improved Efficiency
Faster Practice Growth

1145 Towbin Avenue Lakewood, New Jersey 08701
Visit us on the web at www.sswhitedental.com