




Your Partner in Dentistry



ENDODONTIC

INSTRUMENT
MANAGEMENT

ORTHODONTIC

SURGICAL

DIAGNOSTIC

RESTORATIVE

PERIODONTAL

SUPPLIES &
CONSUMABLES

www.hiteckmedical.com

TABLE OF CONTENTS

ENDO EXPLORERS	03
ENDO EXCAVATORS	04
ENDO ACCESSORIES	05
PLUGGERS	06
Anterior Pluggers Posterior Pluggers ISO Sized Pluggers LUKS Pluggers	
RUBBER DAM CLAMPS.....	12
Anterior Rubber Dam Clamps Lower Molars Rubber Dam Clamp Premolars Molars General Upper Molars	
SPREADERS.....	16
RUBBER DAM FORCEPS, FRAMES & PUNCH	18
Rubber Dam Forceps Rubber Dam Punch Rubber Dam Frame	

ENDO EXPLORERS

ENDODONTIC



DG 16

1119

\$29.99



DG 16/17

1120

\$29.99



DG 16/23

1121

\$29.99

ENDO EXCAVATORS

HiTeck's Endodontic excavators feature extended-length shanks and are meticulously handcrafted from the highest grade stainless steel. They help the clinician access the pulp chamber for easy and efficient removal of dentin. The result is a superior clinical outcome and improved efficiency.

[1.5mm]



2 Glick

2405

\$29.99

[1.0mm]



31L Spoon

2406

\$29.99

[1.3mm]



31W Oval Spoon

2407

\$29.99

[1.5mm]



32L Spoon

2408

\$29.99

[2.0mm]



33L Spoon

2409

\$29.99

ENDO ACCESSORIES

ENDODONTIC

POST & SILVER POINT REMOVAL FORCEPS



Stieglitz

1295
\$48.99



Stieglitz Curved

1296
\$53.99



Stieglitz Angled

1297
\$55.99



Peet Curved

1298
\$39.99



Stainless Steel Ruler

ERULE6
\$25.99

6 Inches



Rubber Dam Clamp Organizing Board


RDCOB8
\$55.99

Store, sterilise and organise with eight pegs that securely hold your clamps. Clamps not included.
Dimensions: 4.62" x 2.4" x .375"
(117 mm x 61 mm x 10 mm). 8 Pieces.

ANTERIOR PLUGGERS

Anterior Root Canal Pluggers are used to compact filling material during vertical condensation. Marked at 5 mm intervals to assess penetration depth. Anterior condensers are 7 mm longer than posterior condensers.


[0.4mm]



8A 29mm **1270**

\$29.99


[0.55mm]



9A 29mm **1271**

\$29.99


[0.75mm]



10A 29mm **1272**

\$29.99


[1.05mm]



11A 29mm **1273**

\$29.99

[0.45mm]





8.5A 29mm **1274**

\$29.99

PLUGGERS




ANTERIOR PLUGGERS

Anterior Root Canal Pluggers are used to compact filling material during vertical condensation. Marked at 5 mm intervals to assess penetration depth. Anterior condensers are 7 mm longer than posterior condensers.



POSTERIOR PLUGGERS




Posterior Root Canal Pluggers used to compact filling material during vertical condensation. Marked at 5 mm intervals to assess penetration depth.

[0.4mm]		8 22mm	1279
			\$29.99
[0.55mm]		9 22mm	1280
			\$29.99
[0.75mm]		10 22mm	1281
			\$29.99
[1.05mm]		11 22mm	1282
			\$29.99
[0.45mm]		8.5 22mm	1283
			\$29.99

PLUGGERS








POSTERIOR PLUGGERS

Posterior Root Canal Pluggers used to compact filling material during vertical condensation. Marked at 5 mm intervals to assess penetration depth.



ISO SIZED PLUGGERS

ISO sized pluggers are used to compact filling material during vertical condensation.


[0.3mm]		ISO Size 30 (24mm)	1262 \$29.99
[0.4mm]		ISO Size 40 (24mm)	1263 \$29.99
[0.5mm]		ISO Size 50 (24mm)	1264 \$29.99
[0.6mm]		ISO Size 60 (24mm)	1265 \$29.99
[0.4/0.45mm]		1/3 Root Canal Plugger (21mm)	2257 \$29.99
[0.5/0.75mm]		5/7 Root Canal Plugger (21mm)	2351 \$29.99
[1.0/1.15mm]		9/11 Root Canal Plugger (21mm)	2352 \$29.99

PLUGGERS

LUKS PLUGGERS

LUKS pluggers are used to compact filling material during vertical condensation.


[0.45mm]



L1 18mm **1266**
\$29.99

The image shows the L1 18mm LUKS plugger, which has a 0.45mm tip diameter. It features a grey handle with a textured grip and a metal shaft with a hook-shaped tip. The HITECK logo is visible on the handle.


[0.50mm]



L2 18mm **1267**
\$29.99

The image shows the L2 18mm LUKS plugger, which has a 0.50mm tip diameter. It features a grey handle with a textured grip and a metal shaft with a hook-shaped tip. The HITECK logo is visible on the handle.


[0.75mm]



L3 18mm **1268**
\$29.99

The image shows the L3 18mm LUKS plugger, which has a 0.75mm tip diameter. It features a grey handle with a textured grip and a metal shaft with a hook-shaped tip. The HITECK logo is visible on the handle.

[1.00mm]









L4 18mm **1269**
\$29.99







The image shows the L4 18mm LUKS plugger, which has a 1.00mm tip diameter. It features a grey handle with a textured grip and a metal shaft with a hook-shaped tip. The HITECK logo is visible on the handle.




RUBBER DAM CLAMPS

ANTERIOR RUBBER DAM CLAMPS

					
0	6	9	9S	212	212SA
1307 \$12.99	1317 \$12.99	1322 \$12.99	1323 \$12.99	1345 \$12.99	1346 \$12.99
Flat jaws and high bow for small lower canines and first premolars	Slightly offset jaws for labial caries on anterior teeth	Offset flat jaws and stiff bow for anterior teeth	Offset flat serrated jaws for anterior teeth	Flat jaws for labial caries on anterior teeth	Slightly offset staggered jaws for odd shaped anteriors

LOWER MOLARS RUBBER DAM CLAMP

					
3	7	7A	12A	13A	26
1314 \$12.99	1318 \$12.99	1319 \$12.99	1326 \$12.99	1327 \$12.99	1333 \$12.99
Flat narrow jaws for small lower molars	Flat jaws for lower molars	Slightly reversed festooned jaws for large lower molars	Serrated offset jaws for lower right molars; paired with RDCM13A	Serrated offset jaws for lower left molars; paired with RCM12A	Slightly reversed festooned jaws for small lower molars

		
28	W3	W7
1337 \$12.99	1348 \$12.99	1350 \$12.99
Flat jaws for lower molars	Flat narrow jaws for small lower molars	Flat jaws for lower molars

RUBBER DAM CLAMPS

ENDODONTIC

PREMOLARS



0

1308
\$12.99

Flat jaws and high bow for longer premolars



1

1309
\$12.99

Slightly festooned jaws for upper premolars



1A

1310
\$12.99

Deeply festooned jaws for canines and premolars



2

1311
\$12.99

Flat jaws for premolars



2A

1312
\$12.99

Flat jaws for larger premolars



2AS

1313
\$12.99

Flat serrated jaws for larger premolars



27N

1336
\$12.99

Slightly festooned jaws for premolars



29

1338
\$12.99

Flat jaws for lower premolars



W2

1347
\$12.99

Flat jaws for premolars



209

2233
\$12.99

Small, flat jaws with stiff bow for lower premolars

RUBBER DAM CLAMPS

MOLARS GENERAL



8A

1321
\$12.99

Deeply festooned four-point jaws for molars



10

1324
\$12.99

Flat jaws designed for small upper left and lower right molars; paired with RDCM11



11

1325
\$12.99

Flat jaws designed for small upper right and lower left molars; paired with RDCM10



14

1328
\$12.99

Deeply festooned jaws for partially erupted or irregularly shaped molars



14A

1329
\$12.99

Deeply festooned jaws for partially erupted or irregularly shaped molars; larger than RDCM14



18

1330
\$12.99

Flat jaws for upper right and lower left molars; paired with RDCM25



24

1331
\$12.99

Flat jaws for upper left and lower right molars; paired RDCM24



25

1332
\$12.99

Flat jaws for large molars



27

1335
\$12.99

Slightly festooned jaws for large molars



56

1341
\$12.99

Slightly festooned serrated jaws for large partially erupted or irregularly shaped molars



56S

1342
\$12.99

Slightly festooned jaws for molars



W8A

1351
\$12.99

Deeply festooned jaws for partially erupted or irregularly shaped molars



W14A

2231
\$12.99

Deeply festooned jaws and extended bow for distal access to molars



8AD

2232
\$12.99

Deeply festooned jaws and extended bow for distal access to molars

RUBBER DAM CLAMPS

UPPER MOLARS



4

1315
\$12.99

Flat narrow jaws for small upper molars



5

1316
\$12.99

Slightly festooned jaws for large upper molars



8

1320
\$12.99

Deeply festooned jaws for upper molars



26N

1334
\$12.99

Flat jaws for upper molars



30

1339
\$12.99

Offset jaws for upper right molars; paired with RDCM31



31

1340
\$12.99

Offset jaws for upper left molars; paired with RDCM30



201

1343
\$12.99

Slightly festooned jaws for upper molars



205

1344
\$12.99

Slightly festooned jaws for large upper molars



W4

1349
\$12.99

Flat narrow jaws for small upper molars

[0.3mm]



1S Wakai 22mm

1288

\$29.99

[0.4mm]



2S Wakai 22mm

1289

\$29.99

[0.4mm]



D11 Standard 23mm

1290

\$29.99

[0.35mm]



D11S Thin 23mm

2228

\$29.99

[0.23mm]



D11T Thinner 23mm

2229

\$29.99

[0.20mm]




D11TS Thinnest 21mm

2230

\$29.99

SPREADERS

[0.25mm]




GP1 20mm

GP1

\$29.99

The image shows a GP1 20mm spreader, which is a curved instrument with a 0.25mm tip. It features a grey handle with a textured grip and the HITECK logo. The instrument is shown in a side view, highlighting its curved tip and the handle's design.

[0.30mm]




ISO Size 30 24mm

ISOS30

\$29.99

The image shows an ISO Size 30 24mm spreader, which is a curved instrument with a 0.30mm tip. It features a grey handle with a textured grip and the HITECK logo. The instrument is shown in a side view, highlighting its curved tip and the handle's design.

[0.40mm]




ISO Size 40 24mm

ISOS40

\$29.99

The image shows an ISO Size 40 24mm spreader, which is a curved instrument with a 0.40mm tip. It features a grey handle with a textured grip and the HITECK logo. The instrument is shown in a side view, highlighting its curved tip and the handle's design.

[0.50mm]



ISO Size 50 28mm

ISOS50

\$29.99

The image shows an ISO Size 50 28mm spreader, which is a curved instrument with a 0.50mm tip. It features a grey handle with a textured grip and the HITECK logo. The instrument is shown in a side view, highlighting its curved tip and the handle's design.

[0.60mm]



ISO Size 60 28mm

ISOS60

\$29.99

The image shows an ISO Size 60 28mm spreader, which is a curved instrument with a 0.60mm tip. It features a grey handle with a textured grip and the HITECK logo. The instrument is shown in a side view, highlighting its curved tip and the handle's design.

RUBBER DAM FORCEPS, FRAMES & PUNCH

RUBBER DAM FORCEPS



Ivory Lightweight

1299
\$65.99
17cm



Brewer

1300
\$75.99
17cm



U of Washington

1301
\$75.99
17cm

RUBBER DAM FORCEPS, FRAMES & PUNCH

RUBBER DAM PUNCH



Ivory
1302
\$95.99



Ainsworth
1303
\$95.99

RUBBER DAM FRAME



Small Child

1304
\$25.99
3.5 inch/9cm

Child

1305
\$25.99
5 inch/13cm

Adult

1306
\$25.99
6 inch/15cm

OUR OTHER PRODUCT CATEGORIES



DIAGNOSTIC



RESTORATIVE




SURGICAL



PERIODONTAL



ORTHODONTIC



INSTRUMENT
MANAGEMENT



ENDODONTIC



SUPPLIES &
CONSUMABLES

To Learn more about HiTeck's Products & Services visit hiteckmedical.com

Manufacturer: 18-215 Traders Blvd East | Mississauga, ON L4Z 3K5 | Canada

Toll Free: 877.899.1861 **Toll Free Fax:** 888.891.9792

Website: www.hiteckmedical.com

Copyright 2020 HiTeck Medical Instruments. All rights reserved. The contents of this product catalog including, but not limited to, all text, photographs and graphics, are the sole and exclusive property of HiTeck Medical Instruments and may not be copied, reproduced or distributed without prior written permission.

HiTeck reserves the right to change features of products, prices, images and text at any time without notice.

For warranty information, please see www.hiteckmedical.com

HiTeck

Your Partner in Dentistry