




PERIODONTAL




ORTHODONTIC



DIAGNOSTIC



INSTRUMENT
MANAGEMENT



ENDODONTIC



RESTORATIVE



SURGICAL

PROCEDURAL KITS


Hand picked By Experts

Compilation of clinicians' most popular procedural set-ups. When you order a Collection, you receive high quality HiTeck instruments and cassettes together in procedural set-ups at a reduced price.

SUMMARY OF KITS

HYGIENE KIT	\$199.99
EXAM KIT	\$74.99
ORTHO ARCHWIRE ADJUSTMENT KIT	\$799.99
SIYAH LUXATING HYBRID KIT	\$499.95
SERRATED LUCAS CURETTE KIT	\$389.99
SIYAH PERIOTOME KIT	\$399.99
ORTHO BANDING/BONDING KIT	\$549.99
ATRAUMATIC EXTRACTION KIT	\$899.99
NDECC EXAM INSTRUMENT KIT	\$1049.99
SUTURE REMOVAL KIT	\$269.99
NEW PATIENT ASSESSMENT KIT	\$189.99
PROPHY STANDARD KIT	\$369.99
ORAL SURGERY EXTRACTION KIT	\$1,799.99

HYGIENE KIT



H6/H7

204S

11/12
GRACEY

13/14
GRACEY

4R/4L
COLUMBIA

13/14
COLUMBIA

ODU 11/12

23/QW
EXPRO

MIRROR
HANDLE

COLLEGE
PLIER

Special Price

\$199.99
/set

Regular Price
\$273.99
/set

Special Offer

Add on a Sterilization
Cassette for only:

\$59.99 **\$94.99**
Daily Line Premium
Series

- Comes with Lifetime Warranty!
- Made from High Grade Surgical Steel
- Ergonomic Handle for user comfort

EXAM KIT



ODU 11/12



23/QW
EXPRO



MIRROR
HANDLE



COLLEGE
PLIER

Special Price
\$74.99
/set

Regular Price
\$96.99
/set

Special Offer

Add on a Sterilization
Cassette for only:

\$52.99 **\$82.99**

Daily Line Premium
Series

- Comes with Lifetime Warranty!
- Made from High Grade Surgical Steel
- Ergonomic Handle for user comfort


ORTHO ARCHWIRE ADJUSTMENT STANDARD KIT

Special Set Price

\$799.99

Regular:

\$1,149.99



LIGATURE
DIRECTOR & BAND
PUSHER



#5 CONE SOCKET
SIYAH MIRROR - HD



MICRO MINI
PIN &
LIGATURE
CUTTER



WEINGART
PLIERS



UNIVERSAL CUT &
HOLD DISTAL END
CUTTER



BIRD BEAK
PLIERS



HARD WIRE
CUTTER 15
DEGREES



THREE JAW
PLIERS




SCALER/
BAND
PUSHER




SIGNATURE
MIRROR
HANDLE



ORTHO SLIM TIP
MATHIU NEEDLE
HOLDER



ORTHO SLIM
CASSETTE



DIRECT BOND
BRACKET
HOLDER

LUXATING HYBRID ELEVATOR SET

Set of 5

\$499.99

Regular:

\$749.95

Special Offer

Buy a Complete set of
13 for:

\$1099.99

Regular:

\$1949.87



- ✓ A unique cross between a periotome and a luxating elevator, Luxating Hybrids are used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.
- ✓ Full-control during apical progression of the sharp working-edge
- ✓ Preserves the socket while gently removing the tooth
- ✓ Tapered, one-piece design drives tip apically
- ✓ Allows for palm pressure to advance the instrument further

SERRATED LUCAS CURETTE KIT

Special Set Price

\$389.99

Regular:

\$499.99




3.5MM

2.8MM

2.5MM

2MM




5
INSTRUMENT
CASSETTE

- ✓ The cutting property is improved by the tothing, which acts as a "saw" for cutting and removal. Sharper than standard LUCAS curettes, serrated LUCAS curettes are used for more invasive applications, which makes it possible to act on tougher fabrics and provides a better deep curettage of the alveolus.

PERIOTOME KIT

Periotomes are used in a procedure that facilitates the removal of teeth with minimal damage to the surrounding alveolar bone. This procedure may be used whenever an extraction is indicated, but is particularly important when considering the placement of a dental implant.



PT6 Anterior Straight **PT5** Anterior Angled **PT4** Posterior Left **PT3** Posterior Right **PT2** Anterior **PT1** Posterior

Set of 6
\$399.99
Regular:
\$509.94

Special Offer

Add on a Sterilization
Cassette for only:

\$59.99 **\$94.99**
Daily Line Premium
Series

- Thin, clinically-designed blades sever periodontal ligament easily
- Blades are permanently fixed to handle and will not rotate under pressure
- Siyah coating for the sharpest and longest lasting cutting edges

ORTHO BANDING/BONDING STANDARD SET-UP




#5 CONE SOCKET
SIYAH MIRROR - HD



LIGATURE
DIRECTOR & BAND
PUSHER

ORTHO SLIM
CASSETTE



MICRO MINI PIN &
LIGATURE CUTTER



UNIVERSAL CUT &
HOLD DISTAL END
CUTTER



COLLEGE
TWEEZER

Special Set Price
\$549.99
Regular:
\$749.99



ORTHO SLIM TIP
MATHIU NEEDLE
HOLDER



DIRECT BOND
BRACKET HOLDER



SIGNATURE
MIRROR HANDLE



1/2 HOLLENBACK
CARVER

ATRAUMATIC EXTRACTION KIT

HiTeck's Atraumatic Extraction Kit comprises of specially chosen instruments to give you everything you need to perform your extraction. These instruments are specially designed to help you perform extraction, atraumatically!

Special Set Price

\$899.99

Regular:

\$1,259.99



Apical Forceps 150 & 151

Thin, tapered and sharp parallel beaks improve root and crown contact and longitudinal serrations provide additional grip. The design of Apical Forceps facilitates an atraumatic tooth extraction.



2/4 Siyah Molt Curette

Designed for curettage, cyst removal and tooth socket debridement with 3.6mm and 7mm working ends.



PUSH Scalpel Handle

Used for cutting gingival tissue and making surgical incisions. Provides greater rotational movement with fingertip control and enhances dexterity. Has a push mechanism to safely remove the blade.



Crile-Wood Needle Holder

Used to grasp and guide the needle when suturing. Recommended suture sizes: 3-0, 4-0, 5-0, 6-0



16 Instrument Sterilization Cassette

Made of long-lasting, low maintenance stainless steel. Soft, colored silicone rails securely hold instruments throughout the processing cycle.



9 Siyah Molt Periosteal

Designed for reflecting and retracting the mucoperiosteum after incisions of the gingival tissue. 4.3mm/8.4mm



Goldman-Fox Scissor, Straight

Used for trimming tissue or cutting sutures. Featuring tapered blade tips, tips that are not sharp on outside and are used as tenotomy scissors. 12.5cm with 1 serrated blade



Small Straight, Small Curved, Straight Spade, Siyah Luxating Hybrid Elevator

A unique cross between a periosteal and a luxating elevator, used apically within the periodontal ligament space to cut periodontal ligaments prior to extraction.

NDECC EXAM INSTRUMENT KIT

Special Set Price

\$1049.99

Regular:

\$1,399.99



#5 CONE SOCKET SIYAH MIRROR - HD

SIGNATURE MIRROR HANDLE

27/29 FOOTBALL BURNISHER

21B ACORN BURNISHER

CR21 KNIFE

15/16 HATCHET

28 DISTAL & 29 MESIAL MARGIN TRIMMERS

1/2 HOLLENBACK CARVER

89/92 DISCOID CLEOID

H1 HOLLENBACK PLUGGER

REGULAR/ LARGE CARRIER

LOCKING TWEEZER

COLLEGE TWEEZER

23QW WILLIAMS EXPRO



GOLDMAN FOX CURVED




KELLY CURVED



MILLER FORCEP

RING PLACEMENT FORCEP



CASSETTE FOR 20 INSTRUMENTS



CASSETTE FOR 16 INSTRUMENTS PLUS ACCESSORY



NON SLIP AMALGAM WELL



32L EXCAVATOR



DG 16 EXPLORER



24 SPATULA



3 LADMORE SIYAH BURNISHER



#3 SIYAH GOLDSTEIN PFI



CVIPC INTERPROXIMAL CARVER

SUTURE REMOVAL STANDARD SET-UP

Special Set Price

\$269.99

Regular:

\$364.94



10 INSTRUMENT CASSETTE



SIGNATURE MIRROR HANDLE




#5 SIYAH EXPLORER



COLLEGE TWEEZER



#5 CONE SOCKET SIYAH MIRROR - HD



IRIS SCISSOR TUNGSTEN CARBIDE

NEW PATIENT ASSESSMENT KIT

Special Set Price


\$189.99

Regular:

\$299.99



#5 CONE SOCKET
SIYAH MIRROR - HD



5 INSTRUMENT
CASSETTE



#5 SIYAH
EXPLORER



SIGNATURE
MIRROR
HANDLE



23/UNC12
EXPRO



ODU 11/12
EXPLORER



2N NABERS
PROBE

PROPHY STANDARD SET-UP

Special Set Price

\$369.99

Regular:

\$479.99



#5 CONE SOCKET
SIYAH MIRROR - HD



10 INSTRUMENT
CASSETTE



13/14
COLUMBIA



13/14
GRACEY



11/12
GRACEY



H6/H7
SCALER



CP 12
PROBE



ODU 11/12
EXPLORER



SIGNATURE
MIRROR
HANDLE

ORAL SURGERY EXTRACTION KIT

Special Set Price

\$1,799.99

Regular:

\$2,799.99



8 CRANE PICK



1C COUPLAND



A30 CRYERS




A31



69W WOODWARD CHISEL



OLSEN HEGAR



ASPIRATING SYRINGE



KELLY CURVED SCISSOR



MEAD MALLET



7.5MM OSTEOTOME STRAIGHT



4MM OSTEOTOME STRAIGHT



2 GARDNER CHISEL



23 SELDIN



#3 BLADE HANDLE




GAUZE PACKER



#11 MILLER CURETTE



#73 LOWER MOLAR



KELLY CURVED SERRATED



ALLIS TISSUE FORCEP



16 PC STERILIZATION CASSETTE

visit us online at hiteckmedical.com



Your Partner in Dentistry

OUR OTHER PRODUCT CATEGORIES



DIAGNOSTIC



RESTORATIVE



SURGICAL



PERIODONTAL



ORTHODONTIC



INSTRUMENT
MANAGEMENT



ENDODONTIC



IMPLANTS

To Learn more about HiTeck's Products & Services, visit hiteckmedical.com

Manufacturer: 18-215 Traders Blvd. East | Mississauga, ON L4Z 3K5 | Canada

Toll Free: 877.899.1861 **Toll Free Fax:** 888.891.9792

Website: www.hiteckmedical.com

Copyright 2022 HiTeck Medical Instruments. All rights reserved. The contents of this product catalog including, but not limited to, all text, photographs and graphics, are the sole and exclusive property of HiTeck Medical Instruments and may not be copied, reproduced or distributed without prior written permission.

HiTeck reserves the right to change features of products, prices, images and text at any time without notice.

For warranty information, please see www.hiteckmedical.com