



Your Partner in Dentistry



MOST POPULAR DENTAL INSTRUMENTS

www.hiteckmedical.com

SURGICAL

APICAL FORCEPS	4
ATRAUMAIR FORCEPS	4
SIYAH DEEP GRIP FORCEPS	4
STANDARD FORCEPS	5
PEDO FORCEPS	6
SURGICAL ELEVATORS	6
LUXATING ELEVATORS	10
LUXATING ELEVATORS - SIYAH SERIES	10
LUXATING HYBRID ELEVATORS	11
SCISSORS	12
NEEDLE DRIVERS	13
ASPIRATORS	14
BONE RONGEURS	14
CASTROVIEJO CALIPERS	15
HEMOSTATS & CLAMPS	15
MALLET	15
DRESSING & TISSUE FORCEPS	16
OSTEOTOMES	16
PERIOSTEALS	17
PERIOTOMES	18
RETRACTORS	19
ROOT TIP PICKS	20
SCALPEL HANDLES	20
BONE FILES	21
BONE CHISELS	21
PERIODONTAL BONE CHISELS	21
CURETTES	22
PERIODONTAL KNIVES	23
SINUS LIFT	23
SYRINGE	23
WIRE	23

DIAGNOSTIC

TWEEZERS	24
MIRRORS & HANDLES	24
SIYAH EXPLORERS	24
EXPLORERS	25
PROBES	25
EXPROS	27

PERIODONTAL

SICKLE SCALERS	28
JACQUETTE SCALERS	28
MCCALL SCALER/CURETTE COMBO	29
GOLDMAN FOX SCALER/CURETTE COMBO	29
GOLDMAN FOX SET	30
LANGER UNIVERSAL CURETTES	30
COLUMBIA UNIVERSAL CURETTES	30
BARNHART UNIVERSAL CURETTES	31
TAHA SCALERS	31
GRACEY CURETTES	32
MINI GRACEY CURETTES	33

ENDODONTIC

ENDO EXPLORERS	34
ENDO EXCAVATORS	34
ANTERIOR PLUGGERS	34
POSTERIOR PLUGGERS	35
ROOT FRAGMENT FORCEPS	35
ENDO ACCESSORIES	35
ISO SIZED PLUGGERS	36
LUKS PLUGGERS	36
RUBBER DAM CLAMPS	37
RUBBER DAM FORCEPS	39
RUBBER DAM PUNCH	39
RUBBER DAM FRAMES	39

RESTORATIVE

GOLDSTEIN SIYAH COMPOSITE	40
SIYAH COMPOSITE INSTRUMENTS	41
SIYAH BURNISHERS	42
GOLDFOGEL FREEHAND INSTRUMENTS	43
STANDARD PFI'S	44
HATCHETS	45
MARGIN TRIMMERS	45
CHISELS	45
PLACEMENT INSTRUMENT	45
BLADE EXCAVATORS	46
ENGLISH PATTERN EXCAVATORS	46
OVAL EXCAVATORS	46
SPOON EXCAVATORS	47
AMALGAM WELL	47
ARTICULATING FORCEP.....	47

PLUGGERS/CONDENSERS	48
CARVERS	50
DISCOID CLEOID	50
ACORN BURNISHERS	51
BALL BURNISHERS	51
BEAVERTAIL BURNISHERS	51
LADMORE BURNISHERS	51
CHRISTENSEN CROWN REMOVERS	52
GOLDSTEIN CROWN REMOVERS	52
TEMPORARY CROWN REMOVERS	52
SPATULAS & KNIFE	53
MEASURING DEVICES	53
P.K THOMAS INSTRUMENTS	54

ORTHODONTIC

CLEAR ALIGNER INSTRUMENTS	55
CUTTERS	55
UTILITY PLIERS	56
WIRE FORMING PLIERS	57
STEEL PHOTOGRAPHY MIRRORS	58
HAND INSTRUMENTS	59

INSTRUMENT MANAGEMENT

PREMIUM SERIES CASSETTES	61
DAILY LINE CASSETTES	66
ACCESSORY CASSETTES.....	67

APICAL FORCEPS



AF1 Standard Apical
Upper anteriors
2524 | \$119.99



AF150 Apical
Upper universal
2525 | \$119.99



AF151 Apical
Lower universal
2526 | \$119.99



AF222 Apical
Lower molars
2527 | \$119.99



AF67A Apical
Upper molars
2528 | \$119.99



AF74N Apical
Lower anteriors
2529 | \$119.99

ATRAUMAIR FORCEPS



67
Upper Molars
F67XS | \$149.99



1 STANDARD
Upper Incisors
F1XS | \$149.99



1 Apical
Upper Incisors
FAF1XS | \$149.99



35
Upper Canines & Premolars Molars
F35XS | \$149.99



36
Lower Incisors
F36XS | \$149.99



36 Apical
Lower Premolars & Incisors
FAF36XS | \$149.99



13
Lower Premolars
F13XS | \$149.99



51 Root Forceps
Upper roots, serrated
FX51XS | \$149.99



79
Lower molars
FX79XS | \$149.99



73
Lower molars
FX73XS | \$149.99



23 Cowhorn
Lower molars
F23XS | \$149.99



45 Roots
Upper roots, serrated
FX45XS | \$149.99

SIYAH DEEP GRIP FORCEPS



Upper Anterior
SEF902 | \$279.99



Upper Universal
SEF903 | \$279.99



Lower Universal
SEF904 | \$279.99



Universal Upper Molars
SEF905 | \$279.99



AF1 Standard Apical
Upper Canines & Premolars Molars
FAF35XS | \$149.99



Lower Anteriors
Lower anteriors
SEF906 | \$279.99



Lower Molars
Lower molars
SEF907 | \$279.99



Lower Roots
Serrated spade tip
SEF908 | \$279.99



Upper Roots
Serrated spade tip
SEF909 | \$279.99

The Siyah Deep Grip Forceps are coated with Aluminium Titanium Nitride that is significantly harder, scratch-resistant and non-reflective for long lasting performance.

Increased hardness ensures long-term integrity of both the beak serrations and forceps overall.

- Extremely smooth surface seals the pores of the stainless steel, minimizing protein adhesion and reducing time to clean.
- Matte finish that reduces distracting glare and reflections mid-procedure.
- Machined beaks give you precision and control over your extractions.



13 Serrated
Lower Premolars
1423 | \$109.99



151 Cryer Universal
Lower Incisors, canines & premolars & roots
1498 | \$109.99



MD3 Mead Serrated
Lower incisors & roots
EF06 | \$109.99



22 Serrated
Lower molars
1424 | \$109.99



16 Cowhorn
1st & 2nd lower molars
1494 | \$109.99



73 Serrated
Lower molars
1435 | \$109.99



79 Serrated
Lower 3rd molars
1440 | \$109.99



23 Cowhorn
1st & 2nd lower molars
1497 | \$109.99



36 Serrated
Special for lower premolars
FX36E | \$109.99



49 Serrated
Upper roots
1410 | \$109.99



51 Serrated
Upper roots
1411 | \$109.99



97 Serrated
Upper anteriors
1419 | \$109.99



46L Serrated
Lower roots
1433 | \$109.99



74N Narrow Beaks Serrated
Lower roots
1437 | \$109.99



86A
Lower molars with carious or broken crowns Forceps
1443 | \$109.99



76S Serrated
Upper roots
2532 | \$109.99



300 Serrated
Upper roots
2530 | \$109.99



301 Serrated
Lower roots
2531 | \$109.99



7 Serrated
Upper premolars
1403 | \$109.99



32 Parmly Alveolar
Upper canines, premolars, & molars
1478 | \$109.99



150 Cryer Universal
Upper incisors & canines, premolars, & roots universal
1487 | \$109.99



35
Upper incisors, canines, premolars & roots
EF19 | \$109.99



286
Upper roots, incisors & premolars
EF17 | \$109.99



17 Serrated
Upper molars; right
1404 | \$109.99



18 Serrated
Upper molars; Left
1405 | \$109.99



18L Harris
Forceps for 1st & upper 2nd molars
1476 | \$109.99



18R Harris
Forceps for 1st & upper 2nd molars
1477 | \$109.99



67A Serrated
Upper 3rd molars
1415 | \$109.99



88L Nevius
Forceps for 1st & upper 2nd molars
1485 | \$109.99



88R Nevius
Forceps for 1st & upper 2nd molars
1486 | \$109.99

PEDO FORCEPS



1C Pedo
Lower incisors
FX1CE | \$109.99



2C Pedo
Lower premolars
FX2CE | \$109.99



3C Pedo
Lower roots
FX3CE | \$109.99



4C Pedo
Upper premolars
FX4CE | \$109.99



5C Pedo
Upper incisors & canines
FX5CE | \$109.99



6C Pedo
Upper molars
FX6CE | \$109.99



150S Pedo
Upper primary teeth & roots universal
F150S | \$109.99



151S Pedo
Lower primary teeth & roots universal
F151S | \$109.99

SURGICAL ELEVATORS



301
E3 | \$74.99



4 Schmeckebeier
Displacement of broken or deep seated root tips.
E28 | \$74.99



5 Schmeckebeier
Displacement of broken or deep seated root tips.
E29 | \$74.99



71 Miller
Displacement of broken or deep seated root tips.
E30 | \$74.99



72 Miller
Displacement of broken or deep seated root tips.
E31 | \$74.99



302
E34 | \$74.99



303
E35 | \$74.99



73 Miller
E32 | \$74.99



74 Miller
E33 | \$74.99



1C
F13XS | \$74.99



2C
F13XS | \$74.99



3C
F13XS | \$74.99



A1C English Pattern
F13XS | \$74.99



A2C English Pattern
F13XS | \$74.99



A3C English Pattern
F13XS | \$74.99



21 Cryer Mini
E5 | \$74.99



22 Cryer Mini
E6 | \$74.99



44 Cryer Small
E7 | \$74.99



45 Cryer Small
E8 | \$74.99



25 Cryer Medium
10.7 mm wide tip
E9 | \$74.99



26 Cryer Medium
10.7 mm wide tip
E10 | \$74.99



31 Cryer Large
E11 | \$74.99



32 Cryer Large
E12 | \$74.99



A30 Cryer
E13 | \$74.99



A31 Cryer
E14 | \$74.99



77R Serrated
E2 | \$74.99



B Cogswell
E4 | \$74.99



6 Potts
E42 | \$74.99



7 Potts
E43 | \$74.99



8 Crane Pick
E44 | \$74.99



11 Stout 3rd Molars
E45 | \$74.99



11A Stout 3rd Molars
E46 | \$74.99



11M MacMillan Gouge
E47 | \$74.99



12M MacMillan Gouge
E48 | \$74.99



13M MacMillan Gouge
E49 | \$74.99



15 Woodward
E50 | \$74.99



16 Woodward
E51 | \$74.99



19
E52 | \$74.99



20
E53 | \$74.99



23
E54 | \$74.99



31F Friedman [Craig]
E55 | \$74.99



32F Friedman [Craig]
E56 | \$74.99



38 Serrated
E57 | \$74.99



39 Serrated
E58 | \$74.99



44C Bradley
E59 | \$74.99



46
E60 | \$74.99



46R
E61 | \$74.99



51 Curtis
E62 | \$74.99



52 Curtis
E63 | \$74.99



63
E64 | \$74.99



69W Woodward
E65 | \$74.99



77 Serrated
E66 | \$74.99



80
E67 | \$74.99



81
E68 | \$74.99



81A
E69 | \$74.99



81E
E70 | \$74.99



92 Serrated
E71 | \$74.99



93 Serrated
E72 | \$74.99



98
E73 | \$74.99



190
E74 | \$74.99



191
E75 | \$74.99



A Cogswell
E76 | \$74.99



H5 Hourigan Serrated
E77 | \$74.99



H6 Hourigan Serrated
E78 | \$74.99



6 Potts T-bar Handle
E79 | \$80.00



7 Potts T-bar Handle
E80 | \$80.00



1 Bein, Large
E124 | \$74.99



2 Bein, Medium
E125 | \$74.99



3 Bein, Small
E126 | \$74.99



20 Elevator
E53 | \$74.99



Lindo Levian Root, Serrated
1528 | \$74.99



1 Heidbrink
E15 | \$74.99



2 Heidbrink
E16 | \$74.99



40 Heidbrink
E17 | \$74.99



41 Heidbrink
E18 | \$74.99



1L Seldin
E19 | \$74.99



1R Seldin
E20 | \$74.99



2S Seldin
E21 | \$74.99



3S Seldin
E22 | \$74.99



34 Seldin
E23 | \$74.99



34S Seldin
E1 | \$74.99



4L Seldin
E24 | \$74.99



4R Seldin
E25 | \$74.99



301W Seldin
E26 | \$74.99



304W Seldin
E27 | \$74.99



Warwick-James Straight
E81 | \$80.00



Warwick-James Left
E82 | \$80.00



Warwick-James Right
E83 | \$80.00



1LX
E84 | \$74.99



1RX
E85 | \$74.99



11LX
E86 | \$74.99



11RX
E87 | \$74.99



12LX
E88 | \$74.99



12RX
E89 | \$74.99



14LX
E90 | \$74.99



14RX
E91 | \$74.99



11 Bernard Straight
E91 | \$80.00



6 Bernard Right
E91 | \$80.00



7 Bernard Left
E91 | \$80.00



Curved, 2 mm
E92 | Regular | \$74.99
E92S | Serrated | \$74.99



Straight, 2 mm
E93 | Regular | \$74.99
E93S | Serrated | \$74.99



Curved, 3 mm
E94 | Regular | \$74.99
E94S | Serrated | \$74.99



Straight, 3 mm
E95 | Regular | \$74.99
E95S | Serrated | \$74.99



Curved, 4 mm
E96 | Regular | \$74.99
E96S | Serrated | \$74.99



Straight, 4 mm
E97 | Regular | \$74.99
E97S | Serrated | \$74.99



Curved, 5 mm
E98 | Regular | \$74.99
E98S | Serrated | \$74.99



Straight, 5 mm
E99 | Regular | \$74.99
E99S | Serrated | \$74.99



Curved, 2 mm
EL2CX | Regular | \$98.99
EL2CXS | Serrated | \$98.99
SEL2CX | Slim Long | \$98.99
SSEL2CX | Slim Short | \$98.99



Straight, 2 mm
EL2SX | Regular | \$98.99
EL2SXS | Serrated | \$98.99
SEL2SX | Slim Long | \$98.99
SSEL2SX | Slim Short | \$98.99



Curved, 3 mm
EL3CX | Regular | \$98.99
EL3CXS | Serrated | \$98.99
SEL3CX | Slim Long | \$98.99
SSEL3CX | Slim Short | \$98.99



Straight, 3 mm
EL3SX | Regular | \$98.99
EL3SXS | Serrated | \$98.99
SEL3SX | Slim Long | \$98.99
SSEL3SX | Slim Short | \$98.99



Curved, 4 mm
EL4CX | Regular | \$98.99
EL4CXS | Serrated | \$98.99
SEL4CX | Slim Long | \$98.99
SSEL4CX | Slim Short | \$98.99



Straight, 4 mm
EL4SX | Regular | \$98.99
EL4SXS | Serrated | \$98.99
SEL4SX | Slim Long | \$98.99
SSEL4SX | Slim Short | \$98.99



Curved, 5 mm
EL5CX | Regular | \$98.99
EL5CXS | Serrated | \$98.99
SEL5CX | Slim Long | \$98.99
SSEL5CX | Slim Short | \$98.99



Straight, 5 mm
EL5SX | Regular | \$98.99
EL5SXS | Serrated | \$98.99
SEL5SX | Slim Long | \$98.99
SSEL5SX | Slim Short | \$98.99



Small, Curved
E100 | \$149.99

Small, Straight
E101 | \$149.99

Small, Distal
E102 | \$149.99

Small, Mesial
E103 | \$149.99

Large, Distal
E104 | \$149.99

Large, Mesial
E105 | \$149.99

Straight Spade
E117 | \$149.99

Distal Spade
E118 | \$149.99

Mesial Spade
E119 | \$149.99

Angled Right Spade
E120 | \$149.99

Angled Left Spade
E121 | \$149.99

Straight, Large
4.5mm
E122 | \$149.99

Curved, Large
4.5mm
E123 | \$149.99

Set of 13 Elevators with Cassette
HYBRIDLUX13 | \$1849



Dean, Serrated
05-1652-17 | 17 CM | \$55.99



Kelly, Serrated
05-1454-16 | 16 CM Straight | \$52.99
05-1455-16 | 18 CM Curved | \$52.99



Iris
05-1484-11 | 11.5 CM Straight | \$32.99
04-1485-11 | 11.5 CM Curved | \$32.99



Iris, Tungsten Carbide
05-1842-11 | 11.5 CM Straight | \$87.99
05-1843-11 | 11.5 CM Curved | \$87.99



Goldman Fox Serrated
05-1840-13 | 13 CM Straight | \$97.99
05-1841-13 | 13 CM Curved | \$97.99



Dean Serrated
200-133 | 16.5 CM | \$109.99



Locklin Curved, Serrated
200-134 | 16.5 CM | \$97.99



Kelly Serrated
200-131 | 16 CM Curved | \$98.99
200-132 | 16 CM Straight | \$98.99



Castroviejo Curved
S31 | 10 CM | \$174.99



Castroviejo Curved
SPV | 18 CM | \$194.99



Iris
05-1843-11C | 11.5 CM Curved | \$129.99
05-1842-11C | 11.5 CM Straight | \$129.99




Mathieu
 1753 | 14 CM | \$44.99
 200-121 | 14 CM Tungsten Carbide | \$94.99



Crile Wood
 09-2515-15 | 15 CM with Groove | \$44.99



Mayo Hegar
 09-2511-16 | 16 CM | \$44.99



Olsen Hegar
 09-2693-14 | 14 CM | \$114.99
 09-2694-17 | 17 CM | \$124.99
 09-2692-11 | 11 CM | \$104.99



Mayo Hegar
 09-2652-16 | 16 CM | \$124.99
 09-2652-18 | 18 CM | \$124.99



Crile Wood
 09-2661-15 | 15 CM | \$89.99
 09-2661-18 | 18 CM | \$89.99



HMT Thin Jaw
 200-118 | 15 CM | \$89.99



Micro Vascular
 200-123 | 18 CM | \$124.99



Barraquer Needle Holder
 2512 | 16 CM, Straight | \$189.99
 2513 | 18 CM, Straight | \$214.99
 2514 | 16 CM, Curved | \$189.99
 2515 | 18 CM, Curved | \$214.99



Yankauer

- 1359** | 20 CM | **\$74.99**
- 1360** | 27 CM | **\$88.99**

Frazier Suction Tube

- 1362** | 1.5 mm \circ | **\$74.99**
- 1363** | 2.0 mm \circ | **\$74.99**
- 1364** | 2.5 mm \circ | **\$74.99**
- 1365** | 3.0 mm \circ | **\$74.99**



15CM Mini Friedman Bone Rongeur

- 1599** | Straight | **\$139.99**
- 1600** | 30° | **\$139.99**
- 1601** | 45° | **\$139.99**
- 1602** | 90° | **\$139.99**



Friedman Bone Rongeur

- 1603** | 30° | **\$139.99**
- 1604** | 45° | **\$139.99**



30° Cleaveland Side Cutting

- 1614** | 14CM | **\$139.99**
- 1615** | 17CM | **\$139.99**



15CM Blumenthal

- 1621** | 30° | **\$139.99**
- 1622** | 45° | **\$139.99**
- 1623** | 90° | **\$139.99**



30° Cleaveland Bone Rongeur

- 1612** | 14CM | **\$139.99**
- 1613** | 17CM | **\$139.99**



90MM Shaft Kerrison Shortnose

- 1618** | 2MM Tip | **\$244.99**
- 1600** | 3 MM Tip | **\$244.99**
- 1620** | 5 MM Tip | **\$244.99**



15CM Bohler Bone Rongeur

- 1609** | Straight | **\$139.99**



12CM Mini Friedman Bone Rongeur

- 1605** | 30° | **\$139.99**



18 CM Weiss Modified
CAL1 | Straight | **\$164.99**
CAL2 | Curved | **\$164.99**



18 CM Standard
CAL3 | Straight | **\$164.99**
CAL4 | Curved | **\$164.99**



11 CM Rounded Castroviejo
CAL5 | Straight | **\$114.99**
CAL6 | Curved | **\$114.99**



18 CM Rounded Castroviejo
CAL7 | Straight | **\$164.99**
CAL8 | Curved | **\$164.99**



16 CM Kelly Rankin Artery Forcep
07-2130-16 | Straight | **\$29.99**
07-2131-16 | Curved | **\$29.99**



12 CM Halstead Mosquito Forcep
07-2113-12 | Straight | **\$28.99**
07-2114-12 | Curved | **\$28.99**



12 CM Halstead Mosquito Forcep
07-2115-12 | Straight | **\$28.99**
07-2116-12 | Curved | **\$28.99**



14 CM Ochsner Kocher Artery Forcep
07-2156-14 | Straight | **\$29.99**
07-2157-14 | Curved | **\$29.99**



16 CM Rochester Pean Fine Artery Forcep
07-2154-16 | Straight | **\$29.99**
07-2155-16 | Curved | **\$29.99**



18 CM Adson Artery Forcep
07-2177-18 | Straight | **\$35.99**
07-2178-18 | Curved | **\$35.99**

#1 Mead Mallet
MAL1 | **\$89.99**





12 CM Adson Forcep

- 06-1919-12 | Regular Dressing Serrated | \$33.99
- 06-1981-12 | Regular Tissue 1x2 Teeth | \$33.99



12 CM Adson Brown

- 06-1984-12 | 5x6 Teeth | \$33.99



18 CM Gerald Tissue

- 06-1970-18 | 1x2 Teeth, Straight | \$38.99
- 06-2047-18 | 1x2 Teeth, Curved | \$38.99
- 06-1923-18 | Serrated, Straight | \$38.99
- 06-1924-18 | Serrated, Curved | \$38.99



Allis Tissue

- 26-13502-15 | 5x6 Teeth | \$34.99



Corn Suture Plier

- 200-110 | 15 CM | \$34.99



Adson Forcep, Tungsten Carbide, 12 CM

- 06-2033-12 | Serrated | \$88.99
- 06-2044-12 | 1x2 Teeth | \$88.99
- 06-2045-12 | 5x6 Teeth | \$88.99



18 CM Gerald Tissue, Tungsten Carbide

- 200-105 | Serrated | \$97.99
- 200-106 | 1x2 Teeth | \$97.99



Cushing Dressing Forcep, Tungsten Carbide, 18 CM

- 06-2037-18 | Serrated, Angled | \$97.99
- 200-107 | Serrated, Straight | \$97.99
- 06-2038-18 | Serrated, Bayonet | \$97.99



Flat Osteotome, 7-10-13-15-18 MM Markings

- 2303 | 4 MM, Straight | \$74.99
- 2304 | 4 MM, Curved | \$74.99
- 2305 | 7.5 MM, Straight | \$74.99
- 2306 | 7.5 MM, Curved | \$74.99



Round Osteotome, 7-8-13-15-18 MM Markings

- 2308 | #1, 2.7 MM | \$74.99
- 2309 | #2, 3.2 MM | \$74.99
- 2310 | #3, 3.7 MM | \$74.99
- 2311 | #4, 4.2 MM | \$74.99



1 Woodson **2268** **\$67.99**

The Woodson Periosteal Elevator is designed to reflect & retract the mucoperiosteum after incisions of the gingival tissue with 2.2 mm & 3.8 mm width

Buser **2275** **\$67.99**

This product is 3.8 mm wide on one side & 4.3 mm wide on other side

9 Molt **1633** **\$67.99**

This product is 4.3 mm on one side & 8.4 mm on other

16 Freer Strongly Curved **2348** **\$67.99**

5 Molt **1638** **\$67.99**

14 Goldman Fox **1640** **\$67.99**

23 Seldin **1641** **\$67.99**

A large flat handle instrument with one sharp blade & one dull blade for reflecting tissue. This product is 10.7 mm on one side & 13.2 mm on other

24G **1642** **\$67.99**

Periosteal with one round blade & other flat angled sharp blade. This product is 4.0mm on one side

3 Prichard **1643** **\$67.99**

Featuring two blades, one large flat blade & one narrow sharp blade on round tip. This product is 4.7mm wide on one side & 11mm on other side

9 Molt **1633S** **\$84.99**

3 Prichard **1643S** **\$84.99**

Allen Anterior **2276S** **\$84.99**

Allen **2277S** **\$84.99**

Buser	2275S	\$84.99
Regular		
152 Kramer Nevins	2272S	\$84.99
Regular		
1 Posterior	2337	\$64.99
Regular		
1 Posterior	PT1S	84.99
Siyah Series		
2 Anterior	1638	\$64.99
Regular		
2 Anterior	PT2S	84.99
Siyah Series		
3 Posterior Right	2339	\$64.99
Regular		
3 Posterior Right	PT3S	84.99
Siyah Series		
4 Posterior Left	2340	\$64.99
Regular		
4 Posterior Left	PT4S	\$84.99
Siyah Series		
5 Anterior Angled	2341	\$64.99
Regular		
5 Anterior Angled	PT5S	\$84.99
Siyah Series		
6 Anterior Straight	2342	\$64.99
Regular		
6 Anterior Straight	PT6S	\$84.99
Siyah Series		





Channel Retractor
21-11194-08 | 16CM | \$74.99



Bishop Cheek Retractor
1366 | 13 CM | \$29.99



University of Minnesota Cheek Retractor
1367 | 14 CM | \$29.99



Columbia Cheek Retractor
1369 | 14 CM | \$29.99



Weider Tongue Retractor
1374 | Small | \$74.99
1375 | Large | \$74.99



Middeldorpf Retractor
MIDLS | 17 x 14 mm | \$94.99



Surgical Retractor, Downward Curve
SR2 | 10 x 42 mm | \$94.99



Surgical Retractor, Upward Curved
SR8 | 80 x 16 mm | \$94.99



Ramus Retractor, V Notch, Curved Down
SR9 | 70 x 11 mm | \$94.99



Implant Retractor
IMPLR | 23 CM | \$154.99



1 Austin Retractor Left
200-103 | \$84.99



Molt Mouth Gag
1381 | Adult/Large | \$124.99
1382 | Child/Medium | \$124.99
1383 | Infant/Small | \$124.99

ROOT TIP PICKS

1 Heidbrink 1587 \$67.99



2 Heidbrink 1588 \$67.99



3 Heidbrink 1589 \$67.99



13/14 Heidbrink 2241 \$67.99



2/3 Howard 2250 \$67.99



Push Scalpel Handle AS9217 \$99.99



#3 1385 \$20.99



#5 Straight 1389 \$34.99



#5A Curved 1390 \$39.99



#5D Double Blade 1391 \$64.99



SCALPEL HANDLES



1X Miller Colburn	2282	\$79.99
The 1X Miller-Colburn cross-cut bone file is used for final smoothing & shaping of bone with 8.1 mm & 4.8 mm wide working ends		
2X Miller Colburn	2283	\$79.99
The 2X Miller-Colburn cross-cut bone file is used for final smoothing & shaping of bone with 6.2 mm wide working ends		
3X Miller Colburn	2284	\$79.99
The 3X Miller-Colburn cross-cut bone file is used for final smoothing & shaping of bone with 6.2 mm & 4.8 mm working ends		
21 Miller	2280	\$79.99
The 21 Miller straight-cut bone file is used for final smoothing & shaping of bone with 8.1 mm & 4.8 mm working ends		
67 Howard	1649	\$79.99
The 67 Miller straight-cut bone file used for final smoothing & shaping of bone with 4.8mm wide, angled working ends		
45 Miller	1650	\$79.99
The 45 Miller straight-cut bone file used for final smoothing & shaping of bone with 6.2mm wide working ends		
12 Howard	1652	\$79.99
The 12 Howard straight-cut bone file is used for final smoothing & shaping of bone with 8.1mm wide working ends		
64 Miller	1654	\$79.99
The 64 Miller straight-cut bone file used for final smoothing & shaping of bone with 6.2 mm & 4.8mm wide working ends		
1 Wahl	1655	\$79.99
The 1W Wahl straight-cut bone file is used for final smoothing & shaping of bone with 8.1mm wide working ends		
1 Gardner Single Bevel	2285	\$79.99
3mm		
2 Gardner Single Bevel	2286	\$79.99
4mm		
1/2 Siyah Buser Modified	2173S	\$84.99
4 mm/5 mm		
2/3 Siyah Buser Modified	2174S	\$84.99
5 mm/6 mm		
1/3 Siyah Buser Modified	2175S	\$84.99
4 mm/6 mm		

84 Lucas

The 84 Lucas Surgical curette is designed for curettage, cyst removal & tooth socket debridement with 2mm working ends.

- 1624** | Regular | **\$64.99**
1624S | Siyah | **\$84.99**
1624SS | Serrated | **\$94.99**

**85 Lucas**

The 85 Lucas Surgical curette is designed for curettage, cyst removal & tooth socket debridement with 2.5mm working ends.

- 1625** | Regular | **\$64.99**
1625S | Siyah | **\$84.99**
1625SS | Serrated | **\$94.99**

**86 Lucas**

The 86 Lucas Surgical curette is designed for curettage, cyst removal & tooth socket debridement with 2.8mm working ends.

- 1626** | Regular | **\$64.99**
1626S | Siyah | **\$84.99**
1626SS | Serrated | **\$94.99**

**87 Lucas**

The 87 Lucas Surgical curette is designed for curettage, cyst removal & tooth socket debridement with 3.5mm working ends.

- 1627** | Regular | **\$64.99**
1627S | Siyah | **\$84.99**
1627SS | Serrated | **\$94.99**

**88 Lucas**

The 88 Lucas Surgical curette is designed for curettage, cyst removal & tooth socket debridement with 4.7mm working ends.

- 1628** | Regular | **\$64.99**
1628S | Siyah | **\$84.99**
1628SS | Serrated | **\$94.99**

**2/4 Molt**

The 2/4 Molt surgical curette is designed for curettage, cyst removal & tooth socket debridement with 3.5mm & 7mm working ends.

- 2251** | Regular | **\$64.99**
2251S | Siyah | **\$84.99**





#1 Siyah Tunneling **TKN1** \$84.99



#2 Siyah Tunneling **TKN2** \$84.99



1/2 Siyah Orban **2168S** \$84.99



15/16 Siyah Kirkland **1241S** \$84.99



11 Siyah Goldman Fox **1246S** \$84.99



Siyah Grisdale Bone Bone Grafting Condenser **2336S** \$84.99

8-11-13-15-18, 4.0/3.35 mm



Siyah LaBanca Bone Grafting Plugger/Condenser **2331S** \$84.99

3 x 1.5 mm



Type A Aspirating Anesthetic Syringe **1399** \$64.99



Type CW Aspirating Anesthetic Syringe **1400** \$64.99



Wire Twister, TC Tips **13-3892-16** \$144.99

16 CM

TWEEZERS



College **1157** \$19.99



Meriam **1158** \$19.99



College, Locking **1160** \$24.99



Signature Mirror Handle Cone Socket **1168** \$21.99



Siyah Mirror Handle - Cone **1168b** \$29.99

MIRRORS & HANDLES



#4 Cone Socket Mirror Double Mirror
1171 | \$7.99



#5 Cone Socket Mirror Double Mirror
1172 | \$7.99



#4 Cone Socket Mirror
1169 | \$47.99/12pc



#5 Cone Socket Mirror
1170 | \$47.99/12pc



#5 Cone Socket Siyah
SMC5 | \$14.99



#5 Cone Socket Double Siyah
SMDC5 | \$14.99

SIYAH EXPLORERS

2 Explorer **1112s** \$39.99



3CH Cowhorn Explorer **1114s** \$39.99



5 Explorer **1115s** \$39.99




TU17/21 Wilkins/Tufts **1116s** \$39.99



23/CP12 Explorer **1150s** \$39.99





3CH Cowhorn Explorer	1114	\$29.99
5 (17/23) Explorer	1115	\$29.99
Combination of #17 end for subgingival calculus detection and #23 end for caries detection		
TU17/23 Explorer	1116	\$29.99
Combination of #TU17 end for subgingival calculus detection and #23 end for caries detection; the #TU17 end features a balanced tip design for ease and comfort		
ODU 11/12 Explorer	1122	\$29.99
Ideal for posterior calculus detection in deep pockets, especially on proximal surfaces. Designed to explore an entire dentition with a single instrument		
2A Pigtail	2210	\$29.99
2R/2L Explorer	2212	\$29.99
3ES Explorer	2213	\$29.99
UNC12 CPN	1132	\$29.99
1-2-3-4-5-6-7-8-9-10-11-12		
Williams CPN	2217	\$29.99
1-2-3-5-7-8-9-10		
CP12 CPN	2219	\$29.99
3-6-9-12		
UNC15 CPN	2221	\$29.99
1-2-3-4-5-6-7-8-9-10-11-12-13-14-15		
GF/Williams	1143	\$29.99
1-2-3-5-7-8-9-10 1-2-3-5-7-8-9-10 Flat		
CP12/Williams	1144	\$29.99
1-2-3-5-7-8-9-10 3-6-9-12		
UT3/4 Cattoni	1145	\$29.99
3-6-9-12 3-6-9-12		

3N Nabers 2222 \$29.99

1-2-3-5-7-8-9-10 | 1-2-3-5-7-8-9-10 Flat

Williams/Williams Offset 2410 \$29.99

1-2-3-5-7-8-9-10 | 1-2-3-5-7-8-9-10

1N Nabers 1140 \$29.99

2N Nabers 1141 \$29.99

Q2N Nabers 1142 \$29.99

3-6-9-12 | 3-6-9-12

UNC15/11.5B Ball 1137 \$29.99

3.5-5.5-8.5-11.5 Ball | 1-2-3-4-5-6-7-8-9-10-11-12-13-14-15

CP12/11.5B Ball 1138 \$29.99

3.5-5.5-8.5-11.5 Ball | 3-6-9-12

PQW/11.5B Ball 1139 \$29.99

3.5-5.5-8.5-11.5 Ball | 1-2-3-5-7-8-9-10

CP8 1124 \$29.99

3-6-8-11

CP12 1127 \$29.99

3-6-9-12

CP15UNC 1128 \$29.99


1-2-3-4-5-6-7-8-9-10-11-12-13-14-15


PQ Williams 1130 \$29.99

1-2-3-5-7-8-9-10

CP11.5B Ball 1131 \$29.99

3.5-5.5-8.5-11.5 Ball





23/QW Williams	1146	\$29.99
[1-2-3-5-7-8-9-10]		
23/CP8	1147	\$29.99
[3-6-8-11]		
23/CP11	1148	\$29.99
[3-6-8-11]		
23/CP11.5B Ball	1149	\$29.99
[3.5-5.5-8.5-11 Ball]		
23/CP12	1150	\$29.99
[3-6-9-12]		
23/UNC15	1151	\$29.99
[1-2-3-4-5-6-7-8-9- 10-11-12-13-14-15]		
TU17/UNC15	1152	\$29.99
[1-2-3-4-5-6-7-8-9- 10-11-12-13-14-15]		
6/QW Williams	1153	\$29.99
[1-2-3-5-7-8-9-10]		
3A/UNC15	1154	\$29.99
[1-2-3-4-5-6-7-8-9- 10-11-12-13-14-15]		
23/Goldman Fox	1155	\$29.99
[1-2-3-5-7-8-9-10 Flat]		
23/O Michigan	1156	\$29.99
[3-6-8]		
23/UNC12	2223	\$29.99
[1-2-3-4-5-6-7-8-9-10-11-12]		
CPNT3/CP12	2224	\$29.99
[3-6-9-12]		

H6/H7 Hygienist **1175** **\$29.99**

Super-thin, contraangle design is for interproximal areas. Shank angulation provides accessibility to all anterior and premolar surfaces

H5/33 Hygienist Jacquette **1176** **\$29.99**

Super-thin sickle scaler combined with a straight sickle blade. For scaling interproximal and cervical areas of anterior teeth

6/7 Goldman **2102** **\$29.99**

Shortened version of the H6/H7 Scaler. Shank angulation provides accessibility to all anterior and premolar surfaces

137 Anterior **2111** **\$29.99**

Small anterior sickle paired with an anterior universal curette

U15/30 Jacquette **1177** **\$29.99**

Designed for the removal of subgingival calculus

U15/33 Towner Jacquette **1178** **\$29.99**

Designed for the removal of subgingival calculus

204 Sickle **1183** **\$29.99**

Designed for the removal of subgingival calculus

204IUFW Scaler **2105** **\$29.99**

Designed for the removal of subgingival calculus

204S Scaler **1184** **\$29.99**

Posterior sickle scaler designed for the removal of supragingival calculus

204SD Scaler **1185** **\$29.99**

Designed with an angle that allows greater access to distal surfaces

135 Posterior **1187** **\$29.99**

Designed for the removal of subgingival calculus

REM 3/4 Remington **2109** **\$29.99**

Designed for the removal of subgingival calculus

30/33 Jacquette **1180** **\$29.99**

Combination of standard and miniature anterior Jacquette

31/32 Jacquette **1181** **\$29.99**

Posterior standard jacquette

34/35 Jacquette Scaler **1182** **\$29.99**

Posterior miniature Jacquette





H3/H4 Jacquette **2103** \$29.99

14/15 Posterior **2104** \$29.99

McCall 11/12 **1191** \$29.99

McCall 13S/14S Pointed Curette **1192** \$29.99

McCall 17S/18S Curette **1195** \$29.99

McCall 19/20 Curette **1196** \$29.99

McCall 11A/12A **2126** \$29.99

McCall 11/12S Modified **2127** \$29.99

Goldman Fox 1 Universal **1197** \$29.99

Ideal for removing heavy supragingival and subgingival calculus, and on proximal surfaces and interproximal areas during periodontal surgical procedures

Goldman Fox 2 Curette **1198** \$29.99

Contour design for removal of calculus on anterior teeth and premolars

Goldman Fox 3 Curette **1199** \$29.99

Has the same blade surface as GF-2, which is specifically designed to work on the proximal surface of the posterior teeth

Goldman Fox 4 Curette **1200** \$29.99

With an additional bend in the shank for easier calculus removal in the molar region

Goldman Fox 5 Hoe Buccal/Lingual **1201** \$29.99

For removing heavy calculus deposits

Goldman Fox 6 Hoe **1202** \$29.99

For mesial/distal surfaces of the posterior teeth or facial/lingual surfaces of the anterior teeth

GOLDMAN FOX SET

Goldman Fox 21 **1203** **\$29.99**

Sharp spoon on one end is for debridement in curettage procedures. Small scaler on the other end can be used with either a circumferential or vertical stroke



Goldman Fox 21S/33 Combination **2128** **\$29.99**

Small end Goldman-Fox scaler with a #33 straight miniature Jacquette scaler



1/2 Langer Curette **1204** **\$29.99**

For mandibular posterior



3/4 Langer Curette **1205** **\$29.99**



5/6 Langer Curette **1206** **\$29.99**



17/18 Langer Curette **1207** **\$29.99**



H5-L5 Langer **H5L5** **\$29.99**

2-in-1 anterior instrument with an H5 straight sickle scaler on one end for deposit removal and interproximal access



H5-L5 Mini Langer **H5L5M** **\$29.99**

Features a mini Langer version of the H5-L5



Nebraska 128/Langer 5 **N128L5** **\$29.99**

2 in 1 anterior instrument with a Nebraska 128 on one end for narrow, interproximal access



Nebraska 128/Langer Mini 5 **N128L5M** **\$29.99**



2R/2L Columbia **1215** **\$29.99**

For removing supragingival and subgingival calculus; useful on proximal surfaces and interproximal areas during periodontal surgical procedures.



4R/4L Columbia **1216** **\$29.99**

Posterior application



13/14 Columbia **1217** **\$29.99**

Universal application



LANGER UNIVERSAL CURETTES

COLUMBIA UNIVERSAL CURETTES



1/2 Barnhart **1213** **\$29.99**

5/6 Barnhart **1214** **\$29.99**

13/14 Indiana University **1208** **\$29.99**
Anterior

17/18 Indiana Univeristy **1209** **\$29.99**
Posterior

1/2 Bunting **1210** **\$29.99**
Curette Chisel Combination

5/6 Bunting **1211** **\$29.99**

7/8 Younger Good **1212** **\$29.99**

3/4 Rule **1218** **\$29.99**

TAHA 1 Anterior **TAHA1** **\$34.99**

An elongated disc end for anterior lingual stain and calculus removal. Paired with a super thin, anterior curved sickle.

TAHA 2 Posterior **TAHA2** **\$34.99**

A super thin posterior curved sickle scaler with a contra-angle design for better interproximal reach and ergonomic hand positioning. Ideal for light to moderate interproximal calculus removal.

TAHA 3 Posterior **TAHA3** **\$34.99**

This ultra thin, dual cutting blade proves effective for reaching tight, narrow contacts and inter-proximal areas. Ideal for effective removal of light to moderate calculus in tight, narrow posterior inter-proximal.

TAHA 4 Posterior **TAHA4** **\$34.99**

Tapered dual cutting blades provide precise access for effective deposit removal and increased shank rigidity provides for effective removal of burnished heavy and tenacious deposits. Ideal for scaling heavy posterior inter-proximal deposits. Largest blade in the TAHA collection.

TAHA 1/H5 Anterior **TAHA1H5** **\$34.99**

An elongated disc end for anterior lingual stain and calculus removal paired with a larger H5 sickle - a thin curved blade for interproximal and cervical anterior areas.

Super Taha Scaler **STAHA** **\$34.99**

Equivalent to R138 (Montana Jack) scaler.

1/2 Gracey **1229** \$29.99
 1/2 Rigid Gracey **G12R** \$29.99

Slight contra-angle for anterior incisors and canines



3/4 Gracey **1230** \$29.99
 1/2 Rigid Gracey **G34R** \$29.99

Short contra-angle for anterior incisors and canines



5/6 Gracey **1231** \$29.99
 5/6 Rigid Gracey **G56R** \$29.99

Medium contra-angle for anteriors and premolars



7/8 Gracey **1232** \$29.99
 7/8 Rigid Gracey **G78R** \$29.99

Medium contra-angle for premolars and molar facial and lingual surfaces



9/10 Gracey **1233** \$29.99
 9/10 Rigid Gracey **G910R** \$29.99

Long contra-angle for molar buccal and lingual surfaces and a more pronounced angle than 7/8 for less accessible root surfaces



11/12 Gracey **1234** \$29.99
 11/12 Rigid Gracey **G1112R** \$29.99

Angulated to reach mesial surfaces of posterior teeth



13/14 Gracey **1235** \$29.99
 13/14 Rigid Gracey **G1314R** \$29.99

Angulated to reach distal surfaces of posterior teeth



15/16 Gracey **1236** \$29.99
 15/16 Rigid Gracey **G1516R** \$29.99

Same shank bends as 13/14, but blade positioned to reach mesial posterior surfaces



17/18 Gracey **1237** \$29.99
 17/18 Rigid Gracey **G1718R** \$29.99

Accentuated angles improve access to distal posterior surfaces. The 17/18 also has smaller blade and slightly longer terminal shank



11/14 Gracey **2143** \$29.99
 11/14 Rigid Gracey **G1114R** \$29.99



12/13 Gracey **2144** \$29.99
 12/13 Rigid Gracey **G1213R** \$29.99





1/2 Mini Gracey **1220** \$29.99
 1/2 Mini Rigid Gracey **G12MR** \$29.99

Mini version of 1/2 Gracey



3/4 Mini Gracey **1221** \$29.99
 3/4 Mini Rigid Gracey **G34MR** \$29.99

Mini Version of 3/4 Gracey



5/6 Mini Gracey **1222** \$29.99
 5/6 Mini Rigid Gracey **G56MR** \$29.99

Mini Version of 5/6 Gracey



7/8 Mini Gracey **1223** \$29.99
 7/8 Mini Rigid Gracey **G78MR** \$29.99

Mini Version of 7/8 Gracey



9/10 Mini Gracey **1224** \$29.99
 9/10 Mini Rigid Gracey **G910MR** \$29.99

Mini Version of 9/10 Gracey



11/12 Mini Gracey **1225** \$29.99
 11/12 Mini Rigid Gracey **G1112R** \$29.99

Mini Version of 11/12 Gracey



13/14 Mini Gracey **1226** \$29.99
 13/14 Mini Rigid Gracey **G1314MR** \$29.99

Mini version of 13/14 Gracey



15/16 Mini Gracey **1227** \$29.99
 15/16 Mini Rigid Gracey **G1516MR** \$29.99

Mini version of 15/16 Gracey



17/18 Mini Gracey **1228** \$29.99
 17/18 Mini Rigid Gracey **G1718MR** \$29.99

Mini Version of 17/18 Gracey



11/14 Mini Gracey **G114M** \$29.99
 11/14 Mini Rigid Gracey **G1114MR** \$29.99



12/13 Mini Gracey **G1213M** \$29.99
 12/13 Mini Rigid Gracey **G1213MR** \$29.99

ENDO EXPLORERS

DG 16 1119 \$29.99



DG 16/17 1120 \$29.99



DG 16/23 1121 \$29.99



ENDO EXCAVATORS

31L Spoon 2406 \$29.99



31W Oval Spoon 2407 \$29.99



32L Spoon 2408 \$29.99



33L Spoon 2409 \$29.99



2 Glick 2405 \$29.99



ANTERIOR PLUGGERS

8A 29mm . 0.40mm 1270 \$29.99

9A 29mm . 0.55mm 1271 \$29.99

10A 29mm . 0.75mm 1272 \$29.99

11A 29mm . 1.05mm 1273 \$29.99

8.5A 29mm . 0.45mm 1274 \$29.99

9.5A 29mm . 0.65mm 1275 \$29.99

10.5A 29mm . 0.90mm 1276 \$29.99

11.5A 29mm . 1.15mm 1277 \$29.99

12A 29mm . 1.30mm 1278 \$29.99





8 22mm	• 0.40mm	1279	\$29.99
9 22mm	• 0.55mm	1280	\$29.99
10 22mm	• 0.75mm	1281	\$29.99
11 22mm	• 1.05mm	1282	\$29.99
8.5 22mm	• 0.45mm	1283	\$29.99
9.5 22mm	• 0.65mm	1284	\$29.99
10.5 22mm	• 0.90mm	1285	\$29.99
11.5 22mm	• 1.15mm	1286	\$29.99
12 22mm	• 1.30mm	1287	\$29.99

POSTERIOR PLUGGERS

ROOT FRAGMENT FORCEPS

ENDO ACCESSORIES



Stieglitz 1295 \$48.99



Stieglitz Curved 1296 \$53.99



Stieglitz Angled 1297 \$55.99



Peet Curved 1298 \$39.99



Stainless Steel Ruler ERULE6 \$29.99
6 Inches



Rubber Dam Clamp Organizing Board RDCOB8 \$55.99

Store, sterilize and organize with eight pegs that securely hold your clamps.
Dimensions: 4.62" x 2.4" x .375"
(117 mm x 61 mm x 10 mm). 8 Pieces.

ISO SIZED PLUGGERS

ISO Size 30 (24MM/0.3MM) **1262** \$29.99

ISO Size 40 (24MM/0.4MM) **1263** \$29.99

ISO Size 50 (24MM/0.5MM) **1264** \$29.99

ISO Size 60 (24mm) **1265** \$29.99

1/3 Root Canal Plugger (21mm) 0.4MM/0.45MM **2257** \$29.99

5/7 Root Canal Plugger (21mm) 0.5MM/0.75MM **2351** \$29.99

9/11 Root Canal Plugger (21mm) 1.0MM/1.15MM **2352** \$29.99

L1 18mm (0.45 MM) **1266** \$29.99

L2 18mm (0.5 MM) **1267** \$29.99

L3 18mm (0.75 MM) **1268** \$29.99

L4 18mm (1.00 MM) **1269** \$29.99

LUKS PLUGGERS





0
Flat jaws and high bow for small lower canines and first premolars
1307 | \$12.99



6
Slightly offset jaws for labial caries on anterior teeth
1317 | \$12.99



9
Offset flat jaws and stiff bow for anterior teeth
1322 | \$12.99



9S
Offset flat serrated jaws for anterior teeth
1323 | \$12.99



212
Flat jaws for labial caries on anterior teeth
1345 | \$12.99



212SA
Slightly offset staggered jaws for odd shaped anteriors
1346 | \$12.99



3
Flat narrow jaws for small lower molars
1314 | \$12.99



7
Flat jaws for lower molars
1318 | \$12.99



7A
Slightly reversed festooned jaws for large lower molars
1319 | \$12.99



12A
Serrated offset jaws for lower right molars; paired with RDCM13A
1326 | \$12.99



13A
Serrated offset jaws for lower left molars; paired with RCM12A
1327 | \$12.99



26
Slightly reversed festooned jaws for small lower molars
1333 | \$12.99



28
Flat jaws for lower molars
1337 | \$12.99



W3
Flat narrow jaws for small lower molars
1348 | \$12.99



W7
Flat jaws for lower molars
1350 | \$12.99



0
Flat jaws and high bow for longer premolars
1308 | \$12.99



1
Slightly festooned jaws for upper premolars
1309 | \$12.99



1A
Deeply festooned jaws for canines and premolars
1310 | \$12.99



2
Flat jaws for premolars
1311 | \$12.99



2A
Flat jaws for larger premolars
1312 | \$12.99



2AS
Flat serrated jaws for larger premolars
1313 | \$12.99



27N
Slightly festooned jaws for premolars
1336 | \$12.99



29
Flat jaws for lower premolars
1338 | \$12.99



W2
Flat jaws for premolars
1347 | \$12.99



209
Small, flat jaws with stiff bow for lower premolars
2233 | \$12.99



8A
Deeply festooned four-point jaws for molars
1321 | \$12.99



10
Flat jaws designed for small upper left and lower right molars; paired with RDCM11
1324 | \$12.99



11
Flat jaws designed for small upper right and lower left molars; paired with RDCM10
1325 | \$12.99



14
Deeply festooned jaws for partially erupted or irregularly shaped molars
1328 | \$12.99



14A
Deeply festooned jaws for partially erupted or irregularly shaped molars; larger than RDCM14
1329 | \$12.99



18
Flat jaws for upper right and lower left molars; paired with RDCM25
1330 | \$12.99



24
Flat jaws for upper left and lower right molars; paired RDCM24
1331 | \$12.99



25
Flat jaws for large molars
1332 | \$12.99



27
Slightly festooned jaws for large molars
1335 | \$12.99



56
Slightly festooned serrated jaws for large partially erupted or irregularly shaped molars
1341 | \$12.99



56S
Slightly festooned jaws for molars
1342 | \$12.99



W8A
Deeply festooned jaws for partially erupted or irregularly shaped molars
1351 | \$12.99



W14A
Deeply festooned jaws and extended bow for distal access to molars
2231 | \$12.99



8AD
Deeply festooned jaws and extended bow for distal access to molars
2232 | \$12.99




4
Flat narrow jaws for small upper molars
1315 | \$12.99



5
Slightly festooned jaws for large upper molars
1316 | \$12.99



8
Deeply festooned jaws for upper molars
1320 | \$12.99



26N
Flat jaws for upper molars
1334 | \$12.99



30
Offset jaws for upper right molars; paired with RDCM31
1339 | \$12.99



31
Offset jaws for upper left molars; paired with RDCM30
1340 | \$12.99



201
Slightly festooned jaws for upper molars
1343 | \$12.99



205
Slightly festooned jaws for large upper molars
1344 | \$12.99



W4
Flat narrow jaws for small upper molars
1349 | \$12.99



Ivory Lightweight
17cm
1299 | \$65.99

Brewer
17cm
1300 | \$75.99



U of Washington
17cm
1301 | \$75.99



Ivory
1302 | \$95.99



Ainsworth
1303 | \$95.99



Small Child
3.5 inch/9cm
1304 | \$29.99

Child
5 inch/13cm
1305 | \$29.99

Adult
6 inch/15cm
1306 | \$29.99

1 Goldstein 2467 \$39.99

Small universal style with rounded plugger tip and a narrow paddle for initial placement and contouring of Class I, II, and III restorations

**2 Goldstein** 2468 \$39.99

Larger universal style for final placement and contouring of Class I, II, and III restorations.

**3 Goldstein Extra Flex** 2469 \$39.99

Flexible, reversed, flared paddle design for shaping and placement of Class III and IV restorations

**4 Goldstein Extra Flex** 2470 \$39.99

Flexible, paired, offset, paddle-shaped blades for placing and shaping material on posterior, mesial and distal surfaces. Reverse angle is also useful for placing and shaping anterior bonded restorations

**5 Goldstein Flexi Thin** 300-101 \$39.99

Small reverse angle tips make it easy to place fissures, grooves and pits creating the ideal occlusal anatomy in hard to reach posterior areas.

**6 Goldstein Flexi Thin** 300-102 \$39.99

Large reverse angle tips make it easy to place fissures, grooves and pits creating the ideal occlusal anatomy in hard to reach posterior areas

**1 Mini Goldstein** 2473 \$39.99

Mini version of 1 Goldstein for small pits and fissures, tunnel preparations or minor tooth defects on lower anteriors.

**3 Mini Goldstein Extra Flex** 2475 \$39.99

Mini version of 2 Goldstein. Can also be used for packing gingival retraction cord

**4 Mini Goldstein Extra Flex** 2476 \$39.99

Mini version of 4 Goldstein for placing and shaping material in difficult to access mesial and distal posterior restorations





- | | | |
|--|-------------|----------------|
| AB1 Bhogolian | 2363 | \$39.99 |
| Unique combination of thin, knife-shaped blade with standard angled blade. Knife blade allows controlled, efficient manipulation of composite even in gingival areas. Application: Class III, IV, V | | |
| AB2 Bhogolian | 2364 | \$39.99 |
| Used for measuring composite layers and shaping occlusal anatomy | | |
| W3 | 2365 | \$39.99 |
| Combination of medium-sized blade with small condenser tip for universal adaptability. Ideal for placement, layering and general contouring. Application: Class I, II, III, IV, V | | |
| 8A | 2366 | \$39.99 |
| Use for packing gingival retraction cord, as well as to place and contour facial aspects | | |
| 4/5 Gregg | 2367 | \$39.99 |
| Off-angled blades allow easy adaptability to mesial and distal surfaces of posterior teeth, providing increased interproximal access and better visibility of the working area Application: Class II, V | | |
| Interproximal Carver | 2368 | \$39.99 |
| Extremely thin flexible blades are opposed for easy handling of composite materials and interproximal contouring. Application: Class III, IV, V | | |
| Interproximal Carver, Long | 2369 | \$39.99 |
| Used for placement of the composite increments against the cavity wall or adjacent tooth surface | | |
| A6 | 2370 | \$39.99 |
| Large, thin blades are opposed for adaptability to any situation, including veneers, where broad contouring or carving strokes are needed. Application: Class II, III, IV and V Mini version of 4 Goldstein for placing and shaping material in difficult to access mesial and distal posterior restorations | | |
| Contouring Instrument | 2371 | \$39.99 |
| Used for shaping of inclines, planes or developmental lobes for anterior and posterior restorations. The instrument has different angles of curvature on each end that provides a buccal and lingual orientation for posterior shaping or a facial and lingual orientation for anterior shaping | | |
| 1 Loesche | 1936 | \$39.99 |
| 2 Loesche | 1937 | \$39.99 |
| 11 Anterior | 1938 | \$39.99 |
| 179 Anterior | 1939 | \$39.99 |

21B Acorn **300-103** **\$39.99**

Acorn-shaped instrument for forming occlusal anatomy in posterior restorations



27/29 **300-104** **\$39.99**

Used to blend material for final contouring, to achieve sculpting of areas like grooves, fissures or pits. Can also be used to form occlusal anatomy.



2 Ladmore **300-105** **\$39.99**

Medium to large rounded tips for condensing composite materials



3 Ladmore **300-106** **\$39.99**

Small to medium slightly rounded tips for condensing composite materials



Ball 18 **300-107** **\$39.99**

Used to smooth and shape composite



Small/Medium Ball **300-132** **\$39.99**

Used to direct and form the composite increments against the cavity wall. The shape conforms to the rounded cavity surfaces and allows ease of access into the rounded corners or junctions of the cavity surfaces to condense and shape the composite against the cavity wall.



Freedman Duckhead **300-108** **\$44.99**

Used to contour the convexity of the cusp ridge, developing the anatomy in a single motion



Freedman Small Contact Forming **300-109** **\$44.99**

Acorn-shaped instrument for forming occlusal anatomy in posterior restorations



Freedman Large Contact Forming **300-110** **\$44.99**

Oval-shaped paired instrument designed to provide improved contact forming for large Class II Restorations



Small/Medium Contact Forming **300-133** **\$44.99**

Rounded cone-shaped paired instrument designed to provide improved contact forming for small/medium Class II restorations



Medium Large Contact Forming **300-134** **\$44.99**

Rounded cone-shaped instrument to provide improved contact forming for medium/ large Class II restorations





A Cosmetic Contouring 300-111 \$44.99

Identical, opposing, large, flexible, ovalshaped blades, straight and angled, for contouring composite material on larger facial surfaces of central incisors.



B Cosmetic Contouring 300-112 \$44.99

Identical, opposing, spear-shaped blades, straight and angled, used for contouring composite material on smaller facial surfaces of central incisors.



C Cosmetic Contouring 300-113 \$44.99

Flexible, oval-shaped blades – one slightly larger – for interproximal contouring on central incisors.



D Cosmetic Contouring 300-114 \$44.99

Used when working near or at interproximal areas. Straight end compacts composite material, while sharp knife edge cuts composite to avoid bonding to adjacent tooth.



E Cosmetic Contouring 300-115 \$44.99

Small and medium curved blades for thinning and shaping composite material at the gingival areas.



F Cosmetic Contouring 300-116 \$44.99

Uniquely-shaped blades with curved and rounded tips for adding and shaping composite material on desired areas of facial incisors.



G Marginal Ridge & Embrasure Shaping Instrument 300-117 \$44.99

Allows formation of marginal ridges along with buccal and lingual embrasures while composite is uncured.



H Occlusal Anatomy Instrument 300-118 \$44.99

Designed to help attain proper occlusal form, function and improve marginal seal.



I Composite Packing Instrument 300-119 \$44.99

Aids in forming a properly filled axial box and occlusal portion.

1 Goldstein 1834 \$29.99

Small universal style with rounded plugger tip and a narrow paddle for initial placement and contouring of Class I, II, and III restorations.

**2 Goldstein 1835 \$29.99**

Larger universal style for final placement and contouring of Class I, II, and III restorations

**3 Goldstein Extra Flex 1836 \$29.99**

Flexible, reversed, flared paddle design for shaping and placement of Class III and IV restorations.

**4 Goldstein Extra Flex 1837 \$29.99**

Flexible, paired, offset, paddle-shaped blades for placing and shaping material on posterior, mesial and distal surfaces.

**1 Mini Goldstein 1838 \$29.99**

Mini version of 1 Goldstein for small pits and fissures, tunnel preparations or minor tooth defects on lower anteriors.


**2 Mini Goldstein 1839 \$29.99****3 Mini Goldstein Extra Flex 1840 \$29.99**

Mini version of 2 Goldstein. Can also be used for packing gingival retraction cord.

**4 Mini Goldstein Extra Flex 1841 \$29.99**

Mini version of 4 Goldstein for placing and shaping material in difficult to access mesial and distal posterior restorations.

**1 Woodson 1940 \$29.99****2 Woodson 1934 \$29.99****3 Woodson 1935 \$29.99**



13/14 [20-9-14]	1855	\$29.99
15/16 [15-8-14]	1856	\$29.99
17/18 [10-6-14]	1857	\$29.99
14/14 Off Angle [15-8-14]	1860	\$29.99
15/15 Off Angle [15-8-14]	1861	\$29.99
26 Distal [13-95-8-14]	1788	\$29.99
27 Mesial [13-80-8-14]	1789	\$29.99
28 Distal [10-95-7-14]	1790	\$29.99
29 Mesial [10-80-7-14]	1791	\$29.99
3/4 Wedelstaedt [11.5-15-3]	1794	\$29.99
5/6 Wedelstaedt [15-15-3]	1795	\$29.99
1/2 Wedelstaedt [20-15-3]	1796	\$29.99

Calcium Hydroxide Placer **2492** **\$29.99**

Calcium hydroxide or glass ionomer base/liner placement instrument. Also useful as a small burnisher.

BLADE EXCAVATORS

36/37 **1819** **\$29.99**

[15-8-14]



40/41 **1820** **\$29.99**

[11-7-14]



44/45 **1821** **\$29.99**

[20-9-14]



65/66 **1824** **\$29.99**

[10-6-12]



125/126 **1825** **\$29.99**

[2.5mm]



127/128 **1826** **\$29.99**

[2.0mm]



129/130 **1827** **\$29.99**

[1.7mm]



131/132 **1828** **\$29.99**

[1.4mm]



243 **1830** **\$29.99**

[1.0mm]



244 **1831** **\$29.99**

[1.4mm]



245 **1832** **\$29.99**

[2.0mm]



ENGLISH PATTERN EXCAVATORS

17W **1816** **\$29.99**

[1.2mm]



18W **1817** **\$29.99**

[1.5mm]




19W **1818** **\$29.99**

[1.8mm]



OVAL EXCAVATORS



E1	1807	\$29.99
[12-9-15]		
E2	1808	\$29.99
[15-9-15]		
E3	1809	\$29.99
[25-9-15]		
6	1810	\$29.99
[0.85mm]		
14	1811	\$29.99
[1.4mm]		
17	1812	\$29.99
[1.2mm]		
18	1813	\$29.99
[1.5mm]		
19	1814	\$29.99
[2.0mm]		
38/39	1815	\$29.99
[1.15mm]		



Non Slip	1871	\$29.99
Stainless steel, with removable non-slip base ring. Designed for easy amalgam handling.		



Miller Articulating Paper Forcep	1786	\$27.99
Straight		

SPOON EXCAVATORS

AMALGAM WELL

ARTICULATING FORCEP

1/2 Andrews **1874** \$29.99

[1.4mm], [2.0mm]



H1 Hollenback **1875** \$29.99

[1.0mm], [1.65mm]



H2 Hollenback **1876** \$29.99

[2.0mm], [2.7mm]



H3 Hollenback **1877** \$29.99

[2.5mm], [2.5mm]



H4 Hollenback **1878** \$29.99

[1.8mm], [2.4mm]



0/1 Marquette **1879** \$29.99

[1.0mm], [1.4mm]



2 Mortonson **1880** \$29.99

[0.6mm], [1.2mm]



1 Oregon **1887** \$29.99

[1.9mm], [1.2mm]



2 Oregon **1888** \$29.99

[1.0mm], [1.6mm]



3 Oregon **1889** \$29.99

[1.3mm], [2.2mm]



4 Oregon **1890** \$29.99

[2.0mm], [3.0mm]



1 Ward **1891** \$29.99

[1.2mm], [2.5mm]



2 Ward **1892** \$29.99

[2.2mm], [3.0mm]



3 Ward **1893** \$29.99

[1.5mm], [1.9mm]





9/10 Back Action **1899** **\$29.99**
[1.2mm], [1.6mm]



11/12 Back Action **1900** **\$29.99**
[1.5mm], [2.0mm]



153 Serrated **1901** **\$29.99**
[1.5mm], [2.5mm]



154 Serrated **1902** **\$29.99**
[1.4mm], [2.4mm]



1/2 Black **2373** **\$29.99**
[1.5mm], [2.5mm]



1/2 Serrated **2374** **\$29.99**
[1.5mm], [2.5mm]



0/1 Marquette Serrated **2375** **\$29.99**
[1.0mm], [1.4mm]



2 Serrated Smith **2376** **\$29.99**
[1.5mm], [1.5mm]



SPO Serrated 76MM **2377** **\$29.99**
[1.4mm], [2.0mm]



1T Tanner **1896** **\$29.99**
[0.6mm], [0.76mm]



2T Tanner **1897** **\$29.99**



3T Tanner **1898** **\$29.99**

Interproximal **1903** \$29.99

Extremely thin, flexible blade; ideal for interproximal contouring.

1/2 Hollenback **1904** \$29.99

Universal adaptability. Ideal for placing, carving and contouring amalgam.

3S Hollenback Carver **1905** \$29.99

Design characteristics similar to 1/2 Hollenback but with slightly larger blades.

3 Hollenback Carver **1906** \$29.99

Design characteristics similar to 1/2 Hollenback but with significantly larger blades.

8 Wiland **1907** \$29.99

Extremely thin curved blade; ideal for adapting to interproximal surfaces.

Interproximal Off Angle **2378** \$29.99

Extremely thin, flexible blade. Off-angle provides better access to posterior areas.

6 Tanner **2382** \$29.99

1 Hollenback **2395** \$29.99

2 Hollenback **2396** \$29.99

L7 Levy Carver **2397** \$29.99

14L Carver **CS10** \$39.99

Also known as Gold Knife Carver


1/2 Discoid Cleoid **1908** \$29.99

4/5 Discoid Cleoid **1909** \$29.99

3/6 Discoid Cleoid **1910** \$29.99

89/92 Discoid Cleoid **1912** \$29.99





21 Acorn **1917** **\$29.99**

Extremely thin, flexible blade; ideal for interproximal contouring.

21B Acorn **2314** **\$29.99**

Extremely thin, flexible blade; ideal for interproximal contouring.

18 **1918** **\$29.99**

Extremely thin, flexible blade; ideal for interproximal contouring.

26/27S **1919** **\$29.99**

Extremely thin, flexible blade; ideal for interproximal contouring.

27/29 **1920** **\$29.99**

Extremely thin, flexible blade; ideal for interproximal contouring.

28/29 **1921** **\$29.99**

Extremely thin, flexible blade; ideal for interproximal contouring.

2 Beavertail **1914** **\$29.99**

Extremely thin, flexible blade; ideal for interproximal contouring.

2 Double Ended Beavertail **1915** **\$29.99**

Extremely thin, flexible blade; ideal for interproximal contouring.

2/29 Beavertail **1916** **\$29.99**

Extremely thin, flexible blade; ideal for interproximal contouring.

2 Ladmore **1922** **\$29.99**

Extremely thin, flexible blade; ideal for interproximal contouring.

3 Ladmore **1923** **\$29.99**

Extremely thin, flexible blade; ideal for interproximal contouring.

CHRISTENSEN CROWN REMOVERS



Right Angled 2604 \$74.99
Posterior



Straight 2603 \$74.99
Anterior



45° Angle 300-124 \$74.99
For cuspids, bicuspid and even first molars.



Straight Interproximal 300-123 \$74.99
For anterior crown removal.



Right Angle 300-126 \$74.99
For molars.



Occlusal 300-125 \$74.99
For occlusal separation especially in hard to remove crowns that have been bonded to the tooth.

TEMPORARY CROWN REMOVERS



Temporary Crown Remover 1805 \$124.99



Crown Removing Pliers, 14cm 1806 \$74.99



Trial Lower Crown Remover 2605 \$109.99



Trial Upper Crown Remover 2606 \$109.99

RESTORATIVE



24 Flexible (44mm) **1924** **\$29.99**

Flexible blade for mixing medium body cements.



324 Rigid (51mm) **1925** **\$29.99**

Rigid blade for mixing heavier or medium body cements.



A6 Rigid (25mm) **1926** **\$29.99**



7 Tapered **1927** **\$29.99**



8R Rigid **300-130** **\$29.99**



5A Knife **300-131** **\$29.99**



Boley Gauge

1785 | **\$44.99**



1 Iwanson Spring Caliper

For metal and porcelain

1783 | **\$29.99**



2 Iwanson Spring Caliper

For wax

1784 | **\$29.99**

PKT 1 **1928** **\$29.99**

Curved tapered tips used to flow on molten wax.

**PKT 2** **1929** **\$29.99**

Curved tapered tips used to flow on molten wax.

**PKT 3** **1930** **\$29.99**

Pointed burnisher used to perfect and enhance the supplemental and developmental grooves.

**PKT 3R** **1931** **\$29.99**

Similar to PKT 3, but with a rounded tip vs. a pointed one.

**PKT 4** **1932** **\$29.99**

Modified carver used to perfect the external contours and remove excess wax at the cavosurface margins.

**PKT 5** **1933** **\$29.99**

Special carver used to remove excess wax as cusp ridges are developed; its contour maintains desired convexity at these ridges.





The Tear Drop
400-101 | \$149.99



The Vertical
400-102 | \$149.99



The Horizontal
400-103 | \$149.99



The Hole Punch
400-104 | \$149.99



The Spot
678-804 | \$149.99



The Petite Punch
678-805 | \$149.99



The Counter
Clockwise Wedge
678-807 | \$149.99



The Clockwise Wedge
678-808 | \$149.99



The Eraser
678-811 | \$149.99



Universal Cut & Hold
Distal End Cutter
400-105 | \$149.99



Universal Cut & Hold
Distal End Cutter, Long
Handle
400-105L | \$149.99



Flush Cut Distal End
Cutter, No Hold
400-106 | \$149.99



Multi Use Pin &
Ligature Cutter
400-107 | \$149.99



Hard Wire Cutter, 15°
400-108 | \$149.99



Mini Pin & Ligature
Cutter
400-109 | \$149.99



Hard Wire Cutter
400-110 | \$149.99



Micro-Mini Pin &
Ligature Cutter
400-111 | \$149.99



Micro-Mini Pin &
Ligature Cutter, Long
Handle
400-111L | \$149.99



Micro-Mini Pin &
Ligature Cutter, 15°
400-113 | \$149.99



Slim Micro Cutter
400-114 | \$149.99



Micro Cutter
400-115 | \$149.99



Lingual Distal End
Cutter
400-116 | \$149.99



Lingual Pin & Lig
Cutter
400-118 | \$149.99



Flush Cut & Hold Distal
End Cutter
678-111 | \$149.99



3 Silicone Replacement
Rings for 678-113
678-113 | \$149.99



Weingart Pliers
400-119 | \$149.99



Weingart Pliers, Long Handle
400-119L | \$134.99



Slim Weingart Pliers
400-120 | \$149.99



How Pliers
400-121 | \$134.99



How Pliers, Offset
400-122 | \$134.99



Band Seating Pliers
400-123 | \$134.99



Adhesive Removing Pliers
400-124 | \$149.99



Posterior Band Removing Pliers, Long
400-125 | \$149.99



Posterior Band Removing Pliers, Short
400-126 | \$134.99



Tip Back Pliers
400-129 | \$134.99



Bracket Removing Pliers
400-130 | \$149.99



Angulated Bracket Removing Pliers
400-131L | \$134.99



Crown & Band Contouring Pliers
400-132 | \$134.99




Slim Crown & Band Contouring Pliers
400-132M | \$134.99



Slim Crown & Band Contouring Pliers, Modified
400-132 MC | \$134.99



Band Crimping Pliers
400-134 | \$134.99



Band Crimping Pliers with Spring
400-134M | \$134.99



Ultra Slim Weingart Pliers
400-135 | \$134.99



Slim Band Removing Pliers
400-136 | \$149.99



Lingual Weingart Pliers
400-137 | \$149.99
400-138 | \$149.99

90 Degree Bend



Lingual Bracket Removing Pliers
400-139 | \$149.99



Force Module Separating Pliers
678-210 | \$149.99



Coon Style Ligature Tying Pliers
678-211 | \$149.99



Steiner Tying Pliers
678-222 | \$149.99



Lingual Debonding Pliers
678-709 | \$149.99



Three Jaw Pliers
400-140 | \$149.99




Omega Loop Forming Pliers
400-141 | \$134.99



Bird Beak Pliers
400-142 | \$134.99



Light Wire Bird Beak Pliers with Groove
400-143 | \$134.99



Tweed Loop Forming Pliers
400-144 | \$134.99



Arch Bending Pliers
400-145 | \$134.99



Tweed Pliers, Slim
400-145S | \$134.99



Rectangular Arch Bending Pliers
400-146 | \$134.99



Lingual Arch Forming Pliers
400-147 | \$134.99



Hollow Chop Contouring Pliers
400-148 | \$134.99



O'Brien Pliers
400-149 | \$134.99



Mini Three Jaw Pliers
400-150 | \$134.99



General Contouring Pliers
400-151 | \$134.99



Marcotte Looping Pliers
400-152 | \$134.99



Hollow Chop Pliers
400-153 | \$149.99



Light Wire Bird Beak Pliers
400-154 | \$149.99



Nance Loop Forming Pliers
400-155 | \$149.99



Adams Pliers
400-156 | \$149.99



V-Stop Pliers
400-157 | \$134.99



NiTi Three Jaw Pliers
400-158 | \$149.99



Optical Pliers
400-159 | \$149.99



Old Style Loop Pliers
400-160 | \$149.99



Bird Beak Pliers with Cutter
400-161 | \$134.99



Long Tapered Bird Beak Pliers
400-162 | \$134.99



Hammerhead NiTi Tie-Back Pliers
400-163 | \$134.99



Torquing Pliers
400-164 | \$134.99
.016"/.018" Key
400-165 | \$134.99
.018"/.022" Key



Step Pliers
400-166 | \$134.99
400-166: 1/2 MM
400-167 | \$134.99
400-167: 1 MM
400-168 | \$149.99
400-168: 3/4 MM

WIRE FORMING PLIERS



Slim Hammerhead NiTi Pliers
400-169 | \$149.99



Lingual Hammerhead NiTi Pliers
400-170 | \$149.99



Lingual Hammerhead NiTi Pliers
400-173 | \$134.99



A-Lingual
42.5x61x160MM
HMSM-105 | \$87.99



B-Buccal
45x44.2x157MM
HMSM-205 | \$87.99




Occlusal-Small with Handle, 60x90x100MM
HMSM-304 | \$87.99



Occlusal-Medium with Handle, 70x90x100MM
HMSM-204 | \$87.99



Occlusal-Large with Handle, 77x90x100MM
HMSM-104 | \$87.99



C-Occlusal Adult
75.5x68x102MM
HMSM-305 | \$87.99



D-Occlusal XL Adult
77.5x71x129.5MM
HMSM-405 | \$87.99



E-Occlusal Medium
69x60.5x129.5MM
HMSM-505 | \$87.99



Lateral Photography
40x90x100MM
HMSM-404 | \$87.99



Adult/Child Rounded Occlusal, 150° angle
HMSM-107 | \$87.99



Adult/Child Squared Occlusal, 150° angle
HMSM-307 | \$87.99



Child Occlusal/Buccal 150° angle
HMSM-707 | \$87.99



Adult/Child Buccal 150° angle
HMSM-607 | \$87.99



Adult Occlusal/Buccal 150° angle
HMSM-407 | \$87.99



Lingual Upper/Lower 150° angle
HMSM-507 | \$87.99




Buccal 150° angle
HMSM-207 | \$87.99

STEEL PHOTOGRAPHY MIRRORS



Crown & Gold Scissor, 12CM

SCGC | Curved | \$54.99
 SCGS | Straight | \$54.99



Ortho Mathieu, Narrow Tip Pliers

400-174 | Narrow Tip | \$84.99
 400-175 | Wide Tip | \$84.99
 400-176 | Tungsten Carbide | \$119.99



Lingual Mathieu Pliers
 400-177 | \$89.99



Crown & Gold Scissors, Pedro
 400-178 | \$69.99



Wire-Cutting Scissors
 400-179 | \$69.99



Direct Bond Bracket Holder
 400-180 | \$54.99




Slim Direct Bond Bracket Holder
 400-181 | \$54.99



Slim Buccal Tube Bonder Tweezers
 400-183 | \$54.99



Band Pusher
 400-197 | \$38.99



Band Pusher with Short tip
 400-198 | \$44.99



Distal Clinch Back Hand Instrument
 400-184 | \$29.99



Lingual Ligature Director
 400-185 | \$29.99



Lingual Ligature Director, U Shape
 400-186 | \$29.99





Scaler/Band Pusher
400-187 | \$29.99



Ligature Director Scaler
400-188 | \$29.99



Dull Scaler/Band Pusher
400-189 | \$29.99



Ligature Director/Band Pusher
400-191 | \$29.99



Large Plugger/Picker
400-192 | \$29.99



Dull Scaler
400-193 | \$29.99



Dull Scaler (H6/7)
400-194 | \$29.99



Picker/Ligature Director
400-196 | \$29.99



Band Pusher Short Tip/Scaler
400-200 | \$29.99



Tucker/Ligature Director
400-201 | \$29.99



Gingrass Ligature Director
400-202 | \$29.99



Lingual Torquing Key, .018"
400-203 | \$29.99



Straight .020" Wire Applicator
400-204 | \$29.99



Lingual Torquing Key/Picker
678-708 | \$29.99




Ligature Director/Plugger
LDP | \$29.99



Standard Ligature Director
PLPK | \$29.99





Exam Cassette for 5 Instruments

HMCT-1605 | \$119.99


Exterior Dimensions 3" x 8" x 1.25"
(200 mm x 80 mm x 32 mm)



Small Cassette for 10 Instruments

HMCT-1610 | \$139.99

Exterior Dimensions 5.5" x 8" x 1.25"
(140 mm x 203 mm x 32 mm)




Large Cassette for 20 Instruments

HMCT-1620 | \$184.99

Exterior Dimensions 8" x 11" x 1.25"
(203 mm x 280 mm x 32 mm)






Exam Cassette for 5 Instruments & Accessory Area

HMCT-2205 | \$149.99


Exterior Dimensions 8" x 5" x 1.25"
(202 mm x 130 mm x 32 mm)



Small Cassette for 10 Instruments & Accessory Area

HMCT-2210 | \$169.99

Exterior Dimensions 8" x 7.5" x 1.25"
(202 mm x 195 mm x 32 mm)




Large Cassette for 16 Instruments & Accessory Area

HMCT-2216 | \$189.99

Exterior Dimensions 8" x 11" x 1.25"
(202 mm x 280 mm x 32 mm)






Large Cassette for 20 Instruments & Accessory Area

HMCT-2212 | \$189.99


Exterior Dimensions 8" x 11" x 1.25"
(202 mm x 280 mm x 32 mm)



Large Cassette for 20 Instruments & Accessory Area

HMCT-2220 | \$299.99

Exterior Dimensions 8" x 14.5" x 1.25"
(230 mm x 368 mm x 32 mm)




Oral Surgery Cassette

HMCT-2214 | \$299.99

Exterior Dimensions 8" x 14.5" x 1.25"
(202 mm x 280 mm x 32 mm)






Small Twin Cassette - Holds 10 Instruments

HMSCT-1910 | \$189.99


Exterior Dimensions 7.9" x 3.2" x 1.55"
(200 mm x 82 mm x 40 mm)



Medium Twin Cassette - Holds 14 Instruments

HMSCT-1914 | \$209.99

Exterior Dimensions 7.9" x 4.5" x 1.55"
(114 mm x 203 mm x 40 mm)




Large Twin Cassette - Holds 20 Instruments

HMSCT-1920 | \$219.99

Exterior Dimensions 7.9" x 5.2" x 1.55"
(200 mm x 132 mm x 40 mm)






Large Thin Orthodontic Cassette

HMMT-7ORTH | \$209.99


Exterior Dimensions 8" x 11" x 0.75"
(203 mm x 280 mm x 19 mm)



Small Thin Orthodontic Cassette

HMMT-4ORTH | \$194.99


Exterior Dimensions 7" x 8" x 0.75"
(178 mm x 203 mm x 19 mm)



Mini Exam Cassette for 3 Instruments

HMCT-1403 | \$62.99

Exterior Dimensions 7" x 2.2" x 0.9"
(180 mm x 55 mm x 22 mm)




Exam Cassette for 5 Instruments

HMCT-1405 | \$69.99

Exterior Dimensions 7" x 2.8" x 0.9"
(180 mm x 70 mm x 22 mm)







Small Cassette for 7 Instruments

HMCT-1407 | \$79.99



Exterior Dimensions 7" x 3.5" x 0.9"
(180 mm x 90 mm x 22 mm)



Small Cassette for 10 Instruments

HMCT-1410 | \$94.99



Exterior Dimensions 7" x 5.1" x 0.9"
(180 mm x 130 mm x 22 mm)



Large Cassette for 20 Instruments

HMCT-1420 | \$112.99


Exterior Dimensions 10.2" x 7" x 0.9"
(260 mm x 180 mm x 22 mm)



Cassette for Multiple Handpieces

CS7 | \$129.99


Exterior Dimensions 7" x 3.5" x 1.25"
(180 mm x 90 mm x 32 mm)



Single Handpiece Cassette

HMMT-4112 | \$69.99


Exterior Dimensions 6" x 1.6" x 1.2"
(150 mm x 40 mm x 30 mm)



Small Parts Cassette - Size A

HMMT-4104 | \$54.99


Exterior Dimensions 3.1" x 1.6" x 0.8"
(80 mm x 40 mm x 20 mm)



Small Parts Cassette - Size B

HMMT-4105 | \$59.99

Exterior Dimensions 3.1" x 2.4" x 1"
(80 mm x 60 mm x 20 mm)



Small Parts Cassette - Size C

HMMT-4106 | \$59.99

Exterior Dimensions 4" x 2.4" x 1"
(100 mm x 60 mm x 25 mm)



Clip for Cassettes

AWCLIP | \$6.99

Used to securely hold air water syringe tip and hinged instruments on the lid of the cassette.



DIAGNOSTIC



RESTORATIVE



SURGICAL



PERIODONTAL



ORTHODONTIC



INSTRUMENT
MANAGEMENT



ENDODONTIC



IMPLANTS

To Learn more about HiTeck's Products & Services visit hiteckmedical.com

Manufacturer: 18-215 Traders Blvd East | Mississauga, ON L4Z 3K5 | Canada

Toll Free: 877.899.1861 **Toll Free Fax:** 888.891.9792

Website: www.hiteckmedical.com

Copyright 2020 HiTeck Medical Instruments. All rights reserved. The contents of this product catalog including, but not limited to, all text, photographs and graphics, are the sole and exclusive property of HiTeck Medical Instruments and may not be copied, reproduced or distributed without prior written permission.

HiTeck reserves the right to change features of products, prices, images and text at any time without notice.

For warranty information, please see www.hiteckmedical.com

