

# NITRILE EXAM GLOVES

REF PM6-200X

CHEMO DRUG TESTED | TEXTURED FINGERTIPS | POWDER-FREE | PURPLE



Fit™

## ABOUT THE GLOVE

The Fit is a high-performance glove, suitable for any application requiring medical-grade hand protection. Made of lightweight and durable nitrile, the Fit dons easily and is 100% latex-free, providing reliable protection without the risk of Type 1 latex allergies.

## PROTECTION STANDARDS

**Exceeds the current ASTM D6319:** Standard Specification for Nitrile Examination Gloves for Medical Application.

**Tested against chemotherapy drugs to the current ASTM D6978:** Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.  
(List of tested chemotherapy drugs available upon request.)

**Tested to ASTM F1671 for Viral Penetration.**

## FEATURES & BENEFITS

### LIGHTWEIGHT NITRILE TECHNOLOGY


Combines strength with tactile sensitivity, allowing intricate and precise movement.

### FLEXIBLE AND COMFORTABLE FIT

Moves with the wearer to facilitate natural dexterity.

### FULLY MICRO-TEXTURED WITH ENHANCED FINGERTIP TEXTURING

Provides reliable grip when handling small instruments and objects.



## AVAILABLE SIZES



## QUANTITIES

GLOVES/BOX	BOXES/CASE	GLOVES/CASE
200	10	2000



PRIMED Medical Products Inc.  
200, 2003-91 Street SW, Edmonton, AB, Canada, T6X 0W8  
TF +1.877.877.4633 | F +1.780.497.7670 | W PRIMED.com | E info@PRIMED.com

NA-REVA3