



# Ambulance Air Sanitisation System

**Product codes**


ME908020  
ME901014  
ME901012

**Quantity**

1 x System  
20L  
5L

- Disinfection time: 5 mins, effective on surfaces: up to 18 hrs
- Retrospectively fitted or comes in a kit for new installs
- Sanitise surfaces & environment without having to return to base
- Custom-designed to fit all ambulance layouts
- Particle sizes in the 10–20-micron range with an air velocity at the nozzle of 100 m/s.





- A precision nozzle provides a fine but powerful mist.
- Positioned at the optimum angle to disperse product evenly for the best coverage.
- No liquid residue, leaks, or wastage of product.

particle sizes in the 10–20–micron range with an air velocity at the nozzle of 100m/s




### System operation.

Operation is simple. No need to set timers, our system has been pre-determined for maximum effectiveness. Everything is automatic, press the start button and wait for the decontamination process to finish.



### Streamlined design.

Our system is purposely designed and installed to be as unobtrusive as possible. We are consciously striving towards streamlining and providing a sanitisation solution that is powerful but intended to fit in.



### 360° circulation.

Two precision air jets are strategically positioned inside the ambulance load space, creating a vortex that distributes product evenly and effectively.

### How it works

**The Medi9 Air Sanitisation System uses both Medi9 Air Sanitisation Solution and compressed air to create a vapour-like fog. The vapour is finely atomised and circulated in a vortex for optimum coverage. It begins by deactivating any airborne viruses or bacteria present before forming a layer on all horizontal or vertical surfaces.**

The decontamination process is automatic and safe, with built-in safety features such as the cab's portable key and stop button on the remote. It is a simple but effective system, optimising the time the vehicle and crew remain fully operational. Our system can be provided as a kit for ambulance installers or fitted by our team on-site, and it comes with a year's warranty and service as standard.

Medi9 is clinically proven and non-corrosive. It is a highly effective infection control technique that leaves no liquid build-up or residue, making it more time efficient and safer than other disinfectants. As Medi9 is positively charged, it can provide residual protection of up to 18 hours.



### Custom Built

Most of our components are designed, printed, and built in-house, allowing us to facilitate retrospective or new installs, whatever the layout of the ambulance. We aim to see each ambulance as a new challenge and problem-solve with ease due to our ability to customise and create bespoke designs.

Watch it in action

Scan the QR to watch our Medi9 Air Sanitisation System in use



[ninegrp.com/pages/videos](https://ninegrp.com/pages/videos)





- Adaptable and solution-based, we work with our clients to create the best system for each ambulance.
- Our system is time-effective, with a recommended waiting time of only 5 minutes.

## Medi9 Air Sanitisation Solution

Medi9 Air Sanitisation Solution is designed and manufactured for use in our Medi9 Sanitation systems. This high-performance solution is independently proven effective against 99.9% of bacteria, spores, viruses, and fungi including Coronavirus, Clostridium Difficile, MRSA and Norovirus.

Use Medi9 regularly and effectively in areas where the transfer of germs is a concern to maintain and uphold a safe working environment. With instant germ protection, Medi9 ensures all working areas have been thoroughly sanitised – providing the ultimate solution to the decontamination process.

### Features:

- RTU - ready to use - very important to avoid human error and water contamination
- Fogger - no need to rinse or wipe down after application
- Provides residual protection - use as prevention as well as cure
- Fragrance-free
- Non-hazardous
- Alcohol, chlorine and triclosan-free
- Safe to use in challenging environments and on sensitive equipment
- REACH and COSHH exempt

More details on the entire Medi9 range can be found at: [ninegrp.com/collections/medi9](https://www.ninegrp.com/collections/medi9)



**To find out more about our Medi9 Air Sanitisation System call 01726 813 979 or scan the QR code.**

Together—We Make Life Better. [ninegrp.com](https://www.ninegrp.com)