

## **Datex/Ohmeda RGM 5250 Anesthesia Monitor**

### **Description:**

This cost-effective stand-alone monitor quantifies any one of five anesthetic agents and gives clinicians detailed intra-operative data on:

- Inspired and expired O<sub>2</sub>, CO<sub>2</sub> and N<sub>2</sub>O
- Minute and tidal volume
- Peak, minimum and mean airway pressures
- SpO<sub>2</sub> and pulse rate
- Respiration rate
- Circuit O<sub>2</sub> concentration

Its electroluminescent touch screen display can be remote mounted up to approximately 2 m/6ft from the main unit for added convenience making the 5250 RGM ideal for respiratory monitoring in a wide range of cases.



Available Refurbished  
6 month Warranty

**For additional information, written quote, warranty information, payment terms or product availability, contact our office Toll Free at: 800-989-4909**

### **Do you need medical equipment Repairs or Service?**

For your convenience MESA now offers our in-house bio-medical capabilities through our new service organization, Code Blue Bio-Medical Services, Inc. Please call our service manager at 800-989-4909 for details or an appointment.