

## Crutches

McKesson's push-button crutches feature durable construction, comfortable, adjustable padded hand grips and contoured underarm pads. The convenient and intuitive push-button design allows you to quickly adjust the height in one-inch increments. Non-skid tips provide excellent traction and stable performance. McKesson crutches are available in child, youth, adult and tall sizes.

### Features and Benefits

#### Push-Button Aluminum Crutches

- Lighter than wood, the aluminum crutch provides stability and durability
- Double-extruded center tube provides additional strength to weight-bearing area (wing nut models only)
- Push-button feature allows height adjustment in 1-inch increments
- Easily adjustable hand grip position with Euro-style clip and wing nut adjustment options (tool-free)
- Comfortable, durable underarm pad and hand grip
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Available in child, youth, adult and tall adult sizes
- Not made with natural rubber latex

#### Push-Button Steel Crutches

- Steel construction provides exceptional durability and strength
- Provides support for up to 500 lbs
- Ideal for bariatric patient use
- Push-button feature allows height adjustment in 1-inch increments
- Padded hand grips and ultra-contoured underarm design for lasting comfort
- Tool-free wing nut adjustments for hand grips
- Non-skid, jumbo-size vinyl contoured tips provide excellent traction
- Not made with natural rubber latex



### Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

### Never Compromise on Care<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

## ORDERING INFORMATION

PUSH-BUTTON ALUMINUM CRUTCHES						
MFR #	HAND ADJUSTMENT	SIZE	PATIENT HEIGHT	WEIGHT CAPACITY	HCPCS	QTY
146-10416-1	Wing Nut	Child	4 ft - 4 ft 6 in	175 lbs.	E0114	1 pair/cs
146-10401-8	Wing Nut	Youth	4 ft 6 in - 5 ft 2 in	350 lbs.	E0114	1 pair, 8 pair/cs
146-10400-8	Wing Nut	Adult	5 ft 2 in - 5 ft 10 in	350 lbs.	E0114	1 pair, 8 pair/cs
146-10402-8	Wing Nut	Adult Tall	5 ft 10 in - 6 ft 6 in	350 lbs.	E0114	1 pair, 8 pair/cs
146-10430-8	Euro-Style Clip	Adult	5 ft 2 in - 5 ft 10 in	350 lbs.	E0114	1 pair, 8 pair/cs
146-10432-8	Euro-Style Clip	Adult Tall	5 ft 10 in - 6 ft 6 in	350 lbs.	E0114	1 pair, 8 pair/cs
PUSH-BUTTON STEEL CRUTCHES						
MFR #	HAND ADJUSTMENT	SIZE	PATIENT HEIGHT	WEIGHT CAPACITY	HCPCS	QTY
146-10406	Wing Nut	Adult	5 ft 2 in - 5 ft 10 in	500 lbs.	E0116	1 pair/cs

McKesson Medical-Surgical  
 9954 Mayland Drive  
 Richmond, VA 23233  
[mms.mckesson.com/mckesson-brands](http://mms.mckesson.com/mckesson-brands)

Post-Acute Care and Home Health  
 800.347.2456

Home Medical Equipment  
 888.822.8111

Primary Care  
 866.625.2679

McKesson Brands Support  
 800.777.4908