

Deluxe Aluminum Shower Chair

HCPCS: E0240




FIGURE A




FIGURE B



FIGURE C



FIGURE D

- Easy, tool-free assembly of back, seat and legs (Figure A - B)
- Aluminum frame is lightweight, durable and corrosion proof
- New pinch-free cover allows for push pins to be depressed without pinching finger
- The audible "click" will be heard when secure (Figure C)
- Legs are height adjustable in 1" increments, and angled with suction-style tips to provide additional stability
- Support collar prevents leg movement and reduces rattle
- Made using recycled materials

PRODUCTS

Description

12202KD-4	Bath Seat w/ Back, KD, 4/cs
-----------	-----------------------------

SPECIFICATIONS

Description

Carton Shipping Weight	29 lbs.
Outside Legs	19.5" (W) x 19.5" (D)
Seat Dimensions	19.25" (W) x 11.5" (D)
Seat Height	15.5" - 19.5"
Warranty	Limited Lifetime
Weight	7 lbs.
Weight Capacity	300 lbs

ACCESSORIES

Description

12002T	Tips for 12002,4/pkg., BX
--------	---------------------------