

## Cane Seat

HCPCS: E0100



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Provides a comfortable seat to rest on when open, and a sturdy support cane when closed (Figure A)
- Manufactured with sturdy, extruded aluminum tubing
- Tripod design with vinyl contoured tipped legs
- 9" diameter seat
- 33" tall when closed (Figure B)
- Weighs 2.5 lbs
- Weight Capacity 250 lbs

### PRODUCTS

	Description
RTL10365-ADJ	Adjustable Height Cane Seat, 1 EA

### SPECIFICATIONS

	Description
Cane Height	34" - 38"
Seat Height	19" - 22"
Seat To Floor Height	21"
Warranty	Limited Lifetime
Weight Capacity	250 lbs