

Gel Grip, Aluminum Offset Canes, Height Adjustable

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- Ergonomically designed handle with soft Gel Grip provides comfort and security
- Manufactured with sturdy, extruded aluminum tubing
- Comes standard with premium Gel Grip and wrist strap
- Manufactured with sturdy, extruded aluminum tubing
- Easy-to-use, one-button height adjustment with locking ring Locking ring prevents rattling
- Available in Blue Crackle (main image), Red Crackle (Figure A) or Black (Figure B)
- Handle height adjusts from 30" to 39"
- Weight Capacity 300 lbs
- Approx. user height 5'-6'5"

PRODUCTS

Description

RTL10336MA	Cane,T-Handle Erg, Marble. 1/EA,6/CV
RTL10336MN	Cane,T-Handle Erg, Mech. Navy. 1/EA,6/CV
RTL10336SF	Cane,T-Handle Erg, Spring. 1/EA,6/CV
RTL10372BC	Cane,Offset,Gel Grip,BL Crackle,1EA,6/CS
RTL10372BK	Cane,Offset,Gel Grip,Black,1/EA,6/CS
RTL10372FL	Cane, Offset Gel Grip, Floral, 1/ea,6/cv
RTL10372GL	Cane, Offset Gel Grip, Leaves, 1/ea,6/cv
RTL10372PL	Cane, Offset Gel Grip, Plaid, 1/ea,6/cv
RTL10372RC	Cane,Offset,Gel Grip,Red Crk,1EA,6/CS

SPECIFICATIONS

Description

Handle Height	30" - 39"
Warranty	Limited Lifetime
Weight Capacity	300 lbs

ACCESSORIES

Description

RTL10322BKB	CaneTip, 3/4 Blk, 1/EA, 10/CS, EA
RTL10384GB	GolfBallCaneTip, 1/ea, 4/CV, EA
RTL10384SB	SoccerBallCaneTip, 1/ea, 4/CV, EA
RTL10384TB	TennisBallCaneTip, 1/ea, 4/CV, EA