

Scout 4

HCPCS: K0800-1, K0800



FIGURE A



FIGURE B




FIGURE C




FIGURE D




FIGURE E

- Quick and easy disassembly (Figure A)
- Color-through color panels hide scuffs and scratches
- Armrests are padded and adjustable
- Quick connect batteries
- Ergonomic throttle control
- Height-adjustable swivel seat with fold-down backrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust tiller
- Anti-tip wheels
- Delta tiller
- Easy free wheel operation
- Available with 12AH or 20AH batteries
- 20AH batteries increase scooter's maximum range

PRODUCTS

Description

SFSCOUT4	Scout 4-Wheel Scooter, Red/Blue, 1/ea
SFSCOUT4-EXT	Scout 4-Wheel 20AH Battery, Red/Blue

SPECIFICATIONS

Description

Base Weight	58 lbs
Batteries (pair)	12V x 12AH, 12V x 20AH (EXT Model Only)
Battery Weight (pair)	20 lbs
Brakes	Electromagnetic
Charger	2A Offboard
Climbing Angle	6°
Controller	CTE 45A
Dimensions	42.5" (L) x 21.5" (W)
Floor to Seat Height	17" - 19"
Front Wheels	8" x 2" Flat-free
Ground Clearance	2.5"
HCPCS	K0800
Heaviest Piece Weight	39 lbs
Maximum Range	9 miles (12AH) 15 miles (20AH)
Motor	24V x 270W x 4700 rpm
Rear Wheels	8" x 2" Flat-free

SPECIFICATIONS	Description
Seat Dimensions	16.5" (W) x 16" (D)
Seat Weight	16 lbs
Top Speed	4.25 mph
Total Weight (excl. hand control)	94 lbs
Turning Radius	53.75"
Warranty on Batteries	12 Months
Warranty on Electronic Controller/Drive Train Components	24 Months
Warranty on Frame	Lifetime
Weight Capacity	300 lbs
Wheel Base Dimensions (L x W)	31.9" x 16.9"

ACCESSORIES	Description
SF8050	Sun Shade for Scooter,1/ea, EA
SF8020	Rear Basket Set, LType, 1/ea, EA
SF8030	Crutch Holder Set, SType, 1/ea, EA
SF8010	Oxygen Holder Set, L Type,1/ea, EA

PARTS	Description
LRM402109	12AH Battery,POV,1/ea, EA
LRM412103-11	Speed Control Knob, EA
BC350103	Connecting Rod Set,Bobcat 4,1set/pkg, PKG
C08-069-01000	Tiller Lock Set, Bobcat X3 / X4, 1/ea, EA
S31026	Anti Tip Wheel,Phoenix/Bobcat,1/ea, EA
S35003-P	PlasticBasketBobcat,1ea, EA
S31052	Armrest Knob POV 1/ea, EA
BC20ACB	Bobcat Breaker 20A with Cap,1ea, EA
BC35051	BasketHook Set,Phoenix,Scout,Ventura,1/e, EA
BC953014	Arm Pad Only,Scout, Bobcat, 1/ea, EA
BC35040L	Bobcat Complete Arm Rest,Left,1ea, EA
BC35040R	Bobcat Complete Arm Rest,Right,1ea, EA
BC35045	Bobcat Seat,Black,1ea, EA
BC31021	Tiller Bearing Set,Bobcat 4,1set/pkg, PKG
LRM412123	21AH Battery,POV, 1/ea, EA
LRM302111	2ah Battery Charger, POV,1/ea, EA
BC31006	Bobcat Ignition Switch,1ea, EA
BC31001	Scout,Bobcat, Speed Dial, 1/ea, EA
BCKEY	Key, Bobcat, Scout, 1ea, EA
BC350112	Bobcat Transaxle Only,1ea, EA