

## Gravity 7 Long Term Care Pressure Redistribution Mattress

HCPCS: E0184-1,E0184





FIGURE A

- High quality therapeutic foam mattress provides pressure redistribution and shear/friction reduction
- Drive's most popular therapeutic foam mattress maximizes value and comfort with 3 full layers of latex free foam
- Full width top layer uses unique die-cutting to optimize zoned pressure redistribution
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress
- High density foam perimeters provide solid edge support and facilitate safe patient transfers
- Available with or without raised side rails
- Visco foam heel section offloads pressure from vulnerable heels
- Removable multi-stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633

PRODUCTS	Description
15770	GRAVITY 7 MATTRESS, 80",1/EA
15777	GRAVITY 7 MATTRESS W/ELEV PERIM,80",1/EA
15784	GRAVITY 7 MATTRESS 84",1/EA
15785	GRAVITY 7 MATTRESS W/ELEV PERIM 84",1/EA
15876	GRAVITY 7 MATTRESS 76", 1/EA
15877	GRAVITY 7 MATTRESS W/ELEV PERIM 76",1/EA

SPECIFICATIONS	Description
Cover Material	Stretch, Fluid-Resistant Nylon Top   Waterproof Vinyl Bottom
Dedicated Heel Section	Yes
Mattress Layers of Foam	3
Side Rails	Yes
Unit Dimensions	15784 and 15785: 36" (W) x 84" (L) x 6" (H)   15770 and 15777: 36" (W) x 80" (L) x 6" (H)   15876 and 15877: 36" (W) x 76" (L) x 6" (H)
Warranty	7 Year Limited Foam   3 Year Cover
Weight Capacity	350 lbs

ACCESSORIES	Description
15770C	Mattress Cover 15770C,1/ea, EA
15777C	Mattress Cover 15777C,1/ea, EA
15784C	Mattress Cover 15784C,1/ea, EA
15785C	Mattress Cover 15785C,1/ea, EA
15877C	Mattress Cover 15877,1/ea, EA
15876C	Mattress Cover 15876,1/ea, EA

ACCESSORIES

Description

14427QF

Quick Fill Pump w/adapters, 1 EA/12 CV, EA