

Power Neb Ultra Nebulizer

HCPCS: E0570, A7003, A7005, A7013




FIGURE A



FIGURE B

- A competitive edge price point for providers as one of the smallest compressor nebs at a mere 3.2 lbs
- Built-in neb kit holder for hands-free use
- Delivers consistent particle size of <5 microns and fast treatment times
- Complete with disposable neb kit, and optional reusable neb kit or pediatric aerosol mask
- Powerful piston pump is efficient, dependable and quiet

PRODUCTS	Description
18080	Power Nebulizer, Ultra, 1/ea, 6/cv
18080P	Power Neb Ultra, w/PedMask 1/ea, 6/cv
18081	PowerNebUltraw/Reuse&DispNeb, 1/ea, 6/cv
18082	PowerNeb, Ultra w/2 Reusable kits,6/cs
18260-PDMASK	Ped Mask for 18260 and Adapter, 1/ea

SPECIFICATIONS	Description
Electrical Requirement	115 VAC 60Hz
Liter Flow	, 8 lpm
Maximum Pressure	30 psi
Method of Nebulization	Piston Compressor
Operating Pressure	, 10 psi
Particle Size	0.5 μm to 5 μm
Unit Dimensions	6.3" (W) x 3.0" (H) x 5.3" (D)
Warranty	Limited 5 Year
Weight	3.2 lbs.

ACCESSORIES	Description
18031	Nebulizer Carry Bag, 1/ea, EA
18254	Reusable Neb Kit 10/case, CS
NEB KIT 600	NebKit,Inc Adlt Mask,Tub 50/cs, CS
MASK 001A	Adult Aerosol mask CS/50, CS
MASK 001P	Pediatric Aerosol mask CS/50, CS
NEB KIT 700	NebKit,IncPediMask,Tub ,50/cs, CS
NEB KIT 500	Neb Kit,Inc Mouth pc&Tub50/cs, CS
18260	JetNeb Plus Breath Enhanced Neb,12/CS, CS
18090F	NebFilter 10/bag, BAG