

## Deluxe 3 Wheel Aluminum Rollator, 7.5" Casters

HCPCS: E0143



FIGURE A



FIGURE B



FIGURE C




FIGURE D



FIGURE E



FIGURE F

- Comes standard with basket, tray and pouch (Figure A)
- Basket conveniently folds for easy storage (Figure B)
- Ships in ultra-compact carton Comes with frame folded and compact Pull lightweight, aluminum trombone-style frame to lock in place and secure with lock knob
- 7.5" casters for indoor or outdoor use
- Handles easily adjust in height with self-threading knob (Figure C)
- Special loop lock made of internal aluminum casting operates easily and ensures safety (Figure D)
- Easy, one-hand folding

### PRODUCTS

#### Description

10289BL	3 Whl Alum Rollator,Blue,1/ea
10289RD	3 Whl Alum Rollator,Red,1/ea

### SPECIFICATIONS

#### Description

Carton Shipping Weight	17 lbs.
Handle Height	31" - 38"
Unit Dimensions	25" (W) x 24" (L)
Warranty	Limited Lifetime
Weight	11 lbs.
Weight Capacity	300 lbs

### PARTS

#### Description

9503E17101-L	Hand Brake with Cable, 10289, left 1/ea, EA
9503E17101-R	Hand Brake with Cable, 10289, Rt, 1/ea, EA
816T	Tray, 10289, 1/ea, EA
9503E1028901-B	Handlebar,BlackFinish,10289,Right/Lft,1e, EA
815B	Folding Basket, RTL199, 1/EA, EA
STDS1060	Wheelchair Fork Bearing,1/ea, EA
808	BLACK TOTE,1/ea, EA
9503E17106	Rear Axle with Nut,10289, 1/ea, EA
9503E17108	Wheel Assembly,10289, 1/ea, EA
9503E17112	Folding Hinge,10289, 1/ea, EA
9503E17103	Height Adjustment Knobs, RTL199,1/PR, PR
9503E17104	Wheel Brake,Left,10289, 1/ea, EA
STDS1011	Wheel Bearing, 8mm x 22mm 1/ea, EA
950330049901	Hand Grip, Right, 3-wheel Rollator, 1/ea, EA

PARTS

Description

950330039901

Hand Grip, Left, 3-wheel Rollator, 1/ea, EA