

Duet Rollator/Transport Chair, 8" Casters

HCPCS: E0143+E0156



FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F

- Combines the features of a rollator and transport chair in one unit. An individual can ambulate independently or be pushed safely by a caregiver.
- To use as a transport chair, the Duet has locking flip-up armrests (Figure A - Figure C) and a flip-down footrest (Figure D) to aid the user. A safety lock release is used to store the arms, and a hook and loop fastener strap is used to store the footrest while using the Duet for ambulation.
- The tool-free two-position contoured padded backrest is reversible, depending on whether you are using the Duet as a rollator or in transport chair mode.
- Comes with a padded seat with a carry pouch for added privacy and security of personal items. Inside the pouch is a handle to aid with the folding.
- The folding frame has a lock-out feature to ensure the unit stays open while in use and ensures the unit will stay folded during transport.
- 8" casters are ideal for indoor and outdoor use.
- Handles are adjustable in height, with a self-threading knob (Figure E).

PRODUCTS

PRODUCTS	Description
795B	Duet Rollator/TransportChair,Blue, 1/ea
795BK	Duet Rollator/TransportChair,B
795BU	Duet Rollator/TransportChair,Burgdy,1/cs

SPECIFICATIONS

SPECIFICATIONS	Description
Back Height	35"
Carton Shipping Weight	25 lbs.
Handle Height	31.5" - 37"
Overall Dimensions	25.5" (W) x 37" (H) x 26" (L)
Seat Dimensions	13.5" (W) x 13" (D)
Seat To Floor Height	21"
Seat Width	13.5"
Warranty	Limited Lifetime
Weight	20 lbs.
Weight Capacity	300 lbs

PARTS

PARTS	Description
9501A79512	Handlebar, Right, 795, 1/ea, EA
9502M7951L	Elbow Rest, Left, 795, 2M Ser.#, 1/ea, EA
9502M7951R	Elbow Rest, Right, 795, 2M Ser.#, 1/ea, EA
750FORK-2M	Front Fork, 795, 2M Ser.#, 1/ea, EA

PARTS	Description
9501A7953BK	FootrestAssm.Black,795 1A/2Y Ser.#,1ea, EA
9501A7953R-2M	FootrestAssembly,Red,795 2M/20hSN#, 1/ea, EA
9501A79511	Handlebar, Left, 795, 1/ea, EA
9501A7953B-2M	FootrestAssembly,Blue,7952M/20HSN#, 1/ea, EA
9501A79510-RP	Roll Pin Fork Fastener,795,1ea, EA
9501A79510-2M	Fork Stem Bolt, 795, 2M Ser.#, 1/ea, EA
9501A79506	RearAxleBolt795,1ea, EA
9501A79506-2M	Rear Axle Bolt, 795, 2M Ser.#, 1/ea, EA
9501A79504	Tote for 795 Rollator, 1/ea, EA
750BCL-2M	HandBrakeAssmby,Left,795,2M Ser.#, 1/ea, EA
750BRAKE-2M	Brake, Right or Left, 795, 2M S#, 1ea, EA
9501A75003	Spring Boot for Cable,1/ea, EA
9500075005	Folding Hinge Bolt,750/795,1/ea, EA
9501A75013	Seat Bracket, 795 Rollator, 2/pr, EA
9502M75013	SeatBrackets&Screws,7952M/20H SN#, 1/pr, PR
9501A75008	3 Hole Wheel,8"Rear,750N,795,1/ea, EA