

vive
health



VIVE 4 Wheel Mobility Scooter

MOB1027

The Vive 4 Wheel mobility scooter provides safe and stable mobility assistance. With a turning radius of only 4.3' the scooter easily maneuvers around corners. 8" flat-free tires traverse over smooth indoor and outdoor surfaces such as concrete, pavement, tile, and carpeting. The padded seat features flip-back armrests and the ability to swivel and lock in place. The front mounted storage basket disconnects for optional use. The scooter has a speed control knob that adjusts up to 3.7mph and has a range of up to 12.4 miles per charge. Quickly disassemble into four pieces for easy transportation. The heaviest piece is only 29lbs.

Body Colors



Seat Color



Best Economy in its class



WHAT'S INCLUDE

- | | |
|-------------------|---------------------------|
| A. Front Section | G. Basket |
| B. Rear Section | H. Battery Securing Plate |
| C. Seat Unit | I. Basket Holder |
| D. Battery | J. 2x Bolts |
| E. Charging Cable | K. Seat Pin |
| F. 4x Screw Knobs | L. Arm Rest |

(239) 220-5772

vivehealth.com

service@vivehealth.com

8955 Fontana Del Sol Way, Naples, FL 34109, USA



Specifications

Overall Dimensions (L x W x H)	37.7" x 19.6" x 33"
Total Weight	88.2 lbs.
Weight of Heaviest Piece	29lbs. (rear base)
Seat Height (from ground)	18.5" - 21.5"
Seat Width	17"
Seat Depth	15"
Armrest Height (above seat)	6.75"
Seat Back Height	13"
Net Weight	88.2 lbs.
Maximum Speed	3.7 mph
Maximum Braking Distance	≤5"
Minimum Turning Radius	≤4.3ft.
Weight Limit	265 lbs.
Maximum Travel Distance (in Theory)	12.4 mi
Static Stability	≥ 9°
Dynamic Stability	≥ 6°
Climbing Ability	≤ 9°
Battery	Lead-acid battery (24V/20AH)
Controller Maximum Output Current	45A
Controller Minimum Output Current	2A
Front Wheel Diameter	7.5" (Solid)

SAFELY NAVIGATE SMOOTH SURFACES:

Easy maneuver over smooth, indoor and outdoor surfaces with comfort and ease with the VIVE 4 Wheel mobility scooter. Designed to safely support up to 265 pounds, the mobility scooter features 8" flat-free, non-marking tires, a maximum speed of 3.7mph and a driving distance up to 12.4 miles per charge.

COMFORTABLE SEATING:

Comfortably padded, the wide, 15.75" seat and backrest provide gentle cushioning and support. The seat rotates and is also height adjustable for a customizable fit. Additionally, the scooter seat is surrounded by two padded armrests that adjust and flip upward for greater convenience when exiting or boarding.

EASY TO TRANSPORT:

Quickly disassemble the VIVE 4 Wheel mobility scooter into four separate pieces for easy portability. Featuring a folding tiller and removable seat, each section easily fits in the trunk of standard size vehicles. The rear section weighs 29 pounds and is the heaviest of the four pieces.

INTUITIVE CONTROL PANEL:

Located on the tiller console, the control panel allows you to set the maximum speed with the variable speed dial, change from a forward to a reverse direction, turn on the headlight as well as produce a warning tone as needed.

The display also includes a color-coded battery life indicator for checking the status of the battery at a glance.