

# DOUBLE safety for users and patients


The guidelines of the European Standard EN 13795 for high performance are more than fulfilled by Foliodrape Protect plus across the entire surface!

**Maximize the safety of your operations.**

Surgical procedures produce varied volumes of blood and irrigating fluids<sup>1</sup>. Based on this knowledge, we from HARTMANN have optimized the fluid management of the new Foliodrape Protect Plus precisely at the point where a lot of fluid is produced. Compared to other products, this surface is very large and ensures a reliable absorption strength – precisely where it matters.

**Minimize the risk of germ transfer.**

One of the important requirements of European Standard EN 13795 is a reliable fluid management. As your competent partner in the operating theatre, we from HARTMANN aspire to assist – with Foliodrape Protect plus: Top safety for users and patients with the best of material and quality.




The performance of the Foliodrape products compared to the European Standard EN 13795.

Standard requirement for High performance, critical area		Foliodrape Protect	Foliodrape Protect plus <sup>3</sup>
Particle release (linting)	≤ 4.0 log <sub>10</sub>	< 3	< 3
Purity of particulate material	≤ 3.5 IPM	< 2.5	< 2.5
Fluid resistance	≥ 100 cm (H <sub>2</sub> O)	> 150	> 150
Tear resistance (dry)	≥ 20 N	> 45	> 110
Tear resistance (wet)	≥ 20 N	> 45	> 110
Bursting strength (dry)	≥ 40 kPa	> 150	> 300
Bursting strength (wet)	≥ 40 kPa	> 150	> 300
Microbial penetration (wet)	6.0 BI	6.0	6.0

<sup>1</sup> Aksoy Y, Altinel L, Kose KC. The comparison of the effects of intraoperative bleeding control and postoperative drain clamping methods on the postoperative blood loss and the need for transfusion following total knee arthroplasty. Acta Orthop Traumatol Turc 2011;45:190-4, DOI:10.3944/AOTT.2011.2398.

<sup>2</sup> Reinforced area

# Foliodrape® Protect Drapes



## Uses

For sterile draping of the patient and equipment in the operating theatre, outpatient department, ward and doctors office.

Product	Size	Package contents	Article number
<b>Foliodrape® Protect Drape</b> Sterile, individually wrapped	37.5 x 45 cm	125 Pieces	277 516
	45 x 75 cm	65 Pieces	277 500
	50 x 50 cm	110 Pieces	277 547
	75 x 75 cm	56 Pieces	277 501
	75 x 90 cm	35 Pieces	277 502
	90 x 100 cm	25 Pieces	277 504
	150 x 175 cm	17 Pieces	277 548
	90 x 150 cm	20 Pieces	277 506
	95 x 120 cm	25 Pieces	277 505
	150 x 240 cm	13 Pieces	938 810



<b>Foliodrape® Protect Drape</b> Sterile, individually wrapped (Adhesive)	45 x 75 cm	60 Pieces	277 507
	50 x 50 cm	90 Pieces	277 546
	75 x 75 cm	45 Pieces	938 811
	75 x 90 cm	40 Pieces	277 508
	90 x 100 cm	30 Pieces	277 542
	150 x 240 cm	12 Pieces	277 551
	200 x 175 cm	13 Pieces	938 812
	150 x 175 cm	13 Pieces	277 543
	170 x 175 cm	12 Pieces	277 544
	170 x 300 cm	7 Pieces	277 538



Product	Size	Package contents	Article number
<b>Foliodrape® Protect Drapes with fenestration</b> Sterile, individually wrapped	45 x 75 cm ø7cm	65 Pieces	277 503
	50 x 60 cm ø5 cm	70 Pieces	277 549
	75 x 90 cm ø10 cm	40 Pieces	277 511



<b>Foliodrape® Protect Drapes with fenestration</b> Sterile, individually wrapped (Adhesive fenestration)	50 x 60 cm ø7 cm	70 Pieces	277 537
	75 x 90 cm ø7 cm	40 Pieces	277 517
	75 x 90 cm ø10 cm	35 Pieces	277 555
	120 x 150 cm ø7 cm	22 Pieces	277 545



Product	Size	Package contents	Article number
<b>Foliodrape® Protect Drapes with fenestration</b> Flexible two-piece drapes (Adhesive Fenestration with maximum extension to 20cm)	45 x 75 cm	65 Pieces	277 509
	75 x 90 cm	40 Pieces	277 510

