

## Indication For Use

The face mask is used as the patient interface to connect with breathing circuit which carries air mixed with special medication during anesthesia or artificial ventilation.

## AnesEasy II Anesthetic Face Mask

SKU	Color Coding	Size	Description	Hook	Case Size	Weight (g)
A0817	Red	7	AnesEasy Mask, Extra large adult, red cushion	Yes	20	-
A0816	Orange	6	AnesEasy Mask, large adult, orange cushion	Yes	50	26.2
A0815	Green	5	AnesEasy Mask, medium adult, green cushion	Yes	50	22.2
A0814	Yellow	4	AnesEasy Mask, small adult, yellow cushion	Yes	50	19.2
A0813	White	3	AnesEasy Mask, pediatric, white cushion	Yes	50	15.2
A0812	Grey	2	AnesEasy Mask, infant, grey cushion	No	20	-
A0811	Blue	1	AnesEasy Mask, neonate, blue cushion	No	20	-

- Single use material to reduce possible cross contamination
- Phthalate and LATEX FREE
- BS EN ISO 5356-1 Anaesthetic and respiratory equipment-Conical connectors-Part 1: cones and sockets

## Materials

Cushion	TPE
Face Plate	Polypropylene



## Physical Characteristics

Color	Color coded cushions Transparent for face plate
-------	--

## Connectors

Face Plate connector	ISO 15 mm Male One on each size: 1 and 2 ISO 22 mm Female One on each size: 3, 4, 5, 6 and 7
Hook Rings	Attached to 22 mm female patient connection One on each: size: 3, 4, 5, 6 and 7

## Environmental Conditions

Storage	Temperature: -4 °F (-20 °C) to 140 °F (60 °C)
---------	---

