


NEW SEAMLESS POLYMER GEL



PPS PATIENT POSITIONING SOLUTIONS, INC.

Medical Gel Positioning Products



Patient Positioning Solutions, Inc.
(844) 365-9467 Toll Free
info@pps-products.com | Email



Proudly Made
In The USA

PATENT PENDING



Patient Positioner

Part #	Dimensions	Description	Weight
SPG-3015P	10 X 4-1/2 X 1-1/2"	Small	1.6 lbs
SPG-3017P	14 X 4-1/2 X 1-1/2"	Large	2.3 lbs



Universal Armrest Pad

Part #	Dimensions	Description	Weight
SPG-8009M	12 X 9 X 1/2"	Flat Pad	1 lb



Regular SPG-9015AX



Sand Bags, Flat Bottom

Part #	Dimensions	Description	Weight
SPG-6003S	10 X 4-1/2 X 1-1/2"	Sandbag, Small	1 lb
SPG-6004S	12 X 4-1/2 X 1-1/2"	Sandbag, Medium	2 lbs
SPG-6011S	14 X 4-1/2 X 1-1/2"	Sandbag, Large	4 lbs

Axillary Rolls

Part #	Dimensions	Description	Weight
SPG-9001AX	8 X 2 X 2"	Axillary Roll	.9 lbs
SPG-9003AX	12 X 3 X 3"	Flat bottom	2.7 lbs
SPG-9005AX	12 X 4 X 3"	Flat bottom	4.2 lbs
SPG-9012AX	12 X 4-1/2 X 1-1/2"	Regular	2 lbs
SPG-9015AX	10 X 4-1/2 X 1-1/2"	Regular	1.6 lbs

Axillary Rolls should be used when positioning patients on their side. The axillary gel pads, when placed in the axilla will protect the brachial plexus during prolonged pressure while in the lateral position.



Chest Roll, Flat Bottom

Part #	Dimensions	Description	Weight
SPG-3001P	8 X 2 X 2"	Chest Roll	.9 lbs
SPG-3004P	12 X 4 X 3"	Chest Roll	3.8 lbs
SPG-3047P	8 X 2 X 2"	Neo/Pedi	.9 lbs
SPG-3071P	12 X 3 X 3"	Chest Roll	2.6 lbs