


CLINIC **ARMOUR**

Contents



x 8 releasable cable ties




x 1 T strap connector

Instructions

STEP 1: Pass all straps under the examination table and through the eyelets. Ensure straps do not block the desired movement of any of your table sections.



STEP 2: Pull cable tie to tension, the straps should not be able to slide- you can trim the cable tie but allow for 5-10 mm of length to enable reuse.



To release cable ties press the tab down and reverse the procedure above.



STEP 3: Thread the side strap nearest the foot end of the table through the T strap as shown and tighten. Next, feed the strap and end of the bed through the single loop and tighten. Use the cable ties as shown above for all fixing.

