



I.V. Injection Leg Canine



With our canine i.v. Injection leg (forelimb) the correct course of an intravenous injection can be trained. The forelimb has an artificial fur which can be clipped with a clipper. Multiple shearing operations are possible by rotating the fur. The cephalic vein can be palpated and punctured. A peripheral venous catheter (PVC) can be placed without any problems. The paw comes with a replacement fur and an additional vein tube.

Ref.no.: VET4610

Details:



*Experts in medical education
since 1950*



Erlor-Zimmer GmbH & Co. KG · Hauptstr. 27 · D-77886 Lauf
Tel.: (+49) 7841 / 6003-0 · info@erler-zimmer.de · www.erler-zimmer.de