


"Buttockmate" type 2

- Intramuscular Injection Simulator

● Weight/About 1.2kg ● Case size/W46×D16×H26cm ● Accessories/Belt, Drain bag, Storage case

Reproducing Clarke's point and Hochstetter's point.
Allowing simulation in the lateral recumbent position.



The epidermis is made of a special material that leaves no needle mark.



Selection of an incorrect site is indicated by blinking of the red lamp and buzzing. Successful puncturing at the correct site is indicated by the green lamp.
(Only at Clarke's point and Hochstetter's point)



Being compact, the device needs a small space for storage.

Practice

Hochstetter's point



Skeletal reference points such as the anterior superior and anterior posterior iliac spines and greater trochanter can be felt for efficient training in finding the injection site.

Clarke's point



Fluid can be injected.
* Do not use a drug solution. Use water (tap water only) for practice.