

Analytical Report Nr.

AR-20-YL-008942-01

Sample code Nr.

560-2020-00009589

Date

21/12/2020

ANALYTICAL REPORT**Client Information**

Eurofins Medical Device Testing
2425 New Holland Pike
Lancaster USA
+1 717 725 8250
aditianand@eurofinsus.com
For the attention of Aditi Anand

Sample Information

Order Code: EUAA70-00009654
Reception Date: 9-Dec-2020
Analysis Starting Date: 9-Dec-2020
Analysis Ending Date: 21-Dec-2020
Sample code Nr. 560-2020-00009589
Sample described as: Masks

Requirements and decision rule

Customer requirements: No requirements
Decision Rule: Not applicable.

Information provided by the customer*

Client Reference: DMF4L
Sample Description: On behalf of: Breathh Inc./Lutema
Purchase Order Number:

Batch BLR139

Analytical Report Nr.

AR-20-YL-008942-01

Sample code Nr.

560-2020-00009589

Date

21/12/2020

SAMPLE PICTURE



Analytical Report Nr.

AR-20-YL-008942-01

Sample code Nr.

560-2020-00009589

Date

21/12/2020

CONCLUSION:

TEST PROPERTY	PASS	FAIL	REMARKS
◆ Determination of face mask materials efficiency to penetration by particulates (latex spheres) ASTM F2299 / F2299M - 03(2017)			
A - Mask			REFER RESULT

Remark: Test has been performed as per application request

Analytical Report Nr.

AR-20-YL-008942-01

Sample code Nr.

560-2020-00009589

Date

21/12/2020

COMPONENT LIST:

COMPONENT ID	COMPONENT NAME	MATERIAL DESCRIPTION	COLOR	REMARKS
CUST 01	A - Mask	Mask	Pink	---

Analytical Report Nr.

AR-20-YL-008942-01

Sample code Nr.

560-2020-00009589

Date

21/12/2020

MASKS TESTING

CAS No.

RESULTS

UNC.

LOQ

GUIDELINES

Analyses on:A - Mask**◆ Determination of face mask materials efficiency to penetration by particulates (latex spheres)**

Analysis Ending Date: 21/12/2020

ASTM F2299 / F2299M - 03(2017)

Filtration Efficiency

99.88 %

-

Particle diameter: 0.1 μm \pm 0.015 μm
Complete test data reported at Annex.

Analytical Report Nr.

AR-20-YL-008942-01

Sample code Nr.

560-2020-00009589

Date

21/12/2020

Signed for and on behalf of Eurofins Textile Testing Spain:

eurofins
Eurofins Textile Testing Spain, S.L.U.
C/ Castellón, 11 (Alicante)
E: 0374099

Report electronically validated by

Maria Jesus Martinez Puig

Chemical Lab manager

EXPLANATORY NOTE

- ◆ Test not covered by ENAC accreditation scope
- Test is subcontracted within Eurofins group and is accredited
- Test is subcontracted within Eurofins group and is not accredited
- Test is subcontracted outside Eurofins group and is accredited
- Test is subcontracted outside Eurofins group and is not accredited

N/A = Not Applicable

Eurofins Textile Testing Spain S.L.U is not responsible of the information supplied by the costumer and reported as section "Information provided by the costumer".

Eurofins General Sales Terms and Conditions Applied.

Results obtained refer only to samples, products or material received in Laboratory, as described in section "Sample information" and tested in conditions shown in present report.

Test uncertainties not reported are at customer disposal, for those tests in which it is possible to evaluate the test uncertainty.

The reported expanded uncertainty is based on a standard uncertainty multiplied by a coverage factor $k = 2$, which for a normal distribution provides a level of confidence of approximately 95%.

Reproduction of this document is only valid if it is done completely and under the written permission. Any printed version of this document will be considered as a copy.

If you happen to have any comments, please do it by sending email to textile_spain@eurofins.com and referring to this report number.

End Of Report

Eurofins Textile Testing Spain, S.L.U.

Calle Germán Bernácer, 4

03203 Elche

SPAIN

Phone+34 966 299 638

www.eurofins.com/tex

DETERMINATION OF THE INITIAL EFFICIENCY OF MATERIALS USED IN MEDICAL FACE MASKS TO PENETRATION BY PARTICULATES USING LATEX SPHERES

Test Method: ASTM F2299 / F2299M - 03(2017)

Number of test specimens: 5

Thickness: 0.59 mm

Unit area weight: 97.3 g/m²

Exposed specimen area: 78.5 cm²

Size of test specimen was facing towards the challenge aerosol: Inner

Face velocity: 6 cm/s

Pressure drop: 0.0 mmH₂O

Particle diameter: 0.1 μm ± 0.015 μm

Conditioning and test conditions: T^a between 18°C and 24°C. HR% between 30 and 50 % Hr

Pretreatment techniques used: No

Test duration: Three one-minute counts per specimen

Controls used: Three one-minute control counts were performed, without testing sample, before and after each 5 specimens. Control counts are averaged.

Testing date: 09/12/2020

Deviation from the stated method: Test aerosol is non-neutralized according to the FDA guidance document on surgical face masks.

Test results:

Specimen	Filtration Efficiency (%)
1	99,88
2	99,89
3	99,89
4	99,87
5	99,87
Average	99,88
Standard Deviation	0,01