

# COVID-19 Antigen Saliva Test Kit SELF TEST



RRP \$24.95



RRP \$62.50



**SELF TEST**  
**ARTG LISTED**  
**372335**



## *Avoid the discomfort of taking nasopharyngeal sample*

The COVID-19 Antigen Rapid Test Pen is an in vitro immunoassay. The assay is for the direct and qualitative detection of SARS-CoV-2 viral nucleoprotein antigens from saliva samples



## Benefits:

- Easy saliva sample – no nasal swab
- Simple 2 step process
- Test takes 2 mins
- Results in 15 mins

# Contents



## 2 Pack

- 2 x Saliva Antigen Tests
- 2 x Disposal bags
- 1 x Instructions for Use

## 5 Pack

- 5 x Saliva Antigen Tests
- 5 x Disposal bags
- 1 x Instructions for Use

## Test procedure



## Clinical Evaluation

Table: COVID-19 Rapid Test vs. RT-PCR

		RT-PCR		Total
		Positive	Negative	
COVID-19 Antigen Saliva Test Kit	Positive	40	2	42
	Negative	5	137	142
Total		45	139	184

Relative Sensitivity: 88.9% (76.5% ~ 95.2%)\*

Relative Specificity: 98.6% (94.9% ~ 99.6%)\*

Overall Agreement: 96.2% (92.4% ~ 98.1%)\*

\*95% Confidence Interval

### SELF TEST BARCODES BY PACK FORMAT

**2 pack box**  
200mm x 55mm x 30mm  
12 x 2 packs per display inner box

9359051000035



**5 pack box**  
200mm x 75mm x 50mm  
6 x 5 packs per display inner box

9359051000066



**2 pack display inner box**  
205mm x 166mm x 122mm  
16 inners per carton

9359051000042



**5 pack display inner box**  
205mm x 151mm x 152mm  
12 inners per carton

9359051000073



**2 pack box carton**  
520mm x 410mm x 36mm  
20 cartons per pallet  
16kg per carton

9359051000059



**5 pack box carton**  
630mm x 370mm x 30mm  
20 cartons per pallet  
15 kg per carton

9359051000080

