



# Ultra Health Medical Pty Ltd | Steeldrill Safety — Product Specification Sheet

<b>Name</b>	Ultra Health Disposable Non Woven Crimped Cap 21 inch (White, Yellow, Red, Green, Blue) 										
<b>Product Code</b>	200010/B, 200010/W, 200010/Y, 200010/R, 200010/G										
<b>Applicable Directive &amp; Product Standards</b>	ISO 13485:2016 Certificate: SX 601277 0001 TGA : 157453 HACCP certificate: PE-610-SWG-05 										
<b>Product Description</b>	Lightweight crimped caps are made from breathable, non-toxic, spunbonded polypropylene, double elastic										
<b>Quality Assurance</b>	EC certificate: DD 60127882 0001										
<b>Typical Uses</b>	Food Processing, Laboratory Use, Hospitality ect.										
<b>Weight</b>	14gsm										
<b>Measurements</b>	21 Inch										
<b>Code/Colour</b>	<table border="1"> <tr> <td>200010/W</td> <td>White</td> </tr> <tr> <td>200010/B</td> <td>Blue</td> </tr> <tr> <td>200010/R</td> <td>Red</td> </tr> <tr> <td>200010/G</td> <td>Green</td> </tr> <tr> <td>200010/Y</td> <td>Yellow</td> </tr> </table>	200010/W	White	200010/B	Blue	200010/R	Red	200010/G	Green	200010/Y	Yellow
200010/W	White										
200010/B	Blue										
200010/R	Red										
200010/G	Green										
200010/Y	Yellow										
<b>Shelf Life</b>	5 Years (valid from manufacturing date)										
<b>Packing</b>	100 pieces/pack, 1000 pieces/carton										
<b>Storage Conditions</b>	Temperature: Between -15°C & 40°C Humidity: Below 80% Sunlight: Do not expose to direct sunlight										
<b>Sterilisation</b>	Nil										

