


Ultra Health Medical Pty Ltd | Steeldrill Safety — Product Specification Sheet

| | | |
|---|--|---|
| Product | Ultra Health Surgical Caps Elastic White |  |
| Product Code | 200070/1000, 200070/500 | |
| Applicable Directive & Product Standards | ISO 13485:2016 Certificate: SX 601277 0001 | |
| Product Description | Disposable white surgical cap with elastic. Polypropylene, non-linting | |
| Colour | White | |
| TGA Number | 293266 | |
| Quality Assurance | EC certificate: DD 60127882 0001 | |
| Composition & Weight | Polypropylene, 20gsm | |
| Pack Quantity | 100 | |
| Carton Quantity | 200070/1000 = 1000 pieces, 200070/500= 500 pieces | |
| Shelf Life | 5 years from date of manufacture | |
| Storage | Keep in a cool, dry place out of direct sunlight | |

