

QSi EZY-Breathe P2 Respirator, Duckbill Style - Size Regular

INFO SHEET

PRODUCT CODE: 65972

PRODUCT NAME: QSi EZYBreathe P2 Respirator, Duckbill Style - Size Regular
(Box = 50 Respirators / Carton = 20 Boxes)

EZYBreathe™



INTENDED USE	P2 respirator.
DESCRIPTION	<p>The QSi® EZYBreathe P2 Respirator is a single-use, disposable 'duckbill' style respirator. It provides respiratory protection for the wearer from the transfer of micro-organisms, body fluids, and particulate material.</p> <p>This P2 respirator is certified to AS/NZS 1716:2012 which is the preferred standard for New Zealand and Australia. This standard is the local equivalent to the American N95 standard.</p>
APPROVALS	<ul style="list-style-type: none"> AS/NZS 1716:2012 BSI Certified Product - BMP 692299 AS 4381:2015 test method ISO 22609 WAND: 211027-WAND-6XUAAH
FEATURES	<ul style="list-style-type: none"> Flat-fold design Tight fitting respirator Able to filter out very fine particles when worn correctly EZYBreathe formula for reduced breathing resistance Fiberglass and latex free Plastic coated metal nose wire Green elastic straps to denote regular size
SPECIFICATIONS	<ul style="list-style-type: none"> High level of particulate filtration (as per standard) Resistance to penetration by synthetic blood to 160mmHg (Level 3 fluid resistance)
SUITABLE APPLICATIONS	<ul style="list-style-type: none"> Hospitals Medical centres Aged care Dentists
LIMITATIONS	<p>Single use only. Discard the respirator after use. Store in a cool dry place out of direct sunlight. Refer to AS/NZS 1715:2009 for selection, use and maintenance. Perform user seal checks every time after donning. Ensure fit testing is performed annually.</p>